

COALITION AGX[®]

Stand-Alone ACDF System



Our mission is to deliver cutting-edge technology, research, and innovative solutions to promote healing in patients with musculoskeletal disorders.

Life moves us 

The Surgical Technique shown is for illustrative purposes only. The technique(s) actually employed in each case always depends on the medical judgment of the surgeon exercised before and during surgery as to the best mode of treatment for each patient. Additionally, as instruments may occasionally be updated, the instruments depicted in this Surgical Technique may not be exactly the same as the instruments currently available. Please consult with your sales representative or contact Globus directly for more information.

SURGICAL TECHNIQUE GUIDE

COALITION AGX[®]

Implant Overview	6
Instrument Overview	8
Surgical Technique	
1. Approach and Disc Preparation	18
2. Distraction	19
3. Discectomy/Endplate Preparation	21
4. Implant Sizing	21
<i>Assembling the Plate and Spacer</i>	22
5. Implant Insertion	23
<i>Assembling the Implant/Drill Guide Holder</i>	24
<i>Assembling Angled Instruments</i>	25
6. Screw Hole Preparation	27
<i>Using the Drill Sleeve Adjuster</i>	27
<i>Using the Lateral Holder</i>	28
<i>Using the Midline Implant Holder</i>	28
<i>Aligning the Self-Centering Sleeve</i>	28
7. Screw Insertion	29
8. Screw Blocking	30
Final Position	31
Final Position (Radiographs)	31
Additional Specifications	32
Screw Depths	33
COALITION AGX [®] Instrument Set	35
COALITION AGX [®] Plate Sets	36
COALITION AGX [®] 12x14 and 14x16mm Radiolucent Spacer Set	37
COALITION AGX [®] 15x18 and 15x20mm Radiolucent Spacer Set	38
COALITION AGX [®] 12x14 and 14x16mm Modular Trial Set	39
COALITION AGX [®] 15x18 and 15x20mm Modular Trial Set	40
COALITION AGX [®] Allograft Set	41
Universal ACDF Instrument Set	42
Universal ACDF Rasp Set	43
Universal ACDF Screw Set	47
COALITION AGX [®] RP Instrument Set	48
COALITION AGX [®] RP Plate Sets	49
COALITION AGX [®] RP Modular Trial Set	50
Important Information	52
Allograft Tissue Safety and Product Information	55


Stand-Alone ACDF System

The design allows for ease of placement and fixation and is ideal for adjacent level fusion.

- Lower anterior profile compared to traditional plate-spacers
- Less disruption to adjacent level
- Requires minimal retraction


- Provides easy assembly in the OR
- Permits insertion/fixation of the spacer and plate
- Allows option of allograft or PEEK spacer

- Creates an osteoconductive scaffold that allows gradual bone incorporation




COALITION AGX®

- Anatomical profile plate, screw, and spacer system
- PEEK or allograft interbody spacer options
- Allograft spacer footprints: 12x14 and 14x16mm
- PEEK spacer footprints: 12x14, 14x16, and 15x18mm
- Eight heights: 5–12mm, in 1mm increments
- 7° lordotic option
- Torsional stabilization features
- Pre-assembled blocking mechanism
- Audible, tactile, and visual confirmation of screw blocking
- Large axial graft chamber




Allograft Spacer



PEEK Spacer


Screw Options

- **Diameter:** 3.6mm and 4.2mm
- **Lengths:** 12–20mm, in 2mm increments
- **Type:** Self-drilling or self-tapping
- **Variable and fixed angle**



Screw Angulation

- **Screw orientation:** 35° cephalad/caudal
- **Variable angle screws:** $\pm 4^\circ$
- **Medial angulation:**
 - 16° for 12x14mm
 - 12° for 14x16mm




Fixed Angle Screw

Variable Angle Screw

IMPLANT OVERVIEW


COALITION AGX® (Allograft with Fixation) Spacer

- Cortical allograft bone
- **Footprints:** 8x14, 10x16, and 10.5x18mm
- **Heights:** 5-10mm, in 1mm increments
- **Sagittal profiles:** 7° lordosis




COALITION AGX® PEEK Spacer

- PEEK polymer material
- **Footprints:** 8x14, 10x16, 10.5x18, and 10.5x20mm
- **Heights:** 5-12mm, in 1mm increments
- **Sagittal profiles:** 0° and 7° lordosis




COALITION AGX® RP Plate

- **Widths:** 14, 16, 18, and 20mm
- **Heights:** 5-12mm, in 1mm increments
- **Screw blocking:** Audible, tactile, and visual confirmation



COALITION AGX® Plate

- **Widths:** 14, 16, 18, and 20mm
- **Heights:** 5-12mm, in 1mm increments
- **Screw blocking:** Audible, tactile, and visual confirmation



INSTRUMENT OVERVIEW

DISTRACTION INSTRUMENTS



Distractor, Left 601.020







Distractor Locking Nuts 665.606



Distractor, Right 601.021

Distractor Pins

	Height	Part No.
	12mm	665.612
	14mm	665.614
	16mm	665.616
	18mm	665.618



Distractor Pin Driver 665.608

COALITION AGX® MODULAR TRIALS - FOR USE WITH STANDARD COALITION AGX® PLATES



Trial Holder, Modular Trial/Rasp Heads - Inner Shaft Assembly



Trial Holder, Modular Trial/Rasp Heads - Outer Sleeve 6147.9001



12x14mm		
Size	0°	7°
5mm	6128.1205	6128.1275
6mm	6128.1206	6128.1276
7mm	6128.1207	6128.1277
8mm	6128.1208	6128.1278
9mm	6128.1209	6128.1279
10mm	6128.1210	6128.1280
11mm	6128.1211	6128.1281
12mm	6128.1212	6128.1282

14x16mm		
Size	0°	7°
5mm	6128.1405	6128.1475
6mm	6128.1406	6128.1476
7mm	6128.1407	6128.1477
8mm	6128.1408	6128.1478
9mm	6128.1409	6128.1479
10mm	6128.1410	6128.1480
11mm	6128.1411	6128.1481
12mm	6128.1412	6128.1482



15x18mm		
Size	0°	7°
5mm	6128.1505	6128.1575
6mm	6128.1506	6128.1576
7mm	6128.1507	6128.1577
8mm	6128.1508	6128.1578
9mm	6128.1509	6128.1579
10mm	6128.1510	6128.1580
11mm	6128.1511	6128.1581
12mm	6128.1512	6128.1582

COALITION AGX® MODULAR TRIALS - FOR USE WITH COALITION AGX® RP PLATE



12x14mm		
Size	0°	7°
5mm	6128.2205	6128.2275
6mm	6128.2206	6128.2276
7mm	6128.2207	6128.2277
8mm	6128.2208	6128.2278
9mm	6128.2209	6128.2279
10mm	6128.2210	6128.2280
11mm	6128.2211	6128.2281
12mm	6128.2212	6128.2282



14x16mm		
Size	0°	7°
5mm	6128.2405	6128.2475
6mm	6128.2406	6128.2476
7mm	6128.2407	6128.2477
8mm	6128.2408	6128.2478
9mm	6128.2409	6128.2479
10mm	6128.2410	6128.2480
11mm	6128.2411	6128.2481
12mm	6128.2412	6128.2482



15x18mm		
Size	0°	7°
5mm	6128.2505	6128.2575
6mm	6128.2506	6128.2576
7mm	6128.2507	6128.2577
8mm	6128.2508	6128.2578
9mm	6128.2509	6128.2579
10mm	6128.2510	6128.2580
11mm	6128.2511	6128.2581
12mm	6128.2512	6128.2582




15x20mm		
Size	0°	7°
5mm	6128.2525	6128.2535
6mm	6128.2526	6128.2536
7mm	6128.2527	6128.2537
8mm	6128.2528	6128.2538
9mm	6128.2529	6128.2539
10mm	6128.2530	6128.2540
11mm	6128.2531	6128.2541
12mm	6128.2532	6128.2542

UNIVERSAL ACDF RASPS



12x14mm		
Size	0°	7°
5mm	6147.1025	-
6mm	6147.1026	6147.1126
7mm	6147.1027	6147.1127
8mm	6147.1028	6147.1128
9mm	6147.1029	6147.1129
10mm	6147.1030	6147.1130
11mm	6147.1031	6147.1131
12mm	6147.1032	6147.1132



14x16mm		
Size	0°	7°
5mm	6147.2025	-
6mm	6147.2026	6147.2126
7mm	6147.2027	6147.2127
8mm	6147.2028	6147.2128
9mm	6147.2029	6147.2129
10mm	6147.2030	6147.2130
11mm	6147.2031	6147.2131
12mm	6147.2032	6147.2132



15x18mm		
Size	0°	7°
5mm	6147.3025	6147.3125
6mm	6147.3026	6147.3126
7mm	6147.3027	6147.3127
8mm	6147.3028	6147.3128
9mm	6147.3029	6147.3129
10mm	6147.3030	6147.3130
11mm	6147.3031	6147.3131
12mm	6147.3032	6147.3132

IMPLANT INSERTION INSTRUMENTS




Inner Shaft



Outer Shaft



Implant/Drill Guide Holder 684.001
(Assembled)



Lateral Holder 684.400



Impactor 665.607

IMPLANT HOLDER TIPS - FOR USE WITH STANDARD COALITION AGX® PLATES


DTS Guide, Pre-Set Angle

Height	12x14mm	14x16mm	15x18mm
5mm	6128.4205	6128.4405	6128.4505
6mm	6128.4206	6128.4406	6128.4506
7mm	6128.4207	6128.4407	6128.4507
8mm	6128.4208	6128.4408	6128.4508
9mm	6128.4209	6128.4409	6128.4509
10mm	6128.4210	6128.4410	6128.4510
11mm	6128.4211	6128.4411	6128.4511
12mm	6128.4212	6128.4412	6128.4512



Midline Holder Tips

Gold (Small) 12x14mm and 14x16mm		Silver (Large) 15x18mm	
Height	Part No.	Height	Part No.
5mm	6128.4005	5mm	6128.4105
6mm	6128.4006	6mm	6128.4106
7mm	6128.4007	7mm	6128.4107
8mm	6128.4008	8mm	6128.4108
9mm	6128.4009	9mm	6128.4109
10mm	6128.4010	10mm	6128.4110
11mm	6128.4011	11mm	6128.4111
12mm	6128.4012	12mm	6128.4112



IMPLANT HOLDER TIPS - FOR USE WITH COALITION AGX® RP PLATE

DTS Guide, Pre-Set Angle

Height	12x14mm	14x16mm	15x18mm
5mm	684.155	684.165	684.175
6mm	684.156	684.166	684.176
7mm	684.157	684.167	684.177
8mm	684.158	684.168	684.178
9mm	684.159	684.169	684.179
10mm	684.160	684.170	684.180
11mm	684.161	684.171	684.181
12mm	684.162	684.172	684.182



Midline Holder Tips

Height	Part No.
5mm	-
6mm	684.406
7mm	684.407
8mm	684.408
9mm	684.409
10mm	684.410
11mm	684.411
12mm	684.412



Free Hand Gripping Tips



12x14mm
684.113



14x16mm
684.123



15x18mm
684.133

STRAIGHT INSTRUMENTS




Quick-Connect Handle, Swivel 636.450



Awl, Straight with Self-Centering Sleeve 684.403



Tap, Straight 684.004



Awl with Retractable Sleeve 684.006



Drills, Straight with Self-Centering Sleeve

Lengths	Part No.
12mm	684.422
14mm	684.424
16mm	684.426
18mm	684.428
20mm	684.430



Screwdriver, 2.5mm Hex, Self-Retaining, with Cap 684.305

ANGLED INSTRUMENTS



Awl, Bent with Self-Centering Sleeve 684.404



Angled Driver (Sleeve) Body 684.415



Drill Guide 6128.3000



Angled Driver Nut 684.416



Counter-Torque 684.421



Angled Driver Shaft 684.417

Angled Tap, Driver, and Drills

Instrument	Part No.
Angled Tap	684.416
Angled Driver	684.418
Angled Drill, 12mm	684.432
Angled Drill, 14mm	684.434
Angled Drill, 16mm	684.436
Angled Drill, 18mm	684.438
Angled Drill, 20mm	684.440



SET SCREW INSTRUMENTS



Set Screw Positioner, 2.0mm Hex, Torque Limiting 650.312



Screwdriver, 2.1mm Hex, QC 671.313

PLATE AND SPACER ASSEMBLY



Assembly Tool 6128.3102 Upper



Assembly Tool 6128.3101 Lower


SURGICAL TECHNIQUE

COALITION AGX[®]

Refer to the device insert (also printed at the back of this manual) for the complete description, indications, contraindications, and warnings.

STEP 1 APPROACH AND DISC PREPARATION


The patient is placed under anesthesia and positioned with support of the posterior cervical spine to maintain cervical lordosis. Traditional cervical retractors may be used. The operative area is carefully cleaned and an incision is made at the appropriate fusion level. COALITION AGX[®] implants may be used in the cervical spine from C2 to T1.



STEP 2 DISTRACTION

Distraction may be accomplished using the **Distractor** (available in the Universal ACDF System or COLONIAL® System).

To use the Distractor, first determine pin placement within the vertebral bodies. Select the appropriate pin length and place the **Distractor Pins** into adjacent vertebral bodies using the **Distractor Pin Driver**. Care should be taken in placing pins to avoid interference with COALITION AGX® implant insertion.



Inserting Distractor Pins



Distractor Pins inserted

DISTRACTION (CONT'D)

Place the Distractor (Right or Left) over pins until seated. Secure the Distractor in place by attaching the **Locking Nuts** and rotating clockwise. Rotate the ratchet handle to distract to the desired amount, being careful not to over-distract the segment. This method may be used throughout the technique to provide visualization and access to the disc and osseous structures.

Note: The Distractor is removed for visual clarity in the remaining images of this technique guide.



Distracting disc space
using Distractor

STEP

3

DISCECTOMY/ENDPLATE PREPARATION

Leaving the lateral annulus intact, remove the intervertebral disc and osteophytes using rongeurs and other preparation instruments. The Box Curette may be used to remove the disc as well as superficial layers of the cartilaginous endplates to expose bleeding bone. Alternatively, **Universal ACDF Rasps** match the trial design and may be used to expose bleeding bone.



Endplate preparation
using Box Curette



Endplate preparation alternatively
using Universal ACDF Rasp

STEP

4


IMPLANT SIZING

Determine the appropriate implant profile for the desired cervical segment. Insert the smallest **COALITION AGX® Trial** into the disc space first, moving to larger trials as needed. Determine which trial best fits the prepared disc space. A secure fit is desirable in order to maintain disc height and stabilize the segment. This can be confirmed using fluoroscopy and tactile feedback.

Plate and Spacer Sizing Matrix		
Titanium Plate	Allograft or PEEK Spacer	Assembled
14mm width	8x14mm	12x14mm footprint
16mm width	10x16mm	14x16mm footprint
18mm width	10.5x18mm	15x18mm footprint
20mm width*	10.5x20mm	15x20mm footprint

*Only PEEK Spacer

IMPLANT SIZING (CONT'D)



Implant sizing using
Universal ACDF Trial




ASSEMBLING THE PLATE AND SPACER

The allograft must be thoroughly hydrated prior to assembly. Hydration in 0.9% saline solution or the patient's own blood for at least 1 minute is recommended.



Choose the appropriate footprint and slide the spacer onto the gold track with the opening facing inward.



Insert the matching plate onto the silver track with the open side facing outward. Press down gently on the plate until it is engaged with the spacer.



Align to the correct size and place the lid onto the base. Press down to assemble the plate and spacer as shown. Visually ensure the plate and spacer are fully assembled.



Implant fully assembled

STEP**5****IMPLANT INSERTION**

The **Midline Implant Holder**, **Implant/Drill Guide Insertion Holder**, **Lateral Implant Holder**, **Free Hand Gripping Tip**, or **DTS Guide, Pre-Set Angle** may be used for insertion.

Assemble the insertion holder and attach the appropriate sized DTS Guide, Pre-Set Angle, Free Hand Gripping Tip, or Midline Implant Holder.

Load the implant on the selected implant holder and insert into the intervertebral space. Light impaction may be used.



Insertion with
Midline Implant Holder



Insertion with
DTS Guide, Pre-Set Angle


⚙️ ASSEMBLING THE IMPLANT/DRILL GUIDE HOLDER

The Implant/Drill Guide Holder has two components. To assemble, slide the outer sleeve over the inner shaft and thread by rotating the inner shaft clockwise, as shown below. Continue to rotate until contact is made between the inner shaft and outer sleeve.


Note: Do not continue to rotate the inner shaft after initial contact.



Connect the Midline Implant Holder by rotating the handle of the inner shaft clockwise until secure.



Prior to loading the implant, confirm that the outer and inner shafts are touching. Rotate the outer shaft counterclockwise to spread the tips. Attach the implant onto the holder. The implant is held lightly but is not yet secured.




Attach implant onto Midline Holding Tip




Implant loading position

To secure the implant to the holder, rotate the outer sleeve clockwise until it stops. The implant is now secured to the holder.



Implant secured into Midline Holding Tip




Secured position

To release the implant, hold the inner shaft and rotate the outer sleeve counterclockwise until the inner shaft and outer sleeve touch. Continue to rotate the outer sleeve counterclockwise until the Midline Holding Tip releases the implant.



Implant releases from Midline Holding Tip



Disengaged position

ASSEMBLING ANGLED INSTRUMENTS

Select the appropriate Angled Driver, Angled Drill, or Angled Tap for assembly with the body and shaft.

Hold the **Angled Driver Body** pointed downward with the distal cutout facing upward. Insert the selected tip into the cutout.




Insert the **Angled Driver Shaft** into the driver body until the gears on the shaft and the selected tip mate.



Place the **Backing Nut** over the shaft. Rotate the threads clockwise until the nut sits flush with the driver body.




Attach a **Quick-Connect Handle, Swivel**. The driver is now ready for use.




ASSEMBLING ANGLED INSTRUMENTS (CONT'D)

For additional control, a **Counter-Torque** handle may be attached to the angled instruments.

Starting from the top, slide the Counter-Torque from the smooth portion of the Angled Driver Body to the knurled portion until fully seated.




Attaching Counter-Torque



Counter-Torque in final position

Rotate the Counter-Torque clockwise to final tighten.



Using the Pre-Set Angled Drill Guide

Insert the **Awl with Self-Centering Sleeve, Bent**, through the screw hole within the plate and break the cortex.

Determine the desired depth and select the appropriate fixed length drill. If desired, a depth gauge may be used to determine the appropriate depth. Drills are color-coded by screw length. Insert the drill into the drill guide and drill to the stop.

Screw holes may be tapped using the Straight or Angled Tap.



Using Awl with Self-Centering Sleeve, Bent



Using Angled Drill with Self-Centering Sleeve

Note: The same screw hole preparation instruments are used for the Implant/Drill Guide Holder, Lateral Holder, and Midline Implant Holder.



USING THE DRILL SLEEVE ADJUSTER

Ensure that the drill sleeve is in the proper position by first gripping the drill sleeve with the **Drill Sleeve Adjuster**, then pulling the sleeve distally.




Gripping drill sleeve




Pulling sleeve distally

⚙️ USING THE LATERAL HOLDER



Fully open the Lateral Holder by rotating the locking nut counterclockwise.



Place the appropriately sized implant on the Lateral Holder, compress the handles, and rotate the locking nut clockwise until finger tight.

⚙️ USING THE MIDLINE IMPLANT HOLDER




Insert the Awl with Self-Centering Sleeve, Bent to break the cortex.




A **Self-Centering Drill** and **Tap** may be used to further prepare the screw hole. Depending on the angle and position, a straight or angled instrument may be used.

⚙️ ALIGNING THE SELF-CENTERING SLEEVE

The Self-Centering Sleeves ensure proper screw trajectory without the use of a drill guide. The sleeve must be properly engaged with the plate before advancing any screw hole preparation instrument. Proceed to screw insertion prior to preparing the remaining screw hole.



INCORRECT



CORRECT

STEP


7

SCREW INSERTION

Load the desired screw from the module using the **Screwdriver, 2.5mm Hex, Self-Retaining, with Cap** or **Angled Driver**. Confirm screw length and diameter using the gauges within the screw module. Insert the screw through the screw hole. As the screw is inserted, the implant lags to the bone. Repeat for the second screw.



Screw insertion with DTS Guide, Pre-Set Angle



Screw insertion with Midline Implant Holder



STEP

8

SCREW BLOCKING

Once the screws are fully seated within the plate, insert the **Set Screw Positioner, 2.0mm Hex, 0.3Nm Torque Limiting** into the blocking set screw and rotate clockwise approximately 90°. The positioner provides audible and tactile confirmation that the screw is blocked from backing out.



If the blocking set screw does not rotate, check screw angulation and ensure that screws are fully seated.



Blocking set
screw unlocked

Initial Unlocked Position

Blocking set screw *facing* bone screw




Blocking set
screw locked


Final Locked Position

Blocking set screw *facing away* from bone screw

FINAL POSITION



FINAL POSITION (RADIOGRAPHS)











AP final position
(allograft spacer)















Lateral final position
(allograft spacer)













ADDITIONAL SPECIFICATIONS

Instrument	Diameter (mm)	Depth (mm)	Design
Straight/Bent Awl	2.4	8.6	
Retractable Awl	2.4	8.6	
Straight/Angled Taps	3.5	10.4	
12mm Drill	2.4	8.6	
14mm Drill	2.4	10.6	
16mm Drill	2.4	12.6	
18mm Drill	2.4	14.6	
20mm Drill	2.4	16.6	

SCREW DEPTHS - STANDARD PLATE

	12mm Screws	14mm Screws	16mm Screws
12x14mm			
	14mm Screws	16mm Screws	18mm Screws
14x16mm			
	16mm Screws	18mm Screws	20mm Screws
15x18mm			
	16mm Screws	18mm Screws	20mm Screws
15x20mm			

SCREW DEPTHS - RP PLATE

	12mm Screws	14mm Screws	16mm Screws
12x14mm			
	14mm Screws	16mm Screws	18mm Screws
14x16mm			
	16mm Screws	18mm Screws	20mm Screws
15x18mm			
	16mm Screws	18mm Screws	20mm Screws
15x20mm			

COALITION AGX®

INSTRUMENT SET 9128.9001

AGX DTS Guides, Pre-Set Angle

Part No.	Description	Qty
6128.4205	14mm W, 5mm H	1
6128.4206	14mm W, 6mm H	1
6128.4207	14mm W, 7mm H	1
6128.4208	14mm W, 8mm H	1
6128.4209	14mm W, 9mm H	1
6128.4210	14mm W, 10mm H	1
6128.4211	14mm W, 11mm H	1
6128.4212	14mm W, 12mm H	1
6128.4405	16mm W, 5mm H	1
6128.4406	16mm W, 6mm H	1
6128.4407	16mm W, 7mm H	1
6128.4408	16mm W, 8mm H	1
6128.4409	16mm W, 9mm H	1
6128.4410	16mm W, 10mm H	1
6128.4411	16mm W, 11mm H	1
6128.4412	16mm W, 12mm H	1
6128.4505	18mm W, 5mm H	1
6128.4506	18mm W, 6mm H	1
6128.4507	18mm W, 7mm H	1
6128.4508	18mm W, 8mm H	1
6128.4509	18mm W, 9mm H	1
6128.4510	18mm W, 10mm H	1
6128.4511	18mm W, 11mm H	1
6128.4512	18mm W, 12mm H	1

Midline Holding Tip

Part No.	Length	Qty
6128.4005	5mm	1
6128.4006	6mm	1
6128.4007	7mm	1
6128.4008	8mm	1
6128.4009	9mm	1
6128.4010	10mm	1
6128.4011	11mm	1
6128.4012	12mm	1

Midline Holding Tip, Large

Part No.	Length	Qty
6128.4105	5mm	1
6128.4106	6mm	1
6128.4107	7mm	1
6128.4108	8mm	1
6128.4109	9mm	1
6128.4110	10mm	1
6128.4111	11mm	1
6128.4112	12mm	1

Other Instruments

Part No.	Description	Qty
650.312	Set Screw Positioner, 2.0mm Hex, Torque-Limiting	1
684.001	Implant/Drill Guide Holder	1
684.113	Free Hand Gripping Tip, 12x14mm	1
684.123	Free Hand Gripping Tip, 14x16mm	1
684.133	Free Hand Gripping Tip, 15x18mm	1
684.400	Lateral Holder	1
6128.3000	Handheld Drill Guide	1
6128.3101	Assembly Press, Lower	1
6128.3102	Assembly Press, Upper	1
9128.0001	COALITION AGX® Graphic Case	

COALITION AGX®

PLATE SETS

COALITION AGX® 12x14 and 14x16mm Plate Set 9128.9002

Part No.	Description	Qty
1128.1205	14mm W, 5mm H	2
1128.1206	14mm W, 6mm H	3
1128.1207	14mm W, 7mm H	3
1128.1208	14mm W, 8mm H	3
1128.1209	14mm W, 9mm H	2
1128.1210	14mm W, 10mm H	2
1128.1211	14mm W, 11mm H	1
1128.1212	14mm W, 12mm H	1
1128.1405	16mm W, 5mm H	2
1128.1406	16mm W, 6mm H	3
1128.1407	16mm W, 7mm H	3
1128.1408	16mm W, 8mm H	3
1128.1409	16mm W, 9mm H	2
1128.1410	16mm W, 10mm H	2
1128.1411	16mm W, 11mm H	1
1128.1412	16mm W, 12mm H	1

9128.0002 COALITION AGX® Plate Module,
12x14 and 14x16mm

COALITION AGX® 15x18 and 15x20mm Plate Set 9128.9003

Part No.	Description	Qty
1128.1505	18mm W, 5mm H	2
1128.1506	18mm W, 6mm H	3
1128.1507	18mm W, 7mm H	3
1128.1508	18mm W, 8mm H	3
1128.1509	18mm W, 9mm H	2
1128.1510	18mm W, 10mm H	2
1128.1511	18mm W, 11mm H	1
1128.1512	18mm W, 12mm H	1
1128.1525	20mm W, 5mm H	1
1128.1526	20mm W, 6mm H	1
1128.1527	20mm W, 7mm H	1
1128.1528	20mm W, 8mm H	1
1128.1529	20mm W, 9mm H	1
1128.1530	20mm W, 10mm H	1
1128.1531	20mm W, 11mm H	1
1128.1532	20mm W, 12mm H	1

9128.0003 COALITION AGX® Plate Module,
15x18 and 15x20mm

COALITION AGX® 12x14 and 14x16mm RADIOLUCENT SPACER SET 9128.9004

Part No.	Description	Qty
3128.1205	COALITION AGX® R Spacer, 8x14mm, 0°, 5mm	1
3128.1206	COALITION AGX® R Spacer, 8x14mm, 0°, 6mm	1
3128.1207	COALITION AGX® R Spacer, 8x14mm, 0°, 7mm	1
3128.1208	COALITION AGX® R Spacer, 8x14mm, 0°, 8mm	1
3128.1209	COALITION AGX® R Spacer, 8x14mm, 0°, 9mm	1
3128.1210	COALITION AGX® R Spacer, 8x14mm, 0°, 10mm	1
3128.1211	COALITION AGX® R Spacer, 8x14mm, 0°, 11mm	1
3128.1212	COALITION AGX® R Spacer, 8x14mm, 0°, 12mm	1
3128.1275	COALITION AGX® R Spacer, 8x14mm, 7°, 5mm	2
3128.1276	COALITION AGX® R Spacer, 8x14mm, 7°, 6mm	3
3128.1277	COALITION AGX® R Spacer, 8x14mm, 7°, 7mm	3
3128.1278	COALITION AGX® R Spacer, 8x14mm, 7°, 8mm	3
3128.1279	COALITION AGX® R Spacer, 8x14mm, 7°, 9mm	2
3128.1280	COALITION AGX® R Spacer, 8x14mm, 7°, 10mm	2
3128.1281	COALITION AGX® R Spacer, 8x14mm, 7°, 11mm	1
3128.1282	COALITION AGX® R Spacer, 8x14mm, 7°, 12mm	1
3128.1405	COALITION AGX® R Spacer, 10x16mm, 0°, 5mm	1
3128.1406	COALITION AGX® R Spacer, 10x16mm, 0°, 6mm	1
3128.1407	COALITION AGX® R Spacer, 10x16mm, 0°, 7mm	1
3128.1408	COALITION AGX® R Spacer, 10x16mm, 0°, 8mm	1
3128.1409	COALITION AGX® R Spacer, 10x16mm, 0°, 9mm	1
3128.1410	COALITION AGX® R Spacer, 10x16mm, 0°, 10mm	1
3128.1411	COALITION AGX® R Spacer, 10x16mm, 0°, 11mm	1
3128.1412	COALITION AGX® R Spacer, 10x16mm, 0°, 12mm	1
3128.1475	COALITION AGX® R Spacer, 10x16mm, 7°, 5mm	2
3128.1476	COALITION AGX® R Spacer, 10x16mm, 7°, 6mm	3
3128.1477	COALITION AGX® R Spacer, 10x16mm, 7°, 7mm	3
3128.1478	COALITION AGX® R Spacer, 10x16mm, 7°, 8mm	3
3128.1479	COALITION AGX® R Spacer, 10x16mm, 7°, 9mm	2
3128.1480	COALITION AGX® R Spacer, 10x16mm, 7°, 10mm	2
3128.1481	COALITION AGX® R Spacer, 10x16mm, 7°, 11mm	1
3128.1482	COALITION AGX® R Spacer, 10x16mm, 7°, 12mm	1
9128.0004	COALITION AGX® Module, R Spacers, 12x14 and 14x16mm	

COALITION AGX® 15x18 and 15x20mm RADIOLUCENT SPACER SET 9128.9005

Part No.	Description	Qty
3128.1505	COALITION AGX® R Spacer, 10.5x18mm, 0°, 5mm	1
3128.1506	COALITION AGX® R Spacer, 10.5x18mm 0°, 6mm	1
3128.1507	COALITION AGX® R Spacer, 10.5x18mm, 0°, 7mm	1
3128.1508	COALITION AGX® R Spacer, 10.5x18mm, 0°, 8mm	1
3128.1509	COALITION AGX® R Spacer, 10.5x18mm, 0°, 9mm	1
3128.1510	COALITION AGX® R Spacer, 10.5x18mm, 0°, 10mm	1
3128.1511	COALITION AGX® R Spacer, 10.5x18mm, 0°, 11mm	1
3128.1512	COALITION AGX® R Spacer, 10.5x18mm, 0°, 12mm	1
3128.1525	COALITION AGX® R Spacer, 10.5x20mm, 0°, 5mm	1
3128.1526	COALITION AGX® R Spacer, 10.5x20mm, 0°, 6mm	1
3128.1527	COALITION AGX® R Spacer, 10.5x20mm, 0°, 7mm	1
3128.1528	COALITION AGX® R Spacer, 10.5x20mm, 0°, 8mm	1
3128.1529	COALITION AGX® R Spacer, 10.5x20mm, 0°, 9mm	1
3128.1530	COALITION AGX® R Spacer, 10.5x20mm, 0°, 10mm	1
3128.1531	COALITION AGX® R Spacer, 10.5x20mm, 0°, 11mm	1
3128.1532	COALITION AGX® R Spacer, 10.5x20mm, 0°, 12mm	1
3128.1535	COALITION AGX® R Spacer, 10.5x20mm, 7°, 5mm	1
3128.1536	COALITION AGX® R Spacer, 10.5x20mm, 7°, 6mm	1
3128.1537	COALITION AGX® R Spacer, 10.5x20mm, 7°, 7mm	1
3128.1538	COALITION AGX® R Spacer, 10.5x20mm, 7°, 8mm	1
3128.1539	COALITION AGX® R Spacer, 10.5x20mm, 7°, 9mm	1
3128.1540	COALITION AGX® R Spacer, 10.5x20mm, 7°, 10mm	1
3128.1541	COALITION AGX® R Spacer, 10.5x20mm, 7°, 11mm	1
3128.1542	COALITION AGX® R Spacer, 10.5x20mm, 7°, 12mm	1
3128.1575	COALITION AGX® R Spacer, 10.5x18mm, 7°, 5mm	2
3128.1576	COALITION AGX® R Spacer, 10.5x18mm, 7°, 6mm	3
3128.1577	COALITION AGX® R Spacer, 10.5x18mm, 7°, 7mm	3
3128.1578	COALITION AGX® R Spacer, 10.5x18mm, 7°, 8mm	3
3128.1579	COALITION AGX® R Spacer, 10.5x18mm, 7°, 9mm	2
3128.1580	COALITION AGX® R Spacer, 10.5x18mm, 7°, 10mm	2
3128.1581	COALITION AGX® R Spacer, 10.5x18mm, 7°, 11mm	1
3128.1582	COALITION AGX® R Spacer, 10.5x18mm, 7°, 12mm	1
9128.0005	COALITION AGX® Module, R Spacers 15x18 and 15x20mm	

COALITION AGX® 12x14 and 14x16mm MODULAR TRIAL SET 9128.9006

COALITION AGX® 12x14mm, 0° Trials

Part No.	Length	Qty
6128.1205	5mm	1
6128.1206	6mm	1
6128.1207	7mm	1
6128.1208	8mm	1
6128.1209	9mm	1
6128.1210	10mm	1
6128.1211	11mm	1
6128.1212	12mm	1

COALITION AGX® 14x16mm, 0° Trials

Part No.	Length	Qty
6128.1405	5mm	1
6128.1406	6mm	1
6128.1407	7mm	1
6128.1408	8mm	1
6128.1409	9mm	1
6128.1410	10mm	1
6128.1411	11mm	1
6128.1412	12mm	1

COALITION AGX® 12x14mm, 7° Trials

Part No.	Length	Qty
6128.1275	5mm	1
6128.1276	6mm	1
6128.1277	7mm	1
6128.1278	8mm	1
6128.1279	9mm	1
6128.1280	10mm	1
6128.1281	11mm	1
6128.1282	12mm	1

COALITION AGX® 14x16mm, 7° Trials

Part No.	Length	Qty
6128.1475	5mm	1
6128.1476	6mm	1
6128.1477	7mm	1
6128.1478	8mm	1
6128.1479	9mm	1
6128.1480	10mm	1
6128.1481	11mm	1
6128.1482	12mm	1

9128.0006 AGX Modular Trial Module, 12x14 and 14x16mm

COALITION AGX® 15x18 and 15x20mm MODULAR TRIAL SET 9128.9007

COALITION AGX® 15x18, 0° Modular Trials

Part No.	Length	Qty
6128.1505	5mm	1
6128.1506	6mm	1
6128.1507	7mm	1
6128.1508	8mm	1
6128.1509	9mm	1
6128.1510	10mm	1
6128.1511	11mm	1
6128.1512	12mm	1

COALITION AGX® 15x20, 0° Modular Trials

Part No.	Length	Qty
6128.1525	5mm	1
6128.1526	6mm	1
6128.1527	7mm	1
6128.1528	8mm	1
6128.1529	9mm	1
6128.1530	10mm	1
6128.1531	11mm	1
6128.1532	12mm	1

COALITION AGX® 15x18, 7° Modular Trials

Part No.	Length	Qty
6128.1575	5mm	1
6128.1576	6mm	1
6128.1577	7mm	1
6128.1578	8mm	1
6128.1579	9mm	1
6128.1580	10mm	1
6128.1581	11mm	1
6128.1582	12mm	1

COALITION AGX® 15x20, 7° Modular Trials

Part No.	Length	Qty
6128.1535	5mm	1
6128.1536	6mm	1
6128.1537	7mm	1
6128.1538	8mm	1
6128.1539	9mm	1
6128.1540	10mm	1
6128.1541	11mm	1
6128.1542	12mm	1

9128.0007 AGX Modular Trial Module, 15x18 and 15x20mm

COALITION AGX®

ALLOGRAFT SET 9128.9010

Part No.	Description	Qty
8128.1275S	COALITION AGX® Allograft, 8x14mm, 7°, 5mm	1
8128.1276S	COALITION AGX® Allograft, 8x14mm, 7°, 6mm	3
8128.1277S	COALITION AGX® Allograft, 8x14mm, 7°, 7mm	3
8128.1278S	COALITION AGX® Allograft, 8x14mm, 7°, 8mm	3
8128.1279S	COALITION AGX® Allograft, 8x14mm, 7°, 9mm	2
8128.1280S	COALITION AGX® Allograft, 8x14mm, 7°, 10mm	1
8128.1475S	COALITION AGX® Allograft, 10x16mm, 7°, 5mm	1
8128.1476S	COALITION AGX® Allograft, 10x16mm, 7°, 6mm	3
8128.1477S	COALITION AGX® Allograft, 10x16mm, 7°, 7mm	3
8128.1478S	COALITION AGX® Allograft, 10x16mm, 7°, 8mm	3
8128.1479S	COALITION AGX® Allograft, 10x16mm, 7°, 9mm	2
8128.1480S	COALITION AGX® Allograft, 10x16mm, 7°, 10mm	1
8128.1575S	COALITION AGX® Allograft, 10.5x18mm, 7°, 5mm	2
8128.1576S	COALITION AGX® Allograft, 10.5x18mm, 7°, 6mm	2
8128.1577S	COALITION AGX® Allograft, 10.5x18mm, 7°, 7mm	2
8128.1578S	COALITION AGX® Allograft, 10.5x18mm, 7°, 8mm	2
8128.1579S	COALITION AGX® Allograft, 10.5x18mm, 7°, 9mm	1
8128.1580S	COALITION AGX® Allograft, 10.5x18mm, 7°, 10mm	1
9128.0010	COALITION AGX® Allograft Soft Case	

UNIVERSAL ACDF INSTRUMENT SET 9147.9001

Part No.	Description	Qty	Part No.	Description	Qty
601.020	Distractor, Left	1	684.421	Counter-Torque, Angled Instrument	2
601.021	Distractor, Right	1	684.422	Straight Drill with Self-Centering Sleeve, 12mm	1
6147.9001	Trial Holder, Modular Trial / Rasp Heads - Outer Sleeve	3	684.424	Straight Drill with Self-Centering Sleeve, 14mm	1
6147.9002	Trial Holder, Modular Trial / Rasp Heads - Inner Shaft Assembly	3	684.425	Angled Driver Tip, Short	0
6147.9003	Packing Block	1	684.426	Straight Drill with Self-Centering Sleeve, 16mm	1
636.450	Quick-Connect Handle, Swivel	2	684.428	Straight Drill with Self-Centering Sleeve, 18mm	1
665.504	Bone Packer	1	684.430	Straight Drill with Self-Centering Sleeve, 20mm	1
665.606	Distractor Locking Nuts	4	684.432	Angled Drill Tip with Self-Centering Sleeve, 12mm	1
665.607	Impactor	1	684.434	Angled Drill Tip with Self-Centering Sleeve, 14mm	1
665.608	Distractor Pin Driver	1	684.436	Angled Drill Tip with Self-Centering Sleeve, 16mm	1
665.610	Distractor Pin, 10mm	0	684.438	Angled Drill Tip with Self-Centering Sleeve, 18mm	1
665.612	Distractor Pin, 12mm	2	684.440	Angled Drill Tip with Self-Centering Sleeve, 20mm	1
665.614	Distractor Pin, 14mm	2	9147.0001	Universal ACDF Instrument Graphic Case	
665.616	Distractor Pin, 16mm	2	984.004	COALITION® Module, Angled Instruments	
665.618	Distractor Pin, 18mm	2			
671.313	VIP® Screwdriver, 2.1mm Hex, QC	1			
684.004	Tap, Straight	1			
684.006	Awl with Retractable Sleeve	1			
684.305	Screwdriver, 2.5mm Hex, Self-Retaining, with Cap	2			
684.309	Drill Sleeve Adjuster	1			
684.401	Self-Centering Sleeve - Short	2			
684.402	Self-Centering Sleeve - Long	2			
684.403	Awl with Self-Centering Sleeve, Straight	1			
684.404	Awl with Self-Centering Sleeve, Bent	1			
684.405	Sleeved Driver	0			
684.415	Angled Sleeve	2			
684.416	Angled Sleeve with Backing Nut	2			
684.417	Angled Driving Shaft	2			
684.418	Hex Driver Assembly	2			
684.419	Angled Tap Tip	1			

UNIVERSAL ACDF 12x14mm RASP SET 9147.9003

Universal ACDF 12x14mm, 0° Modular Rasps

Part No.	Length	Qty
6147.1025	5mm	1
6147.1026	6mm	1
6147.1027	7mm	1
6147.1028	8mm	1
6147.1029	9mm	1
6147.1030	10mm	1
6147.1031	11mm	1
6147.1032	12mm	1

Universal ACDF 12x14mm, 7° Modular Rasps

Part No.	Length	Qty
6147.1126	6mm	1
6147.1127	7mm	1
6147.1128	8mm	1
6147.1129	9mm	1
6147.1130	10mm	1
6147.1131	11mm	1
6147.1132	12mm	1

9147.0003 12x14mm Universal Rasps

UNIVERSAL ACDF 14x16mm RASP SET 9147.9004

Universal ACDF 14x16mm, 0° Modular Rasps

Part No.	Length	Qty
6147.2025	5mm	1
6147.2026	6mm	1
6147.2027	7mm	1
6147.2028	8mm	1
6147.2029	9mm	1
6147.2030	10mm	1
6147.2031	11mm	1
6147.2032	12mm	1

Universal ACDF 14x16mm, 7° Modular Rasps

Part No.	Length	Qty
6147.2126	6mm	1
6147.2127	7mm	1
6147.2128	8mm	1
6147.2129	9mm	1
6147.2130	10mm	1
6147.2131	11mm	1
6147.2132	12mm	1

9147.0004 14x16mm Universal Rasps

UNIVERSAL ACDF 15x18mm RASP SET 9147.9005

Universal ACDF 15x18mm, 0° Modular Rasps

Part No.	Length	Qty
6147.3025	5mm	1
6147.3026	6mm	1
6147.3027	7mm	1
6147.3028	8mm	1
6147.3029	9mm	1
6147.3030	10mm	1
6147.3031	11mm	1
6147.3032	12mm	1

Universal ACDF 15x18mm, 7° Modular Rasps

Part No.	Length	Qty
6147.3126	6mm	1
6147.3127	7mm	1
6147.3128	8mm	1
6147.3129	9mm	1
6147.3130	10mm	1
6147.3131	11mm	1
6147.3132	12mm	1

9147.0005 15x18mm Universal Rasps

UNIVERSAL ACDF 15x20mm RASP SET 9147.9006

Universal ACDF 15x20mm, 0° Modular Rasps

Part No.	Length	Qty
6147.4025	5mm	1
6147.4026	6mm	1
6147.4027	7mm	1
6147.4028	8mm	1
6147.4029	9mm	1
6147.4030	10mm	1
6147.4031	11mm	1
6147.4032	12mm	1

Universal ACDF 15x20mm, 7° Modular Rasps

Part No.	Length	Qty
6147.4126	6mm	1
6147.4127	7mm	1
6147.4128	8mm	1
6147.4129	9mm	1
6147.4130	10mm	1
6147.4131	11mm	1
6147.4132	12mm	1

9147.0006 15x20mm Universal Rasps

UNIVERSAL ACDF SCREW SET 9147.9008

4.2mm Bone Screw Variable, Self-Tapping

Part No.	Length	Qty
184.012	12mm	6
184.014	14mm	6
184.016	16mm	6
184.018	18mm	4
184.020	20mm	4

4.2mm Bone Screw Fixed, Self-Tapping

Part No.	Length	Qty
184.032	12mm	6
184.034	14mm	6
184.036	16mm	6
184.038	18mm	4
184.040	20mm	4

4.2mm Bone Screw Variable, Self-Drilling

Part No.	Length	Qty
184.052	12mm	6
184.054	14mm	6
184.056	16mm	6
184.058	18mm	4
184.060	20mm	4

4.2mm Bone Screw Fixed, Self-Drilling

Part No.	Length	Qty
184.072	12mm	6
184.074	14mm	6
184.076	16mm	6
184.078	18mm	4
184.080	20mm	4

3.6mm Bone Screw Variable, Self-Tapping

Part No.	Length	Qty
184.112	12mm	6
184.114	14mm	6
184.116	16mm	6
184.118	18mm	4
184.120	20mm	4

3.6mm Bone Screw Fixed, Self-Tapping

Part No.	Length	Qty
184.132	12mm	6
184.134	14mm	6
184.136	16mm	6
184.138	18mm	4
184.140	20mm	4

3.6mm Bone Screw Variable, Self-Drilling

Part No.	Length	Qty
184.152	12mm	6
184.154	14mm	6
184.156	16mm	6
184.158	18mm	4
184.160	20mm	4

3.6mm Bone Screw Fixed, Self-Drilling

Part No.	Length	Qty
184.172	12mm	6
184.174	14mm	6
184.176	16mm	6
184.178	18mm	4
184.180	20mm	4

9147.0008 Screw Module, Universal ACDF

COALITION AGX® RP

INSTRUMENT SET 9128.9021

DTS Guides, 12x14mm

Part No.	Length	Qty
684.155	5mm	1
684.156	6mm	1
684.157	7mm	1
684.158	8mm	1
684.159	9mm	1
684.160	10mm	1
684.161	11mm	1
684.162	12mm	1
684.165	5mm	1
684.166	6mm	1
684.167	7mm	1
684.166	6mm	1
684.167	7mm	1
684.168	8mm	1
684.169	9mm	1
684.170	10mm	1
684.171	11mm	1
684.172	12mm	1
684.175	5mm	1
684.176	6mm	1
684.177	7mm	1
684.178	8mm	1
684.179	9mm	1
684.180	10mm	1
684.181	11mm	1
684.182	12mm	1

Midline Holding Tip

Part No.	Length	Qty
684.406	6mm	1
684.407	7mm	1
684.408	8mm	1
684.409	9mm	1
684.410	10mm	1
684.411	11mm	1
684.412	12mm	1

Free Hand Gripping Tip

Part No.	Description	Qty
684.113	12x14mm	1
684.123	14x16mm	1
684.133	15x18mm	1

Other Instruments

Part No.	Description	Qty
6128.3000	Hand Held Drill Guide	1
6128.3101	COALITION AGX® Assembly Press, Lower	1
6128.3102	COALITION AGX® Assembly Press, Upper	1
650.312	Set Screw Positioner, 2.0mm Hex, Torque-Limiting	1
684.001	Implant/Drill Guide Insertion Holder	1
684.400	Lateral Holder	1

9128.0021 COALITION AGX® RP Graphic Case

COALITION AGX® RP PLATE SETS

COALITION AGX® 12x14 and 14x16mm RP Plate Sets 9128.9022

Part No.	Description	Qty
1128.2205	14mm W, 5mm H	2
1128.2206	14mm W, 6mm H	3
1128.2207	14mm W, 7mm H	3
1128.2208	14mm W, 8mm H	3
1128.2209	14mm W, 9mm H	2
1128.2210	14mm W, 10mm H	2
1128.2211	14mm W, 11mm H	1
1128.2212	14mm W, 12mm H	1
1128.2405	16mm W, 5mm H	2
1128.2406	16mm W, 6mm H	3
1128.2407	16mm W, 7mm H	3
1128.2408	16mm W, 8mm H	3
1128.2409	16mm W, 9mm H	2
1128.2410	16mm W, 10mm H	2
1128.2411	16mm W, 11mm H	1
1128.2412	16mm W, 12mm H	1
9128.0022	COALITION AGX® RP Plate Module, 12x14 and 14x16mm	

COALITION AGX® 15x18 and 15x20mm RP Plate Sets 9128.9023

Part No.	Description	Qty
1128.2505	18mm W, 5mm H	2
1128.2506	18mm W, 6mm H	3
1128.2507	18mm W, 7mm H	3
1128.2508	18mm W, 8mm H	3
1128.2509	18mm W, 9mm H	2
1128.2510	18mm W, 10mm H	2
1128.2511	18mm W, 11mm H	1
1128.2512	18mm W, 12mm H	1
1128.2525	20mm W, 5mm H	1
1128.2526	20mm W, 6mm H	1
1128.2527	20mm W, 7mm H	1
1128.2528	20mm W, 8mm H	1
1128.2529	20mm W, 9mm H	1
1128.2530	20mm W, 10mm H	1
1128.2531	20mm W, 11mm H	1
1128.2532	20mm W, 12mm H	1
9128.0023	COALITION AGX® RP Plate Module, 15x18 and 15x20mm	

COALITION AGX® RP 12x14 and 14x16mm MODULAR TRIAL SET 9128.9026

COALITION AGX® RP 12x14, 0° Modular Trials

Part No.	Length	Qty
6128.2205	5mm H	1
6128.2206	6mm H	1
6128.2207	7mm H	1
6128.2208	8mm H	1
6128.2209	9mm H	1
6128.2210	10mm H	1
6128.2211	11mm H	1
6128.2212	12mm H	1

COALITION AGX® RP 14x16, 0° Modular Trials

Part No.	Length	Qty
6128.2405	5mm H	1
6128.2406	6mm H	1
6128.2407	7mm H	1
6128.2408	8mm H	1
6128.2409	9mm H	1
6128.2410	10mm H	1
6128.2411	11mm H	1
6128.2412	12mm H	1

COALITION AGX® RP 12x14, 7° Modular Trials

Part No.	Length	Qty
6128.2275	5mm H	1
6128.2276	6mm H	1
6128.2277	7mm H	1
6128.2278	8mm H	1
6128.2279	9mm H	1
6128.2280	10mm H	1
6128.2281	11mm H	1
6128.2282	12mm H	1

COALITION AGX® RP 14x16, 7° Modular Trials

Part No.	Length	Qty
6128.2475	5mm H	1
6128.2476	6mm H	1
6128.2477	7mm H	1
6128.2478	8mm H	1
6128.2479	9mm H	1
6128.2480	10mm H	1
6128.2481	11mm H	1
6128.2482	12mm H	1

9128.0026 AGX Modular RP Trial Module, 12x14 and 14x16mm

COALITION AGX® RP 15x18 and 15x20mm MODULAR TRIAL SET 9128.9027

COALITION AGX® RP 15x18, 0° Modular Trials

Part No.	Length	Qty
6128.2505	5mm H	1
6128.2506	6mm H	1
6128.2507	7mm H	1
6128.2508	8mm H	1
6128.2509	9mm H	1
6128.2510	10mm H	1
6128.2511	11mm H	1
6128.2512	12mm H	1

COALITION AGX® RP 15x20, 0° Modular Trials

Part No.	Length	Qty
6128.2525	5mm H	1
6128.2526	6mm H	1
6128.2527	7mm H	1
6128.2528	8mm H	1
6128.2529	9mm H	1
6128.2530	10mm H	1
6128.2531	11mm H	1
6128.2532	12mm H	1

COALITION AGX® RP 15x18, 7° Modular Trials

Part No.	Length	Qty
6128.2575	5mm H	1
6128.2576	6mm H	1
6128.2577	7mm H	1
6128.2578	8mm H	1
6128.2579	9mm H	1
6128.2580	10mm H	1
6128.2581	11mm H	1
6128.2582	12mm H	1

COALITION AGX® RP 15x20, 7° Modular Trials

Part No.	Length	Qty
6128.2535	5mm H	1
6128.2536	6mm H	1
6128.2537	7mm H	1
6128.2538	8mm H	1
6128.2539	9mm H	1
6128.2540	10mm H	1
6128.2541	11mm H	1
6128.2542	12mm H	1

9128.0027 AGX Modular RP Trial Module, 15x18 and 15x20mm

IMPORTANT INFORMATION ON THE COALITION® SPACERS

DESCRIPTION

COALITION® Spacers (including COALITION MIS® and COALITION AGX®) are cervical interbody fusion devices used to provide structural stability in skeletally mature individuals following discectomy. COALITION® Spacers are inserted through an anterior cervical approach, and are available in various heights and geometric options to fit the anatomical needs of a wide variety of patients. Protrusions on the superior and inferior surfaces of each device grip the endplates of the adjacent vertebrae to aid in expulsion resistance. The COALITION AGX® Plate is an anterior cervical fixation device that is available in various lengths and widths to fit the anatomical needs of a wide variety of patients. Screws are inserted through the anterior titanium portion of the implants into adjacent vertebral bodies for bony fixation. The COALITION MIS® Spacer may also be used with anchors inserted through the anterior titanium portion of the implant into adjacent vertebral bodies.

COALITION® Spacers are made from radiolucent PEEK polymer and titanium alloy, with titanium alloy or tantalum markers, as specified in ASTM F2026, F136, F1295, and F560. COALITION® Spacers are additionally available in an all-titanium alloy version. COALITION® PEEK spacers are available with commercially pure titanium plasma spray (TPS) coating, as specified in ASTM F67 and F1580. COALITION AGX® Plates are made from titanium alloy, as specified in ASTM F136, F1295, and F1472. The screws and anchors are manufactured from titanium alloy, as specified in ASTM F136 and F1295.

INDICATIONS

All COALITION® Spacers are interbody fusion devices intended for use in skeletally mature patients with degenerative disc disease (DDD) of the cervical spine (C2-T1) for one or two contiguous levels, depending on the system. DDD is defined as discogenic pain with degeneration of the disc confirmed by history and radiographic studies. These patients should be skeletally mature and have had at least six (6) weeks of non-operative treatment. Hyperlordotic implants ($\geq 20^\circ$) must be used with supplemental fixation in addition to the two screws or anchors. These devices are to be filled with autograft bone and/or allogenic bone graft composed of cancellous, cortical and/or corticocancellous bone. All COALITION® TPS coated spacers are indicated for the same use as non-coated PEEK versions.

The COALITION® Spacer and COALITION AGX® Plate and Spacer assembly are stand-alone integrated interbody fusion devices intended for use at one or two levels of the cervical spine (C2-T1) and used with two titanium alloy screws which accompany the implant.

The COALITION MIS® Spacer is an integrated interbody fusion device intended to be used with two titanium alloy screws and/or anchors which accompany the implants. When used with screws, COALITION MIS® Spacers are stand-alone interbody fusion devices intended for use at one or two levels of the cervical spine (C2-T1). When used with anchors, COALITION MIS® Spacers are intended for use at one level of the cervical spine (C2-T1) with additional supplemental fixation such as posterior cervical screw fixation.

COALITION AGX® Spacer is an interbody fusion device intended to be used with supplemental fixation, such as anterior cervical plates or posterior cervical screw fixation, for one or two levels of the cervical spine (C2-T1). When used with the COALITION AGX® Plate, the plate-spacer assembly takes on the indications for use of the COALITION AGX® Spacer, with the COALITION AGX® Plate acting as the supplemental fixation.

The COALITION AGX® Plate is intended for anterior screw fixation to the cervical spine (C2-C7) for the following indications: degenerative disc disease (as defined by neck pain of discogenic origin with degeneration of the disc confirmed by patient history and radiographic studies), trauma (including fractures), tumors, deformity (defined as kyphosis, lordosis, or scoliosis), pseudoarthrosis, failed previous fusion, spondylolisthesis, and spinal stenosis.

WARNINGS

One of the potential risks identified with this system is death. Other potential risks which may require additional surgery, include:

- device component fracture,
- loss of fixation,
- non-union,
- fracture of the vertebrae,
- neurological injury, and
- vascular or visceral injury.

The COALITION AGX® Plates are not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

Plate contouring is not recommended due to the plate's translational components.

Certain degenerative diseases or underlying physiological conditions such as

diabetes, rheumatoid arthritis, or osteoporosis may alter the healing process, thereby increasing the risk of implant breakage or spinal fracture.

Patients with previous spinal surgery at the level(s) to be treated may have different clinical outcomes compared to those without previous surgery.

Components of this system should not be used with components of any other system or manufacturer.

The components of this system are manufactured from PEEK radiolucent polymer, titanium alloy, commercially pure titanium, and tantalum. Mixing of stainless steel implant components with different materials is not recommended for metallurgical, mechanical and functional reasons.

These warnings do not include all adverse effects that could occur with surgery in general, but are important considerations particular to orthopedic implants. General surgical risks should be explained to the patient prior to surgery.

Use this device as supplied and in accordance with the handling and use information provided below.

PRECAUTIONS

The implantation of cervical fixation and intervertebral fusion devices should be performed only by experienced spinal surgeons with specific training in the use of these systems because this is a technically demanding procedure presenting a risk of serious injury to the patient. Preoperative planning and patient anatomy should be considered when selecting implant size.

Surgical implants must never be reused. An explanted implant must never be reimplanted. Even though the device appears undamaged, it may have small defects and internal stress patterns which could lead to breakage.

Adequately instruct the patient. Mental or physical impairment which compromises or prevents a patient's ability to comply with necessary limitations or precautions may place that patient at a particular risk during postoperative rehabilitation.

For optimal implant performance, the surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact the performance of the system.

MRI SAFETY INFORMATION



The COALITION® Spacers are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla only
- Maximum spatial field gradient of 3,000 gauss/cm (30 T/m) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 1 W/kg

Under the scan conditions defined above, the COALITION® Spacers are expected to produce a maximum temperature rise of less than or equal to 3.9°C after 15 minutes of continuous scanning.

The image artifact caused by the device is not expected to extend beyond 35mm from the device when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

CONTRAINDICATIONS

Use of these implants is contraindicated in patients with the following conditions:

1. Active systemic infection, infection localized to the site of the proposed implantation, or when the patient has a suspected or documented allergy, foreign body sensitivity, or known intolerance to any of the implant materials.
2. Signs of local inflammation.
3. Prior fusion at the level(s) to be treated.
4. Severe osteoporosis, which may prevent adequate fixation.
5. Conditions that may place excessive stresses on bone and implants, such as severe obesity or degenerative diseases, are relative contraindications. The decision whether to use these devices in such conditions must be made by the physician taking into account the risks versus the benefits to the patient.
6. Patients whose activity, mental capacity, mental illness, alcoholism, drug abuse, occupation, or lifestyle may interfere with their ability to follow postoperative restrictions and who may place undue stresses on the implant during bony healing and may be at a higher risk of implant failure.
7. Any patient not willing to cooperate with postoperative instructions.
8. Any condition not described in the indications for use.

IMPORTANT INFORMATION ON THE COALITION® SPACERS

9. Fever or leukocytosis.
10. Pregnancy.
11. Any other condition that would preclude the potential benefit of spinal implant surgery, such as the presence of tumors or congenital abnormalities, fracture local to the operating site, elevation of sedimentation rate unexplained by other diseases, elevations of the white blood count (WBC), or a marked left shift in the WBC differential count.
12. Any case not needing a fusion.
13. Patients with a known hereditary or acquired bone friability or calcification problem should not be considered for this type of surgery.
14. These devices must not be used for pediatric cases or where the patient still has general skeletal growth.
15. Any case where the implant components selected for used would be too large or too small to achieve a successful result.
16. Any case that requires the mixing of metals from two different components or systems.
17. Any patient having inadequate tissue coverage at the operative site or inadequate bone stock or quality.
18. Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance.

COMPLICATIONS AND POSSIBLE ADVERSE EFFECTS

Prior to surgery, patients should be made aware of the following possible adverse effects in addition to the potential need for additional surgery to correct these effects:

- Loosening, bending or breakage of components
- Displacement/migration of device components
- Tissue sensitivity to implant material
- Potential for skin breakdown and/or wound complications
- Non-union or delayed union or mal-union
- Infection
- Nerve damage, including loss of neurological function (sensory and/or motor), paralysis, dysesthesia, hyperesthesia, paresthesia, radiculopathy, reflex deficit, cauda equina syndrome
- Dural tears, cerebral spinal fluid leakage
- Fracture of vertebrae
- Foreign body reaction (allergic) to components or debris
- Vascular or visceral injury
- Change in spinal curvature, loss of correction, height and/or reduction
- Urinary retention or loss of bladder control or other types of disorders of the urogenital system
- Ileus, gastritis, bowel obstruction or other types of gastrointestinal system compromise
- Reproductive system compromise including impotence, sterility, loss of consortium and sexual dysfunction
- Pain or discomfort
- Bursitis
- Decrease in bone density due to stress shielding
- Loss of bone or fracture of bone above or below the level of surgery
- Bone graft donor site pain, fracture, and/or delayed wound healing
- Restriction of activities
- Lack of effective treatment of symptoms for which surgery was intended
- Need for additional surgical intervention
- Death

PACKAGING

These implants and instruments may be supplied pre-packaged and sterile, using gamma irradiation. The integrity of the sterile packaging should be checked to ensure that sterility of the contents is not compromised. Packaging should be carefully checked for completeness and all components should be carefully checked to ensure that there is no damage prior to use. Damaged packages or products should not be used, and should be returned to Globus Medical. During surgery, after the correct size has been determined, remove the products from the packaging using aseptic technique.

The instrument sets are provided nonsterile and are steam sterilized prior to use, as described in the STERILIZATION section below. Following use or exposure to soil, instruments must be cleaned, as described in the CLEANING section below.

HANDLING AND USE

All instruments and implants should be treated with care. Improper use or handling may lead to damage and/or possible malfunction. Products should be checked to ensure that they are in working order prior to surgery. All products should be inspected prior to use to ensure that there is no unacceptable deterioration such as corrosion (i.e. rust, pitting), discoloration, excessive scratches, notches, debris, residue, flaking, wear, cracks, cracked seals, etc. Non-working or damaged instruments should not be used, and should be returned to Globus Medical.

Implants are single use devices and should not be cleaned. Re-cleaning of single use implants might lead to mechanical failure and/or material degradation. Discard any implants that may have been accidentally contaminated.

CLEANING

Instruments should be cleaned separately from instrument trays and cases. Lids should be removed from cases for the cleaning process, if applicable. All instruments that can be disassembled must be disassembled for cleaning. All handles must be detached. Instruments may be reassembled following sterilization. The products should be cleaned using neutral cleaners before sterilization and introduction into a sterile surgical field or (if applicable) return of the product to Globus Medical.

Cleaning and disinfecting can be performed with aldehyde-free solvents at higher temperatures. Cleaning and decontamination must include the use of neutral cleaners followed by a deionized water rinse. Note: certain cleaning solutions such as those containing formalin, glutaraldehyde, bleach and/or other alkaline cleaners may damage some devices, particularly instruments; these solutions should not be used.

The following cleaning methods should be observed when cleaning instruments after use or exposure to soil, and prior to sterilization:

1. Immediately following use, ensure that the instruments are wiped down to remove all visible soil and kept from drying by submerging or covering with a wet towel.
2. Disassemble all instruments that can be disassembled.
3. Rinse the instruments under running tap water to remove all visible soil. Flush the lumens a minimum of 3 times, until the lumens flush clean.
4. Prepare Enzo® (or a similar enzymatic detergent) per manufacturer's recommendations.
5. Immerse the instruments in the detergent and allow them to soak for a minimum of 2 minutes.
6. Use a soft bristled brush to thoroughly clean the instruments. Use a pipe cleaner for any lumens. Pay close attention to hard to reach areas.
7. Using a sterile syringe, draw up the enzymatic detergent solution. Flush any lumens and hard to reach areas until no soil is seen exiting the area.
8. Remove the instruments from the detergent and rinse them in running warm tap water.
9. Prepare Enzo® (or a similar enzymatic detergent) per manufacturer's recommendations in an ultrasonic cleaner.
10. Completely immerse the instruments in the ultrasonic cleaner and ensure detergent is in lumens by flushing the lumens. Sonicate for a minimum of 3 minutes.
11. Remove the instruments from the detergent and rinse them in running deionized water or reverse osmosis water for a minimum of 2 minutes.
12. Dry instruments using a clean soft cloth and filtered pressurized air.
13. Visually inspect each instrument for visible soil. If visible soil is present, then repeat cleaning process starting with Step 3.

CONTACT INFORMATION

Globus Medical may be contacted at 1-866-GLOBUS1 (456-2871). A surgical technique manual may be obtained by contacting Globus Medical.

STERILIZATION

These implants and instruments may be available sterile or nonsterile.

Sterile implants and instruments are sterilized by gamma radiation, validated to ensure a Sterility Assurance Level (SAL) of 10⁻⁶. Sterile products are packaged in a heat-sealed double pouch or container/pouch. The expiration date is provided on the package label. These products are considered sterile unless the packaging has been opened or damaged. Sterile implants and

IMPORTANT INFORMATION ON THE COALITION® SPACERS

instruments that become nonsterile or have expired packaging are considered nonsterile and may be sterilized according to instructions for nonsterile implants and instruments below. Sterile implants meet pyrogen limit specifications.

Nonsterile implants and instruments have been validated to ensure an SAL of 10^{-6} . The use of an FDA-cleared wrap is recommended, per the Association for the Advancement of Medical Instrumentation (AAMI) ST79, *Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities*. It is the end user's responsibility to use only sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization cassettes) that have been cleared by the FDA for the selected sterilization cycle specifications (time and temperature).

When using a rigid sterilization container, the following must be taken into consideration for proper sterilization of Globus devices and loaded graphic cases:










- Recommended sterilization parameters are listed in the table below.
- Only FDA-cleared rigid sterilization containers for use with pre-vacuum steam sterilization may be used.
- When selecting a rigid sterilization container, it must have a minimum filter area of 176 in² total, or a minimum of four (4) 7.5in diameter filters.
- No more than one (1) loaded graphic case or its contents can be placed directly into a rigid sterilization container.
- Stand-alone modules/racks or single devices must be placed, without stacking, in a container basket to ensure optimal ventilation.
- The rigid sterilization container manufacturer's instructions for use are to be followed; if questions arise, contact the manufacturer of the specific container for guidance.
- Refer to AAMI ST79 for additional information concerning the use of rigid sterilization containers.

For implants and instruments provided NONSTERILE, sterilization is recommended (wrapped or containerized) as follows:

Method	Cycle Type	Temperature	Exposure Time	Drying Time
Steam	Pre-vacuum	132°C (270°F)	4 Minutes	30 Minutes
Steam	Pre-vacuum	134°C (273°F)	3 Minutes	30 Minutes

These parameters are validated to sterilize only this device. If other products are added to the sterilizer, the recommended parameters are not valid and new cycle parameters must be established by the user. The sterilizer must be properly installed, maintained, and calibrated. Ongoing testing must be performed to confirm inactivation of all forms of viable microorganisms.

CAUTION: Federal (U.S.A.) Law restricts this Device to Sale by or on the Order of a Physician.

SYMBOL TRANSLATION			
	CATALOGUE NUMBER		STERILIZED BY IRRADIATION
	LOT NUMBER		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
	CAUTION		MANUFACTURER
	SINGLE USE ONLY		USE BY (YYYY-MM-DD)
	QUANTITY		

DI140A REV T

ALLOGRAFT TISSUE SAFETY AND PRODUCT INFORMATION

Tissue Bank

Every tissue bank that processes FORGE® Allograft Spacers is accredited by the American Association of Tissue Banks (AATB), complies with FDA regulations and maintains the highest standards in the recovery, processing, storage and distribution of human allograft tissues. These tissue banks assure unsurpassed safety and the highest quality of FORGE® Allograft Spacers.

Stringent Donor Selection, Testing and Donor Release

Each of Globus' partner tissue banks adheres to strict donor criteria and stringent screening procedures. Every donor undergoes an intense review of medical and social history, along with intensive screening and testing. Donors are tested for HIV 1 and 2, Hepatitis C, Hepatitis B, and Syphilis. In addition, numerous microbiologic cultures are performed and evaluated at tissue recovery and allograft packaging. As a final check, a portion of each donor's tissue is destructively tested for the presence of microorganisms. Donors are determined to be suitable for transplantation only after the tissue bank's Medical Director reviews the medical and social history, relevant hospital records, infectious disease testing, physical exam, and autopsy report (if one was performed).

Controlled Processing and Rigorous Cleaning

AATB certified technicians process FORGE® Allograft Spacers in a highly controlled environment to minimize possible contaminants. FORGE® Allograft Spacers are processed, cleaned and disinfected using proprietary technology which is a combination of treatments of antibiotics, alcohol, peroxide and multiple water rinses to reduce bioburden of the processed tissue. Coupled with a removal of bone marrow, lipids, and proteins, this technology is effective in reducing the risk of viral transmission of processed tissues.

Product-Specific Information

Please refer to the package insert for the specific product for details on storage, handling, and preparation of allograft, as well as tissue handling procedures.



Globus Medical
Valley Forge Business Center
2560 General Armistead Avenue
Audubon, PA 19403
www.globusmedical.com

Customer Service:

Phone 1-866-GLOBUS1 (or 1-866-456-2871)

Fax 1-866-GLOBUS3 (or 1-866-456-2873)

©2022 Globus Medical. All rights reserved. Patent www.globusmedical.com/patents.
Life moves us is a registered trademark of Globus Medical. Please refer to package insert for
description, indications, contraindications, warnings, precautions and other important information.



GMTGD145
01.22 Rev E