



# QUARTEX<sup>®</sup>

## Occipito-Cervico-Thoracic Stabilization System



*Our mission is to deliver cutting-edge technology, research, and innovative solutions to promote healing in patients with musculoskeletal disorders.*

***Life Moves Us***

The Surgical Technique shown is for illustrative purposes only. The technique(s) actually employed in each case always depends on the medical judgment of the surgeon exercised before and during surgery as to the best mode of treatment for each patient. Additionally, as instruments may occasionally be updated, the instruments depicted in this Surgical Technique may not be exactly the same as the instruments currently available. Please consult with your sales representative or contact Globus directly for more information.

# SURGICAL TECHNIQUE GUIDE

## QUARTEX®

Overview .....	4
Implant Overview .....	6
Instrument Overview .....	8
Surgical Technique	
1. Approach.....	20
2. Insertion Site Preparation .....	20
3. Screw Insertion.....	21
Optional: Hook Insertion .....	22
4. Rod and Cap Insertion .....	22
Option A: Rocker Fork Reducer .....	23
Option B: Tower Reducer.....	23
Option C: Grip Reducer .....	24
5. Occipital Fixation.....	24
<i>Using the Occipital Plate Bender .....</i>	<i>25</i>
<i>Using the Occipital Plate Tip Bender.....</i>	<i>25</i>
6. Optional Connectors.....	28
Final Construct .....	30
Optional: Implant Removal.....	31
QUARTEX® Instrument Set I.....	32
QUARTEX® Instrument Set II .....	34
QUARTEX® GII Guided instrument Set .....	36
QUARTEX® Implant Set.....	38
QUARTEX® Occipital Set .....	46
QUARTEX® Large Diameter Screw Set .....	50
3.5mm Curved Rod Set .....	52
4.0mm Curved Rod Set.....	53
QUARTEX® Occipital Clamp Set .....	54
QUARTEX® Auxiliary Screw Module.....	55
Important Information.....	56

# QUARTEX<sup>®</sup>

## Occipito-Cervico-Thoracic Stabilization System

The QUARTEX<sup>®</sup> Occipito-Cervico-Thoracic (OCT) Stabilization System offers a variety of solutions to the challenges associated with posterior OCT fusion to aid in easier construct assembly.

With a comprehensive range of instruments, including threaded drivers and streamlined reduction tools, installation of implants is swift and simple in the midst of complex pathologies. The threaded locking cap design allows for single step locking for speed and efficiency.



## Screw Head Capability

The QUARTEX® screw head accepts a 3.5 or 4.0mm rod and offers 90° of conical angulation.

## Threaded Locking Cap

The low-torque, single step locking provides speed and efficiency.

## Refined Instruments

Stable threaded drivers, efficient reduction options and flexible tools facilitate easy construct assembly.

## True Thoracic Pedicle Screw Offerings

QUARTEX® offers true thoracic screws with diameters up to 5.5mm to accommodate patient anatomy.



# IMPLANT OVERVIEW

## Polyaxial Screw

- 90° conical screw angulation of 3.5, 4.0, and 4.5mm screws provides intraoperative versatility for complex pathologies
- Large diameter 5.0 and 5.5mm screws provide 70° conical screw angulation
- Screw head accepts a 3.5 or 4.0mm rod
- Screw lengths from 8–50mm
- Screw diameters from 3.5–5.5mm
- Shoulder screws available



## Screw Thread

- Small diameter screws feature a tapered and fluted tip with a single lead thread
- Large diameter screws feature a tapered tip and a dual lead thread



## Threaded Locking Cap

- Low-torque single step locking simplifies construct assembly



## Rods

- 3.5 and 4.0mm rods in titanium alloy and cobalt chrome alloy
- Curved rods available in lengths from 25–120mm
- Tapered rod transitions from 3.5 or 4.0mm to 4.0, 4.5, 4.75, 5.0, 5.5, 6.0, 6.35 and 6.5mm diameters
- Pre-bent and hinged rod options for occipital fixation



## Hooks


- Threaded low profile hooks in left, right and inline configurations





## Head-to-Head T-Connector

- Allows insertion when screws are directly adjacent
- Low profile curved design accommodates patient anatomy
- Sizes ranging from 21–63mm



## Rod-to-Rod T-Connector

- Low profile design
- Medial-lateral translation allows versatility in placement
- Clamp snaps over rod for controlled insertion
- Sizes ranging from 21–63mm



## Lateral Connector

- 9, 20, 30mm, and angled options



## Parallel Connectors

- Allow attachment of QUARTEX® to 3.5–6.5mm diameter rod systems



## Occipital Plate

- Rod acceptors rotate and translate to accommodate rod position
- Utilizes the same threaded locking cap as polyaxial screws
- Five points of midline fixation
- Variable angle, self-tapping screws in 4.2 and 4.6mm diameters, in lengths from 6–16mm
- Bend zones for contouring
- Design allows plate to lag to the bone
- Multiple plate sizes available



# INSTRUMENT OVERVIEW

## PREPARATION INSTRUMENTS



Awl, 2.4mm 682.113



Pedicule Probe, Straight 6149.0050



Pedicule Probe, Curved 6149.0055

## Drill Bits



Drill Bit 3.5mm (2.4mm diameter) 6149.1035

(S) Sterile Packed 6149.1036S



Drill Bit 3.5, Short 6149.1037

(S) Sterile Packed 6149.1038S



Drill Bit 4.0 (2.9mm diameter) 6149.1040

(S) Sterile Packed 6149.1041S



Drill Bit 4.0, Short 6149.1042

(S) Sterile Packed 6149.1043S



Drill Bit 4.5 (3.4mm diameter) 6149.1045

(S) Sterile Packed 6149.1046S



## PREPARATION INSTRUMENTS (CONT'D)



Drill Guide with 14mm Stop 6149.1014



Drill Guide with Adjustable Stop,  
3.5-4.5 Diameter, 6-50mm 6149.1000

## Taps



3.0mm Screw	6149.1130
3.5mm Screw	6149.1135
4.0mm Screw	6149.1140
4.5mm Screw	6149.1145



5.0mm Screw	6149.1150
5.5mm Screw	6149.1155

## PREPARATION INSTRUMENTS (CONT'D)



Adjustable Depth Tap Sleeve 6149.1110



Ball Tip Probe 682.115

## SCREW INSTRUMENTS



Rigid Driver 6149.2000



Spin Sleeve 6149.2001



Screwdriver, #15, Self-Retaining, Quick-Connect 6149.2100



Quick-Connect Handle, Swivel 636.450



Pushbutton Driver, #15 Hexalobe 6149.2700



Head Turner 6149.3000

## ROD INSTRUMENTS



Rod Template 240mm 615.301




Tube Style Rod Bender 6149.3090



Rod Bender 615.304



*In Situ* Rod Cutter, 4.0mm 6149.3095



Rod Holding Forceps, 3.5mm Rod 682.308



Rod Holding Forceps, 4.0mm Rod 6149.3055

## ROD INSTRUMENTS (CONT'D)



Tower Reducer 6149.3200




Tower Adapter, A/O Quick-Connect 6149.3210



Rocker Fork Reducer 6149.3250



Compressor, 4.0mm Rod 6149.3500



Distractor, 4.0mm Rod 6149.3505

## ROD INSTRUMENTS (CONT'D)



Heavy Duty Tower Reducer, 10mm 6149.3280

Heavy Duty Tower Reducer, 20mm 6149.3290



Top Loading Flip Ratcheting Reducer 6149.3260



Top Loading Ratcheting Reducer 6149.3220

## LOCKING INSTRUMENTS



Threaded Locking Cap Inserter 6149.4000



Screwdriver, #15, Final Tightening, Quick-Connect 6149.4200

## LOCKING INSTRUMENTS (CONT'D)



Torque-Limiting Quick-Connect Handle, 2.0Nm 6149.4220



Counter-Torque 6149.4100

## OTHER INSTRUMENTS



Hook Forceps 6149.3040



Head-to-Head, Bender, Left 6149.5000



Head-to-Head, Bender, Right 6149.5001

## OTHER INSTRUMENTS (CONT'D)



Head-to-Head T-Connector Nut Inserter 682.312



2.5mm Hex Driver, Quick-Connect 610.831

## OCCIPITAL INSTRUMENTS



Occipital Plate Holder 682.603



Occipital Plate Bender 6149.6800



Occipital Plate Tip Bender 682.662



## OCCIPITAL INSTRUMENTS (CONT'D)



Occipital Adjustable Depth Drill Guide 6149.6300



Flexible Drill Bit, 2.9mm 682.622



Occipital U-Joint Sleeve 6149.6500



Occipital U-Joint Screwdriver, #15, Self-Retaining 6149.6555



Occipital U-Joint Caps for Driver 6149.6505

## OCCIPITAL INSTRUMENTS (CONT'D)

### Occipital U-Joint Cap, Stop



6mm Stop	6149.6506
8mm Stop	6149.6508
10mm Stop	6149.6510
12mm Stop	6149.6512
14mm Stop	6149.6514*
16mm Stop	6149.6516*



Occipital U-Joint Tap, 4.2mm Screws 6149.6540

Occipital U-Joint Tap, 4.6mm Screws 6149.6545



Socket Driver 682.602



Offset Handle, Occipital CV Joint 682.620



Offset Handle, Locking Ring 682.619

\*Additionally Available

## OCCIPITAL INSTRUMENTS (CONT'D)



Occipital Rod Slot Positioner 6149.6010



Occipital Rocker Fork Reducer 6149.6250



Jointed Rod Holder 6149.3070



Occipital Counter-Torque 6149.6600

## OCCIPITAL INSTRUMENTS (CONT'D)



Torque-Limiting 2.5mm Hex Driver, 2.5Nm 615.203



Screw Removal Tool 6149.6900



Angled Screw Removal Tool 6149.6910



Angled Screw Removal Tool Counter-Torque 682.657

# SURGICAL TECHNIQUE

## QUARTEX®

QUARTEX® Navigation Instruments, including drills, taps, pedicle probes and screwdrivers, may be used in conjunction with the Medtronic StealthStation®. Refer to the QUARTEX® Guided Instruments (GI) Technique Guide (GMTGD228) for detailed instructions.

### STEP 1 APPROACH

The patient is placed under general anesthesia and positioned prone. The operative area is carefully cleaned and an incision is made at the appropriate level(s). While there are various techniques for implant insertion, a standard midline approach and the assembly of an occipito-cervico-thoracic (OC-T3) construct is described for the purposes of this surgical technique guide. Lateral C-arm fluoroscopy or other radiographic methods may be utilized throughout surgery to ensure correct implant placement.

### STEP 2 INSERTION SITE PREPARATION

Locate pedicles or lateral masses (C1-T3) and remove bone and/or soft tissue as needed using standard instruments.

Align the **Awl, 2.4mm** in the trajectory for the screw and perforate the cortex to create a pilot hole. The **Pedicle Probe, Straight** or **Pedicle Probe, Curved** may be used to open the pedicle or lateral mass pathway. The **Ball Tip Probe** may be used to verify that the walls of the prepared pathway are not violated.


Determine the appropriate screw diameter and length. Choose the corresponding drill bit for the selected screw diameter and length. The **Drill Guide with Adjustable Stop, 6-50mm** allows drill depths from 6-50mm, in 2mm increments. Rotate the gold ring and pull down the tapered sleeve to release the ratchet. Adjust the drill stop until the appropriate depth is indicated and release the sleeve. Rotate the gold ring until secure to lock the drill guide at the appropriate depth. If the required depth is 14mm, the **Drill Guide, Fixed Stop 14mm** may be used.

Attach the drill bit to the **Quick-Connect Handle, Swivel** and insert the assembly through the drill guide and drill to the stop. The **Depth Gauge** may be used to verify depth.

In dense bone, the hole should be tapped to the full length of the intended screw trajectory with the appropriate diameter tap for ease of insertion. If preparing the pedicles with a **Tap**, the **Adjustable Depth Tap Sleeve** may be used to indicate depth. Insert the tap into the sleeve to the desired depth and rotate the knob to lock in place.



Diameter (mm)				
Instrument	3.0	3.5	4.0	4.5
Drill	N/A	2.4	2.9	3.4
Tap	2.9	3.4	3.9	4.4



Diameter (mm)		
Instrument	5.0	5.5
Large Diameter Tap	4.9	5.4

## STEP 3 SCREW INSERTION

To assemble the **Rigid Driver**, connect the Quick-Connect Handle, Swivel onto the flat portion of the proximal end of the driver. Push the gold button on the knurled knob and retract towards the handle. If desired, the **Spin Sleeve** may be loaded onto the driver by inserting the driver into the cannulated portion of the sleeve and rotating the sleeve until it is fully captured. Snap the Spin Sleeve towards the knurled knob to position for loading the screw.

Prior to inserting, screw size can be verified by checking the length and diameter markings on the screw head, or by using the gauges provided in the screw module.


To load a QUARTEX® polyaxial screw onto the driver, press the gold button and retract the outer sleeve toward the handle and fully insert the hexalobe tip into the screw head. Advance the sleeve forward, toward the screw head, and rotate the knurled knob to thread the driver into the screw head until fully seated.



Inserting screw with Rigid Driver



Rigid Driver fully inserted into screw



Rigid Driver

Insert the screw into the prepared pedicle canal or lateral mass. Ensure that the screw head is free of bone obstruction. If required, a **Reamer** is additionally available to prepare the bone at the insertion site for optimal screw head mobility.

To disengage the driver, rotate the outer knurled knob counterclockwise until it releases the screw head.

Alternately, the **Screwdriver, #15 Self-Retaining** or **Pushbutton Driver** may be used for screw insertion. Once the screws are fully inserted, the screw heads can be oriented to receive the rod using the **Head Turner**.

To reposition a screw *in situ*, use the Screwdriver, #15, Final Tightening.



Polyaxial Screw loaded onto Screwdriver, #15, Self-Retaining



Polyaxial Screw loaded onto Pushbutton Driver

### OPTIONAL: HOOK INSERTION

Hook Forceps can be used for controlled insertion of hooks. The construct is not fully locked until final tightening is completed.




## STEP 4 ROD AND CAP INSERTION

### Rod Preparation

Determine the appropriate length and contour of the rod using the **Rod Template, 240mm**. Rods are available in a variety of lengths, and may be cut using the **In Situ Rod Cutter**. If necessary, bend or contour the rod using the **Rod Bender** by placing the rod into the bender and compressing the handles until the desired curvature is achieved. The **In Situ Benders, 4.0** can be used to bend 3.5 or 4.0mm diameter rods.

### Rod Insertion

Insert the rod into the screw head or hook using the **Rod Holding Forceps**. If the rod cannot be inserted fully, the rod may need to be contoured or rod reduction may be necessary.



Rod insertion

### Locking Cap Insertion

Thread the locking cap into the screw head to capture the rod using the **Threaded Locking Cap Inserter**. If necessary, compress or distract the screws or hooks along the rod using the **Compressor** or **Distractor** prior to final tightening the locking cap.

To final tighten the construct, position the **Counter-Torque** over the screw head and attach the **Screwdriver, #15, Final Tightening, Quick-Connect** to the **Torque-Limiting Quick-Connect Handle, 2.0Nm**. Rotate the screwdriver clockwise until it reaches its torque limit.



Threaded locking cap inserter with Counter-Torque




## Rod Reduction

The QUARTEX® system includes multiple options for rod reduction. Please note that rod reduction instruments are designed to seat the rod into the screw, not to bend the rod. Ensure the rod is properly contoured prior to reduction.

### Option A: Rocker Fork Reducer


Tilt the **Rocker Fork Reducer** upwards so that the rocker arms are orientated vertically and slide the pins of the rocker fork into the chevrons of the screw head. Once engaged, rotate the handle downward to push the rod into the screw head as needed.




Inserting Threaded Locking Cap with Rocker Fork

### Option B: Tower Reducer

To use the **Tower Reducer**, rotate the knob counterclockwise to ensure that it has been fully opened. Line up the rod with the slot on the reducer and click into the screw head. To reduce, rotate the knob clockwise until the rod is reduced by the desired amount. The etching on the reducer aligns when fully reduced. To disengage the instrument from the screw head, rotate the knob counterclockwise until it stops and then slightly rotate the instrument while pulling away from the implant. To aid in rotating the knob of the Tower Reducer, the **Tower Adapter**, **AO Quick-Connect** may be attached to a handle and placed over the rotating knob.



Inserting Threaded Locking Cap with Tower Reducer




Rotate axially to disengage reducer from screw head

## ROD AND CAP INSERTION (CONT'D)

### Option C: Grip Reducer

To prepare the **Top Loading Flip Ratcheting Reducer** for use, release the hinge on the handle and ensure that it is fully open. Line up the rod with the slot on the reducer and click onto the top of the screw head. To reduce, slowly compress the handle until the rod is reduced the desired amount. To disengage the instrument from the screw head, release the hinge on the handle to fully open, and then slightly rotate the instrument while pulling away from the implant.

*Note: For tactile feedback during reduction, open the ratcheting lever while compressing the handle.*



Inserting Threaded Locking Cap with Top Loading Flip Ratcheting Reducer

## STEP

## 5

## OCCIPITAL FIXATION

Occipital plates are secured to the occiput using bone screws. Screws up to 12mm in length may be safely inserted into the occipital bone, close to the external occipital protuberance (EOP). Longer screws up to 16mm may be used if they are angled medially close to the EOP. Care must be taken not to penetrate the far cortex of the occipital bone.

Use the **Occipital Plate Holder** to place the appropriate sized plate on the occiput. If it is necessary to increase the plate contour, the **Occipital Plate Bender** may be used to bend the plate within the indicated bend zone.

To use an occipital clamp, slide the rod into the rod acceptor of the clamp and hold in position by provisionally tightening the set screw with the Screwdriver, #15, Self-Retaining and Quick-Connect Handle, Swivel assembly. Prepare the screw hole(s) as described on page 26 and insert the appropriate length screw through the clamp hole(s). The set screw on the clamp must be final tightened using the Screwdriver, #15, Final Tightening, Quick-Connect with the Torque-Limiting Quick-Connect Handle, 2.0Nm.




## USING THE OCCIPITAL PLATE BENDER

Rotate the handle clockwise so that the metal insert of the plate bender elevates to receive the plate. Insert the plate into the bender so that the two locator pins fit into the two lateral holes on the plate.

For more convexity, install the plate's top surface facing down so the rod acceptors are within the instrument, as shown below. For more concavity, install the plate with the top surface facing up.

To securely hold the plate while bending, compress the handle grip. Rotate the rotating handle counterclockwise until the desired plate contour is achieved. The handle stops when the maximum contour is achieved.



## USING THE OCCIPITAL PLATE TIP BENDER

To bend the tip of the plate, ensure that the plate is inserted with the top surface of the plate facing down so the rod acceptors are within the instrument as shown above. The metal insert of the plate bender should be elevated. Compress the handle grip to securely hold the plate.

Slide the **Occipital Plate Tip Bender** completely over the tip of the plate and compress until the tip bender reaches the handle grip.

The tip bender should only be compressed in line with the handle grip, as shown at right. Avoid twisting or tilting the tip bender while compressing. Follow these steps to ensure that bending does not compromise the structural integrity of the plate.



# OCCIPITAL FIXATION (CONT'D)


## Screw Hole Preparation

The Awl 2.4mm may be used to perforate the cortex of the occiput. The **Occipital Adjustable Depth Drill Guide** allows drill depth from 6 to 16mm, in 2mm increments. To adjust the drill stop, rotate the gold ring and pull down the tapered sleeve to release the ratchet and adjust to desired depth. Release the sleeve and rotate the gold ring until secure to lock the drill guide at the appropriate depth.

The **Flexible Drill Bit, 2.9mm** allows up to 90° insertion angle when needed in the occiput. Insert the drill bit and drill through the Occipital Adjustable Depth Drill Guide.

Choose the **Occipital Joint Tap** of appropriate diameter and assemble with the **Occipital U-Joint Sleeve**. Select the desired depth for the **Occipital U-Joint Cap Stop** and thread onto the U-Joint sleeve. Attach the tap assembly to a Quick-Connect handle and tap to the desired depth. Tapping is recommended in occipital bone due to its density. QUARTEX® 6 and 8mm occipital screws do not have self-tapping flutes; 10 to 20mm occipital screws have self-tapping flutes.

An additional **Offset Handle** may be used to gain leverage during tapping or screw insertion. Assemble the offset handle to the U-Joint sleeve and thread the **Offset Handle, Locking Ring** to secure the offset handle.




## Screw Insertion and Locking

Select the desired occipital screw and insert using the **Occipital U-Joint Screwdriver, #15, Self-Retaining**. A joint sleeve is available to maintain the angle of the screw during insertion.

QUARTEX® occipital screws are designed to automatically lock to the occipital plate during insertion. The screw provides a tactile feel and four divots are visible around the screw when fully seated, as shown below. No additional tightening is necessary.



Screw fully seated



Four divots visible when fully seated

## Locking Cap Insertion

Use the **Occipital Rod Slot Positioner** to orient the rod acceptors to receive the rod. Place the rod in the acceptors. Thread the locking caps into the occipital tulips using the Threaded Locking Cap Inserter or hexalobe Occipital U-Joint Driver. The locking caps may be inserted with or without reduction instruments.

Two reduction tools are provided for use with the Occipital plate. The **Occipital Counter-Torque** and the **Occipital Rocker Fork Reducer** are used in the same manner as the QUARTEX® Rod Counter-Torque and QUARTEX® Rocker Fork Reducer, respectively.

The construct is not completely locked until final tightening. The Torque-Limiting Quick-Connect Handle, 2.0Nm, may be attached to the Screwdriver, #15, Final Tightening, Quick-Connect or Occipital U-Joint Driver.

## Final Tightening Occipital Jointed Rods

The **Jointed Rod Holder** may be used for rod placement. Place the holder around the locking set screw of the jointed rod and insert the rod into position. With the holder still in place, use the **Torque-Limiting 2.5mm Hex Driver, 2.5Nm** to tighten the set screw of the jointed rod until it reaches the torque limit.



Offset Handle on Occipital U-Joint Screwdriver



Occipital Rocker Fork



Using Jointed Rod Holder and Torque-Limiting Driver



Various connectors may be used to connect rods or lengthen a construct, as described below. Stainless steel implants should not be connected to titanium alloy or cobalt chrome alloy implants due to the risk of galvanic corrosion.

### Parallel Connectors

Parallel Connectors can be used to aid in transitioning between screw or hook locations and to minimize rod contouring. They also enable connections to pre-existing titanium rods systems of various diameters. The 3.5 and 4.0mm QUARTEX® rods may be connected to a 3.5, 4.0, 4.5, 4.75, 5.0, 5.5, 6.0, 6.35 or 6.5mm diameter rod system using the corresponding connectors.

Slide the Parallel Connectors onto the rods to be connected. Once the Parallel Connector is attached to both rods, the set screws should be final tightened using the Screwdriver, #15, Final Tightening, Quick-Connect, with the Torque-Limiting Quick-Connect Handle, 2.0Nm.



### Lateral Connectors

Lateral Connectors are used to aid in transitioning between screw or hook locations to minimize rod contouring.

Slide the clamp of the lateral connector onto the rod. Introduce the rod portion of the connector into the opening of the polyaxial screw or hook. Insert the threaded locking cap onto the rod with the Threaded Locking Cap Inserter. The Lateral Connector is now connected to the construct, but is not completely locked until final tightening of both the locking cap and set screw.

*Note: Lateral Connectors may be cut to length using the In Situ Rod Cutter. Ensure that the desired rod length passes completely through the screw or hook head.*



### Rod-to-Rod T-Connectors

The Rod-to-Rod T-Connectors are used to securely link the rods of the construct together.

Select the Rod-to-Rod T-Connector that corresponds with the rod diameter that is being implanted. Gently snap one side of the connector onto the rod. Adjust the sliding assembly to the desired width and snap the other side of the connector onto the contralateral rod. To final tighten the connector, install the **Head-to-Head T-Connector Nut Inserter** onto the 2.0Nm Torque Limiting Handle and final tighten the locking nut on each side of the connector.




## Head-to-Head T-Connectors

Head-to-Head T-Connectors attach to contralateral polyaxial screw heads only and are useful when screw heads are in close proximity and a traditional Rod-to-Rod connector cannot be used. These connectors use a special Head-to-Head Locking Cap and Head-to-Head Locking Nut.


Install the Head-to-Head Threaded Locking Cap using the Threaded Locking Cap Inserter. Connect the Screwdriver, #15, Final Tightening, Quick-Connect to the Torque-Limiting Quick-Connect. Final tighten the Head-to-Head Locking Cap. Choose the appropriate length connector and place it onto the screw heads using a pair of forceps. Connect the **Head-to-Head T-Connector Nut Inserter** to the Torque-Limiting Quick-Connect Handle, 2.0Nm. Install a Head-to-Head Locking Nut onto the threaded portion of each locking cap until it reaches its torque limit (2.0Nm).

The Head-to-Head T-Connectors are pre-curved; however, if more curvature is desired, the **Head-to-Head Benders** (right and left) may be used. Place the lateral ends of the Head-To-Head T-Connector into the mating part of the Right and Left Benders and manually bend the T-Connector to the desired curvature.

*Note: The Head-to-Head Locking Cap may be used with the Rocker Fork Reducer, but not with the Top Loading Ratcheting Reducer or Tower Reducer.*



Head-to-Head T-Connector Benders



Placement of Head-to-Head T-Connector

## QUARTEX® ADDITION®-C Implants

Refer to the QUARTEX® ADDITION®-C Technique Guide (GMTGD265) for instructions.



Headed-Headed



Headed-Headed



Open-Closed



Open-Open



Zero-Run-On-Rod



Swivel



CREO AMP  
QUAD®




# FINAL CONSTRUCT




## OPTIONAL: IMPLANT REMOVAL

For revision or removal of QUARTEX® implants, reverse the insertion steps until the desired implants are removed. Loosen and remove all threaded locking caps using the threaded locking cap driver. Once the threaded locking caps are loosened and removed, grasp the rods and remove from the screws or hooks. Remove all screws using a hexalobe driver. Connectors may remain connected on the rods for removal, or may be removed separately.

To remove occipital screws from an occipital plate or clamp that may have stripped in the bone, use either the **Screw Removal Tool** or the **Angled Screw Removal Tool**. For the straight option, prepare the screw removal tool by rotating the silver knob counter-clockwise, exposing the hexalobe tip. Attach the screw removal tool to the intended screw by aligning the drive feature and then rotating the gold knob clockwise to retain the screw's threaded connection. While holding the silver triangular knob stationary, rotate the blue handle counterclockwise to back the screw out of the plate or clamp.



If a steep approach angle is required to remove occipital screws, use the **Angled Screw Removal Tool** with the **Angled Screw Removal Tool Counter-Torque**. To assemble, insert the Angled Screw Removal Tool into the Angled Screw Removal Tool Counter-Torque and rotate the silver wheel until the wheel is finger-tight against the removal tool. Attach the Quick-Connect Handle, Swivel. Attach the removal tool to the intended screw by aligning the drive feature and then rotating the gold knob clockwise to retain the screw's threaded connection. Depress the button on the counter-torque and continuously hold while rotating the blue handle counterclockwise to back the screw out of the plate.






# QUARTEX®

## INSTRUMENT SET I 9149.9001

Instruments	Qty	Instruments	Qty
1 6149.0050 Pedicle Probe, Straight	1	6149.4020 Threaded Locking Cap Inserter, Dual Sided	
2 6149.0055 Pedicle Probe, Curved	1	20 6149.4100 Counter-Torque	1
3 6149.1000 Drill Guide with Adjustable Stop, 3.5-4.5mm, 6-50mm	1	21 6149.4200 Screwdriver, #15, Final Tightening, Quick-Connect	2
6149.1010 Drill Guide with 10mm Stop		6149.4020 Threaded Locking Cap Inserter, Dual Sided	
6149.1012 Drill Guide with 12mm Stop		22 6149.4220 Torque-Limiting Quick-Connect Handle, 2.0Nm	2
4 6149.1014 Drill Guide with 14mm Stop	1	23 6149.5000 Head-to-Head Bender, Left	1
6149.1016 Drill Guide with 16mm Stop		24 6149.5001 Head-to-Head Bender, Right	1
6149.1018 Drill Guide with 18mm Stop		615.111 Lamina Finder	
6149.1020 Drill Guide with 20mm Stop		25 615.301 Rod Template, 240mm	1
5 6149.1035 Drill Bit for 3.5mm Screw	1	26 636.450 Quick-Connect Handle, Swivel	2
6149.1036S Drill Bit for 3.5mm Screw		27 682.112 Depth Gauge	1
6 6149.1037 Drill Bit for 3.5mm Screw, Short	1	28 682.113 Awl, 2.4mm	1
6149.1038S Drill Bit for 3.5mm Screw, Short		29 682.115 Ball Tip Probe	1
7 6149.1040 Drill Bit for 4.0mm Screw	1	682.123 Awl, 2.4mm, Long	
6149.1041S Drill Bit for 4.0mm Screw		30 682.312 Head-to-Head T-Connector Nut Inserter	2
8 6149.1042 Drill Bit for 4.0mm Screw, Short	1		
6149.1043S Drill Bit for 4.0mm Screw, Short		9149.0001 QUARTEX® Instruments Graphics Case	
9 6149.1045 Drill Bit for 4.5mm Screw	1		
6149.1046S Drill Bit for 4.5mm Screw			
10 6149.1110 Adjustable Depth Tap Sleeve	1		
11 6149.1135 Tap for 3.5mm Screw	1		
12 6149.1140 Tap for 4.0mm Screw	1		
13 6149.1145 Tap for 4.5mm Screw	1		
6149.1500 Reamer			
14 6149.2000 Rigid Driver	2		
15 6149.2001 Spin Sleeve	2		
16 6149.2100 Screwdriver, #15, Self-Retaining, Quick-Connect	2		
6149.2105 Screwdriver, #15, Self-Retaining, Small Handle			
6149.2500 Quick-Connect Handle, Small, Ratcheting			
17 6149.2700 Pushbutton Driver, #15 Hexalobe	2		
18 6149.3000 Head Turner	1		
19 6149.4000 Threaded Locking Cap Inserter	4		



**QUARTEX®**  
INSTRUMENT SET I 9149.9001 (CONT'D)





# QUARTEX®

## INSTRUMENT SET II 9149.9011

	Instruments	Qty
1	610.831 2.5mm Hex Driver, Quick-Connect	1
2	6149.3040 Hook Forceps	1
3	6149.3055 Rod Holding Forceps, 4.0mm Rod	1
	6149.3065 Rod Gripper, 4.0mm	
4	6149.3095 <i>In Situ</i> Rod Cutter, 4.0mm	1
5	6149.3100 <i>In Situ</i> Bender, 4.0mm Rod, Left	1
6	6149.3105 <i>In Situ</i> Bender, 4.0mm Rod, Right	1
	6149.3210 Tower Adapter, A/O Quick-Connect	
7	6149.3200 Tower Reducer	2
	6149.3220 Top Loading Ratcheting Reducer	
	6149.3230 Lateral Reducer	
	6149.3240 Inline Reducer	
8	6149.3250 Rocker Fork Reducer	1
9	6149.3260 Grip Reducer	1
	6149.3280 Heavy Duty Tower Reducer, 10mm	
	6149.3290 Heavy Duty Tower Reducer, 20mm	
10	6149.3500 Compressor, 4.0mm Rod	1
11	6149.3505 Distractor, 4.0mm Rod	1
12	615.304 Rod Bender	1
	682.306 <i>In Situ</i> bender, Left. 3.5mm Rod	
	682.307 <i>In Situ</i> Bender, Right, 3.5mm Rod	
13	682.308 Rod Holding Forceps, 3.5mm Rod	1
	682.310 3.5mm Rod Gripper	
	9149.0011 QUARTEX® Instruments II Graphic Case	



# QUARTEX® INSTRUMENT SET II 9149.9011 (CONT'D)



# QUARTEX® GI1


## GUIDED INSTRUMENT SET 9149.9012

	Instruments	Qty
1	6067.0010      Straight Handle, Ratcheting, 1/4" Quick-Connect	2
	6067.0030      Palm Handle, Ratcheting, 1/4" Quick-Connect	
2	6149.9122      Pedicle Probe, CT, GI1	1
	6149.9123      Pedicle Probe, Large, CT, GI1	
3	6149.9124      Tap for 3.5mm Screw, CT, GI1	1
4	6149.9125      Tap for 4.0mm Screw, CT, GI1	1
5	6149.9126      Tap for 4.5mm Screw, CT, GI1	1
6	6149.9127      Tap for 5.0mm Screw, CT, GI1	1
7	6149.9128      Tap for 5.5mm Screw, CT, GI1	1
8	6149.9135      Drill Bit for 3.5mm Screw, CT, GI1	1
9	6149.9136      Drill Bit for 3.5mm Screw, Short, CT, GI1	1
10	6149.9140      Drill Bit for 4.0mm Screw, CT, GI1	1
11	6149.9141      Drill Bit for 4.0mm Screw, Short, CT, GI1	1
12	6149.9145      Drill Bit for 4.5mm Screw, CT, GI1	1
	6149.9150      Drill Bit for 5.0mm Screw, CT, GI1	
	6149.9155      Drill Bit for 5.5mm Screw, CT, GI1	
13	6149.9200      Rigid Driver, CT, GI1	1
14	6149.9210      Screwdriver, #15, Self-Retaining, CT, GI1	1
15	639.407      1/4" Quick-Connect Adapter	1
	9149.0012      QUARTEX® Guided Instruments, GI1 Graphic Case	1
	999.040      QUARTEX® Guided Instruments Usage Code	



# QUARTEX® GII

## GUIDED INSTRUMENT SET 9149.9012 (CONT'D)



# QUARTEX®

## IMPLANT SET 9149.9002

QUARTEX® Polyaxial Screws			Qty	QUARTEX® Polyaxial Screws (Cont'd)			Qty
1149.3508	3.5mm Polyaxial Screw, 8mm		8	1149.4032	4.0mm Polyaxial Screw, 32mm		
1149.3510	3.5mm Polyaxial Screw, 10mm		16	1149.4034	4.0mm Polyaxial Screw, 34mm		
1149.3512	3.5mm Polyaxial Screw, 12mm		16	1149.4036	4.0mm Polyaxial Screw, 36mm		
1149.3514	3.5mm Polyaxial Screw, 14mm		16	1149.4038	4.0mm Polyaxial Screw, 38mm		
1149.3516	3.5mm Polyaxial Screw, 16mm		16	1149.4040	4.0mm Polyaxial Screw, 40mm		
1149.3518	3.5mm Polyaxial Screw, 18mm		8	1149.4042	4.0mm Polyaxial Screw, 42mm		
1149.3520	3.5mm Polyaxial Screw, 20mm		8	1149.4044	4.0mm Polyaxial Screw, 44mm		
1149.3522	3.5mm Polyaxial Screw, 22mm		8	1149.4046	4.0mm Polyaxial Screw, 46mm		
1149.3524	3.5mm Polyaxial Screw, 24mm		8	1149.4048	4.0mm Polyaxial Screw, 48mm		
1149.3526	3.5mm Polyaxial Screw, 26mm		8	1149.4050	4.0mm Polyaxial Screw, 50mm		
1149.3528	3.5mm Polyaxial Screw, 28mm			1149.4510	4.5mm Polyaxial Screw, 10mm		
1149.3530	3.5mm Polyaxial Screw, 30mm			1149.4512	4.5mm Polyaxial Screw, 12mm		
1149.3532	3.5mm Polyaxial Screw, 32mm			1149.4514	4.5mm Polyaxial Screw, 14mm		
1149.3534	3.5mm Polyaxial Screw, 34mm			1149.4516	4.5mm Polyaxial Screw, 16mm		
1149.3536	3.5mm Polyaxial Screw, 36mm			1149.4518	4.5mm Polyaxial Screw, 18mm		
1149.3538	3.5mm Polyaxial Screw, 38mm			1149.4520	4.5mm Polyaxial Screw, 20mm	4	
1149.3540	3.5mm Polyaxial Screw, 40mm			1149.4522	4.5mm Polyaxial Screw, 22mm	4	
1149.3542	3.5mm Polyaxial Screw, 42mm			1149.4524	4.5mm Polyaxial Screw, 24mm	4	
1149.3544	3.5mm Polyaxial Screw, 44mm			1149.4526	4.5mm Polyaxial Screw, 26mm	4	
1149.3546	3.5mm Polyaxial Screw, 46mm			1149.4528	4.5mm Polyaxial Screw, 28mm	4	
1149.3548	3.5mm Polyaxial Screw, 48mm			1149.4530	4.5mm Polyaxial Screw, 30mm	4	
1149.3550	3.5mm Polyaxial Screw, 50mm			1149.4532	4.5mm Polyaxial Screw, 32mm	4	
1149.4010	4.0mm Polyaxial Screw, 10mm	8		1149.4534	4.5mm Polyaxial Screw, 34mm	4	
1149.4012	4.0mm Polyaxial Screw, 12mm	8		1149.4536	4.5mm Polyaxial Screw, 36mm	4	
1149.4014	4.0mm Polyaxial Screw, 14mm	8		1149.4538	4.5mm Polyaxial Screw, 38mm	4	
1149.4016	4.0mm Polyaxial Screw, 16mm	8		1149.4540	4.5mm Polyaxial Screw, 40mm		
1149.4018	4.0mm Polyaxial Screw, 18mm	8		1149.4542	4.5mm Polyaxial Screw, 42mm		
1149.4020	4.0mm Polyaxial Screw, 20mm	8		1149.4544	4.5mm Polyaxial Screw, 44mm		
1149.4022	4.0mm Polyaxial Screw, 22mm	8		1149.4546	4.5mm Polyaxial Screw, 46mm		
1149.4024	4.0mm Polyaxial Screw, 24mm	4		1149.4548	4.5mm Polyaxial Screw, 48mm		
1149.4026	4.0mm Polyaxial Screw, 26mm	4		1149.4550	4.5mm Polyaxial Screw, 50mm		
1149.4028	4.0mm Polyaxial Screw, 28mm	4					
1149.4030	4.0mm Polyaxial Screw, 30mm	4					


# QUARTEX®

## IMPLANT SET 9149.9002 (CONT'D)

QUARTEX® Polyaxial Shoulder Screws			Qty	QUARTEX® Polyaxial Shoulder Screws			Qty
1149.3722	3.5mm Polyaxial Shoulder Screw, 22mm	4		1149.4228	4.0mm Polyaxial Shoulder Screw, 28mm	4	
1149.3724	3.5mm Polyaxial Shoulder Screw, 24mm	4		1149.4230	4.0mm Polyaxial Shoulder Screw, 30mm	4	
1149.3726	3.5mm Polyaxial Shoulder Screw, 26mm	4		1149.4232	4.0mm Polyaxial Shoulder Screw, 32mm	4	
1149.3728	3.5mm Polyaxial Shoulder Screw, 28mm	4		1149.4234	4.0mm Polyaxial Shoulder Screw, 34mm	4	
1149.3730	3.5mm Polyaxial Shoulder Screw, 30mm	4		1149.4236	4.0mm Polyaxial Shoulder Screw, 36mm	4	
1149.3732	3.5mm Polyaxial Shoulder Screw, 32mm	4		1149.4238	4.0mm Polyaxial Shoulder Screw, 38mm		
1149.3734	3.5mm Polyaxial Shoulder Screw, 34mm	4		1149.4240	4.0mm Polyaxial Shoulder Screw, 40mm		
1149.3736	3.5mm Polyaxial Shoulder Screw, 36mm	4		1149.4242	4.0mm Polyaxial Shoulder Screw, 42mm		
1149.3738	3.5mm Polyaxial Shoulder Screw, 38mm			1149.4244	4.0mm Polyaxial Shoulder Screw, 44mm		
1149.3740	3.5mm Polyaxial Shoulder Screw, 40mm			1149.4246	4.0mm Polyaxial Shoulder Screw, 46mm		
1149.3742	3.5mm Polyaxial Shoulder Screw, 42mm			1149.4248	4.0mm Polyaxial Shoulder Screw, 48mm		
1149.3744	3.5mm Polyaxial Shoulder Screw, 44mm			1149.4250	4.0mm Polyaxial Shoulder Screw, 50mm		
1149.3746	3.5mm Polyaxial Shoulder Screw, 46mm						
1149.3748	3.5mm Polyaxial Shoulder Screw, 48mm						
1149.3750	3.5mm Polyaxial Shoulder Screw, 50mm						
1149.4222	4.0mm Polyaxial Shoulder Screw, 22mm	4					
1149.4224	4.0mm Polyaxial Shoulder Screw, 24mm	4					
1149.4226	4.0mm Polyaxial Shoulder Screw, 26mm	4					

### Locking Cap and Nuts

		Qty
1149.0001	QUARTEX® Threaded Locking Cap	24
1149.0002	QUARTEX® Threaded Locking Cap for Head-to-Head T-Connector	4
1149.0003	QUARTEX® Locking Nut for Head-to-Head T-Connector	4



# QUARTEX<sup>®</sup>


## IMPLANT SET 9149.9002 (CONT'D)

Curved Rods		Qty	Curved Rods		Qty
111.825	3.5mm Curved Rod, 25mm	2	1149.7600	4.0mm Curved Rod, 100mm	
111.830	3.5mm Curved Rod, 30mm		1149.7605	4.0mm Curved Rod, 105mm	
111.835	3.5mm Curved Rod, 35mm		1149.7610	4.0mm Curved Rod, 110mm	
111.840	3.5mm Curved Rod, 40mm		1149.7615	4.0mm Curved Rod, 115mm	
111.845	3.5mm Curved Rod, 45mm		1149.7620	4.0mm Curved Rod, 120mm	
111.850	3.5mm Curved Rod, 50mm		Straight Rods		Qty
111.855	3.5mm Curved Rod, 55mm		1149.7040	4.0mm Rod, Straight, 40mm	4
111.860	3.5mm Curved Rod, 60mm		1149.7080	4.0mm Rod, Straight, 80mm	4
111.865	3.5mm Curved Rod, 65mm		1149.7120	4.0mm Rod, Straight, 120mm	2
111.870	3.5mm Curved Rod, 70mm		1149.7240	4.0mm Rod, Straight, 240mm	2
111.875	3.5mm Curved Rod, 75mm		1149.7250	4.0mm Rod, Straight, 500mm	
111.880	3.5mm Curved Rod, 80mm		182.805	3.5mm Rod, 500mm	
111.885	3.5mm Curved Rod, 85mm		3.5mm Cobalt Chrome Straight Rods		
111.890	3.5mm Curved Rod, 90mm		782.801	3.5mm CoCr Rod, 80mm	
111.895	3.5mm Curved Rod, 95mm		782.802	3.5mm CoCr Rod, 120mm	
111.900	3.5mm Curved Rod, 100mm		782.803	3.5mm CoCr Rod, 240mm	
111.905	3.5mm Curved Rod, 105mm		782.804	3.5mm CoCr Rod, 40mm	
111.910	3.5mm Curved Rod, 110mm	782.805	3.5mm CoCr Rod, 500mm		
111.915	3.5mm Curved Rod, 115mm	2	4.0mm Cobalt Chrome Straight Rods		
111.920	3.5mm Curved Rod, 120mm		7149.7040	4.0mm Straight Rod, 40mm, CoCr	
1149.7525	4.0mm Curved Rod, 25mm		7149.7080	4.0mm Straight Rod, 80mm, CoCr	
1149.7530	4.0mm Curved Rod, 30mm		7149.7120	4.0mm Straight Rod, 120mm, CoCr	
1149.7535	4.0mm Curved Rod, 35mm		7149.7240	4.0mm Straight Rod, 240mm, CoCr	
1149.7540	4.0mm Curved Rod, 40mm		7149.7250	4.0mm Straight Rod, 500mm, CoCr	
1149.7545	4.0mm Curved Rod, 45mm				
1149.7550	4.0mm Curved Rod, 50mm				
1149.7555	4.0mm Curved Rod, 55mm				
1149.7560	4.0mm Curved Rod, 60mm				
1149.7565	4.0mm Curved Rod, 65mm				
1149.7570	4.0mm Curved Rod, 70mm				
1149.7575	4.0mm Curved Rod, 75mm				
1149.7580	4.0mm Curved Rod, 80mm				
1149.7585	4.0mm Curved Rod, 85mm				
1149.7590	4.0mm Curved Rod, 90mm				
1149.7595	4.0mm Curved Rod, 95mm				



# QUARTEX®

## IMPLANT SET 9149.9002 (CONT'D)



# QUARTEX®

## IMPLANT SET 9149.9002 (CONT'D)

### QUARTEX® In-Line Hooks

		Qty
1149.1900	In-Line Hook	2
1149.1901	In-Line Hook, Left	1
1149.1902	In-Line Hook, Right	1

### Titanium Tapered Rods - 350mm

		Qty
182.980	3.5mm (175mm)-4.0mm (175mm) Tapered Rod, Titanium, 350mm Length	
182.981	3.5mm (175mm)-4.5mm (175mm) Tapered Rod, Titanium, 350mm Length	
182.982	3.5mm (175mm)-5.0mm (175mm) Tapered Rod, Titanium, 350mm Length	
182.983	3.5mm (175mm)-5.5mm (175mm) Tapered Rod, Titanium, 350mm Length	2
182.984	3.5mm (175mm)-6.0mm (175mm) Tapered Rod, Titanium, 350mm Length	
182.985	3.5mm (175mm)-6.35mm (175mm) Tapered Rod, Titanium, 350mm Length	
182.986	3.5mm (175mm)-6.5mm (175mm) Tapered Rod, Titanium, 350mm Length	
182.987	Tapered Rod, 3.5mm (175mm) - 4.75mm (175mm), Titanium, 350mm Length	
1149.7460	4.0mm (175mm)-4.5mm (175mm) Tapered Rod, Titanium, 350mm Length	
1149.7461	4.0mm (175mm)-4.75mm (175mm) Tapered Rod, Titanium, 350mm Length	
1149.7462	4.0mm (175mm)-5.0mm (175mm) Tapered Rod Titanium, 350mm Length	
1149.7463	4.0mm (175mm)-5.5mm (175mm) Tapered Rod Titanium, 350mm Length	2
1149.7465	4.0mm (175mm)-6.0mm (175mm) Tapered Rod Titanium, 350mm Length	
1149.7466	4.0mm (175mm)-6.35mm (175mm) Tapered Rod Titanium, 350mm Length	
1149.7467	4.0mm (175mm)-6.5mm (175mm) Tapered Rod Titanium, 350mm Length	

### Cobalt Chrome Tapered Rods - 350mm

782.979	3.5mm (175mm)-3.7mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length
782.980	3.5mm (175mm)-4.0mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length
782.981	3.5mm (175mm)-4.5mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length

### Cobalt Chrome Tapered Rods - 350mm (Cont'd)

782.982	3.5mm (175mm)-5.0mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length
782.983	3.5mm (175mm)-5.5mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length
782.984	3.5mm (175mm)-6.0mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length
782.985	3.5mm (175mm)-6.35mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length
782.986	3.5mm (175mm)-6.5mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length
7149.7460	4.0mm (175mm)-4.5mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length
7149.7461	4.0mm (175mm)-4.75mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length
7149.7462	4.0mm (175mm)-5.0mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length
7149.7463	4.0mm (175mm)-5.5mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length
7149.7465	4.0mm (175mm)-6.0mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length
7149.7466	4.0mm (175mm)-6.35mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length
7149.7467	4.0mm (175mm)-6.5mm (175mm) Tapered Rod, Cobalt Chrome, 350mm Length

### Tapered Rods - 500mm

182.989	3.5mm (250mm)-5.5mm (250mm) Tapered Rod, Titanium, 500mm Length
782.989	3.5mm (250mm)- 5.5mm (250mm) Tapered Rod, Cobalt Chrome, 500mm Length
1149.7464	4.0mm (250mm)-5.5mm (250mm) Tapered Rod Titanium, 500mm Length
7149.7464	4.0mm (250mm)-5.5mm (250mm) Tapered Rod, Cobalt Chrome, 500mm Length

### Tapered Rods - 550mm

182.988	Tapered Rod, 3.5mm (50mm) - 5.5mm (500mm), Titanium, 550mm Length
---------	----------------------------------------------------------------------

### Tapered Rods - 600mm

182.978	Tapered Rod, 3.5mm (100mm)-4.5mm (500mm), Titanium, 600mm Length
---------	---------------------------------------------------------------------

### Tapered Rods - 700mm

182.977	Tapered Rod, 3.5mm (200mm)-5.5mm (500mm), Titanium, 700mm Length
782.977	CoCr Tapered Rod, 3.5mm to 5.5mm, 700mm

# QUARTEX®

## IMPLANT SET 9149.9002 (CONT'D)

### Titanium Tapered Rods - 800mm

182.940	Tapered Rod, 3.5mm (200mm)-4.5mm (600mm), Titanium, 800mm Length
182.941	Tapered Rod, 3.5mm (200mm)-4.75mm (600mm), Titanium, 800mm Length
182.942	Tapered Rod, 3.5mm (200mm)-5.0mm (600mm), Titanium, 800mm Length
182.943	Tapered Rod, 3.5mm (200mm)-5.5mm (600mm), Titanium, 800mm Length
182.944	Tapered Rod, 3.5mm (200mm)-6.0mm (600mm), Titanium, 800mm Length
182.945	Tapered Rod, 3.5mm (200mm)-6.35mm (600mm), Titanium, 800mm Length
182.946	Tapered Rod, 3.5mm (200mm)-6.5mm (600mm), Titanium, 800mm Length
1149.7440	Tapered Rod, 4.0mm (200mm)-4.5mm (600mm), Titanium, 800mm Length
1149.7441	Tapered Rod, 4.0mm (200mm)-4.75mm (600mm), Titanium, 800mm Length
1149.7442	Tapered Rod, 4.0mm (200mm)-5.0mm (600mm), Titanium, 800mm Length
1149.7443	Tapered Rod, 4.0mm (200mm)-5.5mm (600mm), Titanium, 800mm Length
1149.7444	Tapered Rod, 4.0mm (200mm)-6.0mm (600mm), Titanium, 800mm Length
1149.7445	Tapered Rod, 4.0mm (200mm)-6.35mm (600mm), Titanium, 800mm Length
1149.7446	Tapered Rod, 4.0mm (200mm)-6.5mm (600mm), Titanium, 800mm Length

### Cobalt Chrome Tapered Rods - 800mm

782.940	Tapered Rod, 3.5mm (200mm)-4.5mm (600mm), Cobalt Chrome, 800mm Length
782.941	Tapered Rod, 3.5mm (200mm)-4.75mm (600mm), Cobalt Chrome, 800mm Length
782.942	Tapered Rod, 3.5mm (200mm)-5.0mm (600mm), Cobalt Chrome, 800mm Length
782.943	Tapered Rod, 3.5mm (200mm)-5.5mm (600mm), Cobalt Chrome, 800mm Length
782.944	Tapered Rod, 3.5mm (200mm)-6.0mm (600mm), Cobalt Chrome, 800mm Length
782.945	Tapered Rod, 3.5mm (200mm)-6.35mm (600mm), Cobalt Chrome, 800mm Length
782.946	Tapered Rod, 3.5mm (200mm)-6.5mm (600mm), Cobalt Chrome, 800mm Length
7149.7440	Tapered Rod, 4.0mm (200mm)-4.5mm (600mm), Cobalt Chrome, 800mm Length
7149.7441	Tapered Rod, 4.0mm (200mm)-4.75mm (600mm), Cobalt Chrome, 800mm Length
7149.7442	Tapered Rod, 4.0mm (200mm)-5.0mm (600mm), Cobalt Chrome, 800mm Length
7149.7443	Tapered Rod, 4.0mm (200mm)-5.5mm (600mm), Cobalt Chrome, 800mm Length
7149.7444	Tapered Rod, 4.0mm (200mm)-6.0mm (600mm), Cobalt Chrome, 800mm Length
7149.7445	Tapered Rod, 4.0mm (200mm)-6.35mm (600mm), Cobalt Chrome, 800mm Length
7149.7446	Tapered Rod, 4.0mm (200mm)-6.5mm (600mm), Cobalt Chrome, 800mm Length





# QUARTEX®

## IMPLANT SET 9149.9002 (CONT'D)

### QUARTEX® Parallel Connectors Qty

1149.8000	Parallel Connector 4.0mm Rod to 4.0mm Rod	2
1149.8001	Parallel Connector 4.0mm Rod to 4.5mm Rod	
1149.8002	Parallel Connector 4.0mm Rod to 4.75mm Rod	
1149.8003	Parallel Connector 4.0mm Rod to 5.0mm Rod	
1149.8004	Parallel Connector 4.0mm Rod to 5.5mm Rod	2
1149.8005	Parallel Connector 4.0mm Rod to 6.0mm Rod	
1149.8006	Parallel Connector 4.0mm Rod to 6.35mm Rod	
1149.8007	Parallel Connector 4.0mm Rod to 6.5mm Rod	

### QUARTEX® Lateral Connectors Qty

1149.8050	Lateral Connector, 4.0mm Rod	2
1149.8051	Lateral Connector, Long, 4.0mm Rod	2
1149.8052	Lateral Connector, Short, 4.0mm Rod	2
1149.8053	Lateral Connector, Angled, 4.0mm Rod	
1149.8054	Lateral Connector, 3.5mm Rod	2
1149.8055	Lateral Connector, Long, 3.5mm Rod	2
1149.8056	Lateral Connector, Short, 3.5mm Rod	2
1149.8057	Lateral Connector, Angled, 3.5mm Rod	
1149.8058	Lateral Connector, Extra Long, 4.0mm Rod	
1149.8059	Lateral Connector, Extra Long, 3.5mm Rod	
1149.8060	Lateral Connector, Closed, Short, 3.5mm Rod	
1149.8061	Lateral Connector, Closed, Long, 3.5mm Rod	
1149.8062	Lateral Connector, Closed, Short, 4.0mm Rod	
1149.8063	Lateral Connector, Closed, Long, 4.0mm Rod	

### QUARTEX® In-Line Connectors Qty

1149.8070	In-Line Connector, Single, 4.0-4.0mm	2
1149.8071	In-Line Connector, Double, 4.0-4.0mm	
1149.8072	In-Line Connector, Single, 4.0-6.5mm	2
1149.8073	In-Line Connector, Double, 4.0-6.5mm	

### QUARTEX® Rod-to-Rod T-Connectors Qty

1149.8100	3.5mm Rod-to-Rod T-Connector, 21-27mm	2
1149.8101	3.5mm Rod-to-Rod T-Connector, 27-33mm	1
1149.8102	3.5mm Rod-to-Rod T-Connector, 33-39mm	2
1149.8103	3.5mm Rod-to-Rod T-Connector, 39-45mm	1

### QUARTEX® Rod-to-Rod T-Connectors Qty

1149.8104	3.5mm Rod-to-Rod T-Connector, 45-51mm	1
1149.8105	3.5mm Rod-to-Rod T-Connector, 51-57mm	
1149.8106	3.5mm Rod-to-Rod T-Connector, 57-63mm	
1149.8200	4.0mm Rod-to-Rod T-Connector, 21-27mm	2
1149.8201	4.0mm Rod-to-Rod T-Connector, 27-33mm	1
1149.8202	4.0mm Rod-to-Rod T-Connector, 33-39mm	2
1149.8203	4.0mm Rod-to-Rod T-Connector, 39-45mm	1
1149.8204	4.0mm Rod-to-Rod T-Connector, 45-51mm	1
1149.8205	4.0mm Rod-to-Rod T-Connector, 51-57mm	
1149.8206	4.0mm Rod-to-Rod T-Connector, 57-63mm	


# QUARTEX<sup>®</sup>

## IMPLANT SET 9149.9002 (CONT'D)

QUARTEX <sup>®</sup> Head-to-Head T-Connectors		Qty
1149.8300	Head-to-Head T-Connector, 21-27mm	1
1149.8301	Head-to-Head T-Connector, 27-33mm	1
1149.8302	Head-to-Head T-Connector, 33-39mm	1
1149.8303	Head-to-Head T-Connector, 39-45mm	1
1149.8304	Head-to-Head T-Connector, 45-51mm	1
1149.8305	Head-to-Head T-Connector, 51-57mm	
1149.8306	Head-to-Head T-Connector, 57-63mm	

### Graphic Cases

9149.0002	QUARTEX <sup>®</sup> Implants Graphics Case
9149.1200	QUARTEX <sup>®</sup> Screw Module
9149.1201	QUARTEX <sup>®</sup> Locking Cap Module
9149.1202	QUARTEX <sup>®</sup> Rod-to-Rod T-Connector Module
9149.1203	QUARTEX <sup>®</sup> Head-to-Head T-Connector Module
9149.1204	QUARTEX <sup>®</sup> Hooks and Connectors Module
9149.1205	3.5mm Curved Rod Module
9149.1206	4.0mm Curved Rod Module



# QUARTEX®




## OCCIPITAL SET 9149.9003

Occipital Instruments				QUARTEX® Occipital Screws		
			Qty			Qty
1	6149.3070	Jointed Rod Holder	1	1149.9006	4.2mm Occipital Screw, 6mm	8
2	6149.6010	Occipital Rod Slot Positioner	1	1149.9008	4.2mm Occipital Screw, 8mm	8
3	6149.6250	Occipital Rocker Fork Reducer	1	1149.9010	4.2mm Occipital Screw, 10mm	8
4	6149.6300	Occipital Adjustable Depth Drill Guide	1	1149.9012	4.2mm Occipital Screw, 12mm	8
5	6149.6500	Occipital U-Joint Sleeve	2	1149.9014	4.2mm Occipital Screw, 14mm	8
6	6149.6505	Occipital U-Joint Cap for Driver	1	1149.9016	4.2mm Occipital Screw, 16mm	8
7	6149.6506	Occipital U-Joint Cap for 6mm Stop	1	1149.9106	4.6mm Occipital Screw, 6mm	4
8	6149.6508	Occipital U-Joint Cap for 8mm Stop	1	1149.9108	4.6mm Occipital Screw, 8mm	4
9	6149.6510	Occipital U-Joint Cap for 10mm Stop	1	1149.9110	4.6mm Occipital Screw, 10mm	4
10	6149.6512	Occipital U-Joint Cap for 12mm Stop	1	1149.9112	4.6mm Occipital Screw, 12mm	4
	6149.6514	Occipital U-Joint Cap for 14mm Stop		1149.9114	4.6mm Occipital Screw, 14mm	4
	6149.6516	Occipital U-Joint Cap for 16mm Stop		1149.9116	4.6mm Occipital Screw, 16mm	4
11	6149.6540	Occipital U-Joint Tap, 4.2mm Screws	1	QUARTEX® Occipital Plates		
12	6149.6545	Occipital U-Joint Tap, 4.6mm Screws	1			
13	6149.6555	Occipital U-Joint Screwdriver, #15, Self-Retaining	1	1149.9200	Occipital Plate, Small	2
14	6149.6600	Occipital Counter-Torque	1	1149.9210	Occipital Plate, Large	2
15	6149.6800	Occipital Plate Bender	1			
16	6149.6900	Screw Removal Tool	1			
17	6149.6910	Angled Screw Removal Tool	1			
18	615.203	Torque-Limiting 2.5mm Hex Driver, 2.5Nm	1			
	636.450	Quick-Connect Handle, Swivel	1			
19	682.602	Socket Driver	1			
20	682.603	Occipital Plate Holder	1			
21	682.619	Offset Handle Locking Ring	2			
22	682.620	Offset Handle, Occipital CV Joint	2			
23	682.622	Flexible Drill Bit, 2.9mm	1			
24	682.657	Angled Screw Remover Tool				
		Counter-Torque	1			
25	682.662	Occipital Plate Tip Bender	1			



# QUARTEX®

## OCCIPITAL SET 9149.9003 (CONT'D)



# QUARTEX®

## OCCIPITAL SET 9149.9003 (CONT'D)

### Occipital Screws

		Qty
1149.9300	QUARTEX® 4.0mm Occipital Rod 240mm, 100	2
1149.9315	QUARTEX® 4.0mm Occipital Rod 240mm, 115	2
1149.9330	QUARTEX® 4.0mm Occipital Rod 240mm, 130	2
1149.9450	QUARTEX® 4.0mm Occipital Jointed Rod, Right	2
1149.9455	QUARTEX® 4.0mm Occipital Jointed Rod, Left	2
182.953	3.5mm ELLIPSE® Occipital Rod, 240mm, 100°	2
182.954	3.5mm ELLIPSE® Occipital Rod, 240mm, 115°	2
182.955	3.5mm ELLIPSE® Occipital Rod, 240mm, 130°	2
182.956	Occipital Jointed Rod, Right	2
182.957	Occipital Jointed Rod, Left	2
7149.9300	4.0mm Occipital Rod, 240mm, 100, CoCr	
7149.9315	4.0mm Occipital Rod, 240mm, 115, CoCr	
7149.9330	4.0mm Occipital Rod, 240mm, 130, CoCr	
782.953	3.5mm CoCr Occipital Rod 240mm, 100°	
782.954	3.5mm CoCr Occipital Rod 240mm, 115°	
782.955	3.5mm CoCr Occipital Rod 240mm, 130°	

### Occipital Clamps

1149.9511	QUARTEX® 3.5mm Occipital Clamp, Single, Short
1149.9512	QUARTEX® 3.5mm Occipital Clamp, Single, Long
1149.9521	QUARTEX® 3.5mm Occipital Clamp, Double, Short
1149.9522	QUARTEX® 3.5mm Occipital Clamp, Double, Long
1149.9531	QUARTEX® 3.5mm Occipital Clamp, Triple, Short
1149.9532	QUARTEX® 3.5mm Occipital Clamp, Triple, Long
1149.9611	QUARTEX® 4.0mm Occipital Clamp, Single, Short
1149.9612	QUARTEX® 4.0mm Occipital Clamp, Single, Long
1149.9621	QUARTEX® 4.0mm Occipital Clamp, Double, Short
1149.9622	QUARTEX® 4.0mm Occipital Clamp, Double, Long
1149.9631	QUARTEX® 4.0mm Occipital Clamp, Triple, Short
1149.9632	QUARTEX® 4.0mm Occipital Clamp, Triple, Long


### Graphic Cases

9149.0003	QUARTEX® Occipital System Graphics Case
9149.1300	QUARTEX® Occipital Screw Module
9149.1301	QUARTEX® Occipital Plate Module
9149.1302	QUARTEX® Occipital Clamp Module
9149.2300	QUARTEX® U-Joint Accessory Module



# QUARTEX®

## OCCIPITAL SET 9149.9003 (CONT'D)




# QUARTEX® LARGE DIAMETER SCREW SET 9149.9004

QUARTEX® Large Diameter Screws		Qty	Instruments		Qty
1149.5020	5.0mm Polyaxial Screw, 20mm	6	1	6067.0001 Pedicle Awl	1
1149.5022	5.0mm Polyaxial Screw, 22mm	6	2	6067.0010 Straight Handle, Ratcheting, 1/4" Quick-Connect	1
1149.5024	5.0mm Polyaxial Screw, 24mm	6	3	6149.1150 Tap for 5.0mm Screw	1
1149.5026	5.0mm Polyaxial Screw, 26mm	6	4	6149.1155 Tap for 5.5mm Screw	1
1149.5028	5.0mm Polyaxial Screw, 28mm	6	5	624.113 Pedicle Probe, Thoracic, Curved	1
1149.5030	5.0mm Polyaxial Screw, 30mm	6	6	624.114 Pedicle Probe, Thoracic, Straight	1
1149.5032	5.0mm Polyaxial Screw, 32mm	6			
1149.5034	5.0mm Polyaxial Screw, 34mm	6			
1149.5036	5.0mm Polyaxial Screw, 36mm	6		9149.0004 QUARTEX® Large Diameter Graphics Case	
1149.5038	5.0mm Polyaxial Screw, 38mm	6			
1149.5040	5.0mm Polyaxial Screw, 40mm	6		9149.1400 QUARTEX® Large Diameter Screw Module	
1149.5042	5.0mm Polyaxial Screw, 42mm				
1149.5044	5.0mm Polyaxial Screw, 44mm				
1149.5046	5.0mm Polyaxial Screw, 46mm				
1149.5048	5.0mm Polyaxial Screw, 48mm				
1149.5050	5.0mm Polyaxial Screw, 50mm				
1149.5520	5.5mm Polyaxial Screw, 20mm	6			
1149.5522	5.5mm Polyaxial Screw, 22mm	6			
1149.5524	5.5mm Polyaxial Screw, 24mm	6			
1149.5526	5.5mm Polyaxial Screw, 26mm	6			
1149.5528	5.5mm Polyaxial Screw, 28mm	6			
1149.5530	5.5mm Polyaxial Screw, 30mm	6			
1149.5532	5.5mm Polyaxial Screw, 32mm	6			
1149.5534	5.5mm Polyaxial Screw, 34mm	6			
1149.5536	5.5mm Polyaxial Screw, 36mm	6			
1149.5538	5.5mm Polyaxial Screw, 38mm	6			
1149.5540	5.5mm Polyaxial Screw, 40mm	6			
1149.5542	5.5mm Polyaxial Screw, 42mm				
1149.5544	5.5mm Polyaxial Screw, 44mm				
1149.5546	5.5mm Polyaxial Screw, 46mm				
1149.5548	5.5mm Polyaxial Screw, 48mm				
1149.5550	5.5mm Polyaxial Screw, 50mm				



# QUARTEX® LARGE DIAMETER SCREW SET 9149.9004 (CONT'D)



# 3.5mm CURVED ROD SET 9149.9005

Part No.	Description	Qty
111.825	3.5mm Curved Rod, 25mm	2
111.830	3.5mm Curved Rod, 30mm	2
111.835	3.5mm Curved Rod, 35mm	2
111.840	3.5mm Curved Rod, 40mm	2
111.845	3.5mm Curved Rod, 45mm	2
111.850	3.5mm Curved Rod, 50mm	2
111.855	3.5mm Curved Rod, 55mm	2
111.860	3.5mm Curved Rod, 60mm	2
111.865	3.5mm Curved Rod, 65mm	2
111.870	3.5mm Curved Rod, 70mm	2
111.875	3.5mm Curved Rod, 75mm	2
111.880	3.5mm Curved Rod, 80mm	2
111.885	3.5mm Curved Rod, 85mm	2
111.890	3.5mm Curved Rod, 90mm	2
111.895	3.5mm Curved Rod, 95mm	2
111.900	3.5mm Curved Rod, 100mm	2
111.905	3.5mm Curved Rod, 105mm	2
111.910	3.5mm Curved Rod, 110mm	2
111.915	3.5mm Curved Rod, 115mm	2
111.920	3.5mm Curved Rod, 120mm	2
9149.1205	QUARTEX® 3.5mm Curved Rod Module	

# 4.0mm CURVED ROD SET 9149.9006

Part No.	Description	Qty
1149.7525	4.0mm Curved Rod, 25mm	2
1149.7530	4.0mm Curved Rod, 30mm	2
1149.7535	4.0mm Curved Rod, 35mm	2
1149.7540	4.0mm Curved Rod, 40mm	2
1149.7545	4.0mm Curved Rod, 45mm	2
1149.7550	4.0mm Curved Rod, 50mm	2
1149.7555	4.0mm Curved Rod, 55mm	2
1149.7560	4.0mm Curved Rod, 60mm	2
1149.7565	4.0mm Curved Rod, 65mm	2
1149.7570	4.0mm Curved Rod, 70mm	2
1149.7575	4.0mm Curved Rod, 75mm	2
1149.7580	4.0mm Curved Rod, 80mm	2
1149.7585	4.0mm Curved Rod, 85mm	2
1149.7590	4.0mm Curved Rod, 90mm	2
1149.7595	4.0mm Curved Rod, 95mm	2
1149.7600	4.0mm Curved Rod, 100mm	2
1149.7605	4.0mm Curved Rod, 105mm	2
1149.7610	4.0mm Curved Rod, 110mm	2
1149.7615	4.0mm Curved Rod, 115mm	2
1149.7620	4.0mm Curved Rod, 120mm	2
9149.1206	QUARTEX® 4.0mm Curved Rod Module	

# QUARTEX®

## OCCIPITAL CLAMP SET 9149.9007

Part No.	Description	Qty
1149.9511	QUARTEX® 3.5mm Occipital Clamp, Single, Short	4
1149.9512	QUARTEX® 3.5mm Occipital Clamp, Single, Long	4
1149.9521	QUARTEX® 3.5mm Occipital Clamp, Double, Short	4
1149.9522	QUARTEX® 3.5mm Occipital Clamp, Double, Long	4
1149.9531	QUARTEX® 3.5mm Occipital Clamp, Triple, Short	4
1149.9532	QUARTEX® 3.5mm Occipital Clamp, Triple, Long	4
1149.9611	QUARTEX® 4.0mm Occipital Clamp, Single, Short	4
1149.9612	QUARTEX® 4.0mm Occipital Clamp, Single, Long	4
1149.9621	QUARTEX® 4.0mm Occipital Clamp, Double, Short	4
1149.9622	QUARTEX® 4.0mm Occipital Clamp, Double, Long	4
1149.9631	QUARTEX® 4.0mm Occipital Clamp, Triple, Short	4
1149.9632	QUARTEX® 4.0mm Occipital Clamp, Triple, Long	4
9149.1302	QUARTEX® Occipital Clamp Module	

# QUARTEX®

## AUXILIARY SCREW MODULE 9149.9008

### QUARTEX® 3.5mm Polyaxial Screw

		Qty
1149.3528	28mm	4
1149.3530	30mm	4
1149.3532	32mm	4
1149.3534	34mm	4
1149.3536	36mm	4
1149.3538	38mm	4
1149.3540	40mm	4
1149.3542	42mm	2
1149.3544	44mm	2
1149.3546	46mm	2
1149.3548	48mm	2
1149.3550	50mm	2

### QUARTEX® 3.5mm Polyaxial Shoulder Screw Qty

1149.3738	38mm	4
1149.3740	40mm	4
1149.3742	42mm	2
1149.3744	44mm	2
1149.3746	46mm	2
1149.3748	48mm	2
1149.3750	50mm	2

### QUARTEX® 4.0mm Polyaxial Screw

		Qty
1149.4032	32mm	4
1149.4034	34mm	4
1149.4036	36mm	4
1149.4038	38mm	4
1149.4040	40mm	4
1149.4042	42mm	2
1149.4044	44mm	2
1149.4046	46mm	2
1149.4048	48mm	2
1149.4050	50mm	2

### QUARTEX® 4.0mm Polyaxial Shoulder Screw Qty

1149.4238	38mm	4
1149.4240	40mm	4
1149.4242	42mm	2
1149.4244	44mm	2
1149.4246	46mm	2
1149.4248	48mm	2
1149.4250	50mm	2

### QUARTEX® 4.5mm Polyaxial Screw Qty

1149.4510	10mm	4
1149.4512	12mm	4
1149.4514	14mm	4
1149.4516	16mm	4
1149.4518	18mm	4
1149.4540	40mm	2
1149.4542	42mm	2
1149.4544	44mm	2
1149.4546	46mm	2
1149.4548	48mm	2
1149.4550	50mm	2

9149.1500 QUARTEX® Auxiliary Screw Module

# IMPORTANT INFORMATION ON QUARTEX® OCCIPITO-CERVICO-THORACIC SPINAL SYSTEM

## DESCRIPTION

The QUARTEX® Occipito-Cervico-Thoracic Spinal System consists of 3.5mm-4.0mm jointed, straight and pre-bent rods, tapered rods, polyaxial screws, hooks, locking caps, t-connectors, lateral connectors, parallel connectors, in-line connectors, rod-to-rod connectors, rod extension clamps and occipital plates. The implants are composed of titanium alloy (per ASTM F136, F1472, or F1295), stainless steel (per ASTM F138) or cobalt chromium molybdenum alloy (CoCr) (per ASTM F1537).

QUARTEX® constructs may be connected to stabilization systems including ELLIPSE®, PROTEX® CT, PROTEX®, CREO®, REVERE®, or BEACON® Systems using corresponding connectors.

The QUARTEX® system include manual surgical instruments manufactured from stainless steel, as specified in ASTM F899. Navigation Instruments are nonsterile, reusable instruments that can be operated manually or under power using a power drill such as POWERASE™ and are intended to be used with the Medtronic StealthStation® System.

## INDICATIONS

The QUARTEX® Occipito-Cervico-Thoracic Spinal System implants are intended to provide immobilization and stabilization of spinal segments as an adjunct to fusion for the following acute and chronic instabilities of the craniocervical junction, the cervical spine (C1-C7) and the thoracic spine (T1-T3): traumatic spinal fractures and/or traumatic dislocations; instability or deformity; failed previous fusions (e.g. pseudoarthrosis); tumors involving the cervical/thoracic spine; and degenerative disease, including intractable radiculopathy and/or myelopathy, neck and/or arm pain of discogenic origin as confirmed by radiographic studies, and degenerative disease of the facets with instability. These implants are also intended to restore the integrity of the spinal column even in the absence of fusion for a limited time period in patients with advanced stage tumors involving the cervical spine in whom life expectancy is of insufficient duration to permit achievement of fusion. In order to achieve additional levels of fixation, rods may be connected to occipital cervical thoracic or thoracolumbar stabilization systems ranging in diameter from 3.2mm to 6.5mm, using corresponding connectors.

Globus Navigation Instruments are intended to be used during the preparation and placement of QUARTEX® screws during spinal surgery to assist the surgeon in precisely locating anatomical structures in either open or minimally invasive procedures. These instruments are designed for use with the Medtronic StealthStation® System, which is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure, such as a skull, a long bone, or vertebra, can be identified relative to a CT or MR based model, fluoroscopy images, or digitized landmarks of the anatomy.

## WARNINGS

The safety and effectiveness of pedicle screw spinal systems have been established only for spinal conditions with significant mechanical instability or deformity requiring fusion with instrumentation. These conditions are significant mechanical instability or deformity of the thoracic spine secondary to degenerative spondylolisthesis with objective evidence of neurological impairment, fracture, dislocation, spinal tumor, and failed previous fusion (pseudoarthrosis). The safety and effectiveness of these devices for any other conditions are unknown.

Possible adverse effects which may occur and may require additional surgery include: failed fusion or pseudoarthrosis leading to implant breakage; allergic reaction to implant materials; device fracture or failure; device migration or loosening; loss of fixation; vertebral fracture; decrease in bone density; pain, discomfort, or abnormal sensations due to the presence of the device; injury to nerves, vessels, and organs; venous thrombosis, lung embolism and cardiac arrest; and death.

The components of this system are manufactured from titanium alloy, stainless steel or cobalt chromium alloy. Dissimilar metals in contact with each other can accelerate the corrosion process due to galvanic corrosion effects. Mixing of implant components of titanium or cobalt chromium with stainless steel is not recommended, for metallurgical, mechanical and functional reasons.

Certain degenerative diseases or underlying physiological conditions such as diabetes, rheumatoid arthritis, or osteoporosis may alter the healing process, thereby increasing the risk of implant breakage or spinal fracture.

These warnings do not include all adverse effects which could occur with surgery in general, but are important considerations particular to orthopedic implants. General surgical risks should be explained to the patient prior to surgery.

Use this device as supplied and in accordance with the handling and use information provided below.

## PRECAUTIONS

The implantation of posterior screw spinal systems should be performed only by experienced spinal surgeons with specific training in the use of this posterior screw spinal system because this is a technically demanding procedure presenting a risk of serious injury to the patient.

Preoperative planning prior to implantation of posterior cervical screw systems should include review of cross-sectional imaging studies (e.g., CT and/or MRI) to evaluate the patient's cervical anatomy including the transverse foramen, neurologic structures, and the course of the vertebral arteries. If any findings would compromise the placement of these screws, other surgical methods should be considered. In addition, use of intraoperative imaging should be considered to guide and/or verify device placement, as necessary.

Use of posterior cervical pedicle screw fixation at the C3 through C6 spinal levels requires careful consideration and planning beyond that required for lateral mass screws placed at these spinal levels, given the proximity of the vertebral arteries and neurologic structures in relation to the cervical pedicles at these levels.

The implants are for single use only. Surgical implants must never be reused. An explanted metal implant must never be reimplanted. Even though the device appears undamaged, it may have small defects and internal stress patterns which could lead to breakage.

Correct handling of the implant is extremely important. Contouring of metal implants should be avoided where possible. If contouring is necessary, or allowed by design, the surgeon should avoid sharp bends, reverse bends, or bending the device at a screw hole. The operating surgeon should avoid any notching or scratching of the device when contouring it. These factors may produce internal stresses which may become the focal point for eventual breakage of the implant.

Metallic implants can loosen, fracture, corrode, migrate, cause pain, or stress shield bone even after a fracture has healed, particularly in young, active patients. While the surgeon must have the final decision on implant removal, we recommend that whenever possible and practical for the individual patient, fixation devices should be removed once their service as an aid to healing is accomplished. Implant removal should be followed by adequate postoperative management.

Adequately instruct the patient. Mental or physical impairment which compromises or prevents a patient's ability to comply with necessary limitations or precautions may place that patient at a particular risk during postoperative rehabilitation.

For optimal implant performance, the surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact the performance of the system.

## MRI SAFETY INFORMATION

The QUARTEX® Occipito-Cervico-Thoracic Spinal System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of QUARTEX® Occipito-Cervico-Thoracic Spinal System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

## CONTRAINDICATIONS

Certain degenerative diseases or underlying physiological conditions such as diabetes or rheumatoid arthritis may alter the healing process, thereby increasing the risk of implant breakage.

Mental or physical impairment which compromises a patient's ability to comply with necessary limitations or precautions may place that patient at a particular risk during postoperative rehabilitation.

Factors such as the patient's weight, activity level, and adherence to weight bearing or load bearing instructions have an effect on the stresses to which the implant is subjected.

## PACKAGING

These implants and instruments may be supplied pre-packaged and sterile, using gamma irradiation. The integrity of the sterile packaging should be checked to ensure that sterility of the contents is not compromised. Packaging should be carefully checked for completeness and all components should be carefully checked to ensure that there is no damage prior to use. Damaged packages or products should not be used, and should be returned to Globus Medical. During surgery, after the correct size has been determined, remove the products from the packaging using aseptic technique.

The instrument sets are provided nonsterile and are steam sterilized prior to use, as described in the STERILIZATION section below. Following use



# IMPORTANT INFORMATION ON QUARTEX® OCCIPITO-CERVICO-THORACIC SPINAL SYSTEM

or exposure to soil, instruments and instrument trays and cases must be cleaned, as described in the CLEANING section below.

## HANDLING

All implants, instruments, and instrument trays and cases should be treated with care. Improper use or handling may lead to damage and/or possible malfunction. Products should be checked to ensure that they are in working order prior to surgery. All products should be inspected prior to use to ensure that there is no unacceptable deterioration such as corrosion (i.e. rust, pitting), discoloration, excessive scratches, notches, debris, residue, flaking, wear, cracks, cracked seals, etc. Non-working or damaged instruments should not be used, and should be returned to Globus Medical.

Implants are single use devices and should not be cleaned. Re-cleaning of single use implants might lead to mechanical failure and/or material degradation. Discard any implants that may have been accidentally contaminated.

## CLEANING

Instruments should be cleaned separately from instrument trays and cases. Lids should be removed from cases for the cleaning process, if applicable. All instruments that can be disassembled must be disassembled for cleaning. All handles must be detached. Instruments may be reassembled following sterilization. The products should be cleaned using neutral cleaners before sterilization and introduction into a sterile surgical field or (if applicable) return of the product to Globus Medical.

Cleaning and disinfecting can be performed with aldehyde-free solvents at higher temperatures. Cleaning and decontamination must include the use of neutral cleaners followed by a deionized water rinse. Note: certain cleaning solutions such as those containing formalin, glutaraldehyde, bleach and/or other alkaline cleaners may damage some devices, particularly instruments; these solutions should not be used.

The following cleaning methods should be observed when cleaning instruments and instrument trays and cases after use or exposure to soil, and prior to sterilization:

1. Immediately following use, ensure that the instruments are wiped down to remove all visible soil and kept from drying by submerging or covering with a wet towel.
2. Disassemble all instruments that can be disassembled.
3. Rinse the instruments under running tap water to remove all visible soil. Flush the lumens a minimum of 3 times, until the lumens flush clean.
4. Prepare Enzo® (or a similar enzymatic detergent) per manufacturer's recommendations.
5. Immerse the instruments in the detergent and allow them to soak for a minimum of 2 minutes.
6. Use a soft bristled brush to thoroughly clean the instruments. Use a pipe cleaner for any lumens. Pay close attention to hard to reach areas.
7. Using a sterile syringe, draw up the enzymatic detergent solution. Flush any lumens and hard to reach areas until no soil is seen exiting the area.
8. Remove the instruments from the detergent and rinse them in running warm tap water.
9. Prepare Enzo® (or a similar enzymatic detergent) per manufacturer's recommendations in an ultrasonic cleaner.
10. Completely immerse the instruments in the ultrasonic cleaner and ensure detergent is in lumens by flushing the lumens. Sonicate for a minimum of 3 minutes.
11. Remove the instruments from the detergent and rinse them in running deionized water or reverse osmosis water for a minimum of 2 minutes.
12. Dry instruments using a clean soft cloth and filtered pressurized air.
13. Visually inspect each instrument for visible soil. If visible soil is present, then repeat cleaning process starting with Step 3.

## CONTACT INFORMATION

Globus Medical may be contacted at 1-866-GLOBUS1 (456-2871). A surgical technique manual may be obtained by contacting Globus Medical.

## STERILIZATION

These implants and instruments may be available sterile or nonsterile.

Sterile implants and instruments are sterilized by gamma radiation, validated to ensure a Sterility Assurance Level (SAL) of  $10^{-6}$ . Sterile products are packaged in a heat sealed, double foil pouch or container/pouch. The expiration date is provided in the package label. These products are

considered sterile unless the packaging has been opened or damaged.

Nonsterile implants and instruments have been validated to ensure an SAL of  $10^{-6}$ . The use of an FDA-cleared wrap is recommended, per the Association for the Advancement of Medical Instrumentation (AAMI) ST79, *Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities*. It is the end user's responsibility to use only sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization cassettes) that have been cleared by the FDA for the selected sterilization cycle specifications (time and temperature). Sterile implants meet pyrogen limit specifications.

When using a rigid sterilization container, the following must be taken into consideration for proper sterilization of Globus devices and loaded graphic cases:










- Recommended sterilization parameters are listed in the table below.
- Only FDA-cleared rigid sterilization containers for use with pre-vacuum steam sterilization may be used.
- When selecting a rigid sterilization container, it must have a minimum filter area of 176 in<sup>2</sup> total, or a minimum of four (4) 7.5in diameter filters.
- No more than one (1) loaded graphic case or its contents can be placed directly into a rigid sterilization container.
- Stand-alone modules/racks or single devices must be placed, without stacking, in a container basket to ensure optimal ventilation.
- The rigid sterilization container manufacturer's instructions for use are to be followed; if questions arise, contact the manufacturer of the specific container for guidance.
- Refer to AAMI ST79 for additional information concerning the use of rigid sterilization containers.

For implants and instruments provided NONSTERILE, sterilization is recommended (wrapped or containerized) as follows:

Method	Cycle Type	Temperature	Exposure Time	Drying Time
Steam	Pre-vacuum	132°C (270°F)	4 Minutes	30 Minutes

*Do not stack trays during sterilization. These parameters are validated to sterilize only this device. If other products are added to the sterilizer, the recommended parameters are not valid and new cycle parameters must be established by the user. The sterilizer must be properly installed, maintained, and calibrated. Ongoing testing must be performed to confirm inactivation of all forms of viable microorganisms.*

**CAUTION:** Federal (USA) Law Restricts this Device to Sale by or on the order of a Physician.

SYMBOL TRANSLATION			
	CATALOGUE NUMBER		STERILIZED BY IRRADIATION
	LOT NUMBER		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
	CAUTION		MANUFACTURER
	SINGLE USE ONLY		USE BY (YYYY-MM-DD)
	QUANTITY		

DI197 REV D



[illegible]

## NOTES

[illegible]



Globus Medical  
Valley Forge Business Center  
2560 General Armistead Avenue  
Audubon, PA 19403  
[www.globusmedical.com](http://www.globusmedical.com)

Customer Service:

Phone 1-866-GLOBUS1 (or 1-866-456-2871)

Fax 1-866-GLOBUS3 (or 1-866-456-2873)

©2021 Globus Medical. All rights reserved. Patent [www.globusmedical.com/patents](http://www.globusmedical.com/patents).  
Life moves us is a registered trademark of Globus Medical. Please refer to package insert for  
description, indications, contraindications, warnings, precautions and other important information.

CE 0297

GMTGD172  
09.21 Rev D