

RISE[®] TLIF

Expandable Lumbar Interbody Spacer



Our mission is to deliver cutting-edge technology, research, and innovative solutions to promote healing in patients with musculoskeletal disorders.

Life moves us 

The Surgical Technique shown is for illustrative purposes only. The technique(s) actually employed in each case always depends on the medical judgment of the surgeon exercised before and during surgery as to the best mode of treatment for each patient. Additionally, as instruments may occasionally be updated, the instruments depicted in this Surgical Technique may not be exactly the same as the instruments currently available. Please consult with your sales representative or contact Globus directly for more information.


SURGICAL TECHNIQUE GUIDE

RISE[®] TLIF

Overview	4
Implant Overview	6
Instrument Overview	7
Surgical Technique	
1. Transforaminal Access and Approach	11
2. Retraction	12
<i>Table Arm Attachment</i>	12
3. Creating Transforaminal Access	13
4. Discectomy/Endplate Preparation	14
5. Distraction and Implant Sizing	14
<i>Assembling the Adjustable Trial Assembly</i>	14
6. Implant Insertion	17
<i>Assembling the Implant Inserter</i>	17
<i>Implant Attachment</i>	18
7. Implant Expansion	19
8. Final Positioning	21
Optional: Implant Removal	21
Supplemental Fixation	22
Final Constructs	23
RISE [®] Implant Sets	24
RISE [®] Instrument Set	26
MARS [™] 3V Retractor Instrument Set	28
MARS [™] Instrument II Set	30
MARS [™] Instrument III Set	32
Posterior Disc Prep Instruments I Set	34
Posterior Disc Prep Instruments II Set	36
MIS Lumbar Discectomy Instrument Set	38
PRESERVE [®] Posterior Unilateral Instrument Set	40
Important Information	42

RISE® TLIF



Expandable Lumbar Interbody Spacer



The RISE® Transforaminal Lumbar Interbody Fusion (TLIF) spacer is an innovative expandable lumbar fusion device that represents the future of minimally invasive surgery by optimizing endplate-to-endplate fit and minimizing insertion force.

Insertion of RISE® TLIF is performed at a contracted height to help reduce the amount of nerve root retraction and musculoskeletal resection required.

The continuous expansion of RISE® TLIF is designed to restore disc height, while controlled distraction helps to properly tension the annulus and surrounding ligaments.

Contracted	Partially Expanded	Fully Expanded
		

Minimal Impaction

Contracted insertion height eases insertion into the disc space, which may help reduce tissue disruption and the amount of nerve root retraction required.




Controlled Disc Height Restoration

Controlled continuous expansion and distraction is designed to restore disc height and reduce the risk of over-distraction.








Optimized Fit

When expanded in the disc space, the implant optimizes endplate-to-endplate fit.



IMPLANT OVERVIEW












- Contracted insertion height helps minimize impaction force
- Controlled continuous expansion is designed to restore disc height
- Axial graft opening to help promote fusion
- Automatic locking for simple operation
- Convex profile designed to fit anatomy
- Multiple footprint options:
 - 8mm: 8x22mm
 - 10mm: 10x22, 10x26, and 10x30mm
 - 12mm: 12x26 and 12x30mm
- Height expansion options from 7–17mm
- Up to 7mm of expansion
- 4°, 10°, and 15° lordotic options
- One instrument for both insertion and expansion
- Slim design for minimally invasive surgery applications



INSTRUMENT OVERVIEW

SIZERS/SHIVERS



	Height	Part No.
	7mm	668.507
	8mm	668.508
	9mm	668.509
	10mm	668.510
	11mm	668.511
	12mm	668.512
	13mm	668.513
	14mm	668.514
	15mm	668.515
	16mm	668.516
	17mm	668.517



T-Handle 601.800

TRIALING INSTRUMENTS



MIS Handle 673.003



Adjustable Trial 10x26mm, 7-14mm 693.212



Adjustable Trial 10x26mm, 8-15mm, 15° 693.215



Torque-Limiting Palm Handle 694.002



Removable Drive, Right Hand 694.218



MIS Handle 673.003
Adjustable Trial 10x26mm, 7-14mm 693.212
Torque-Limiting Palm Handle 694.002
Removable Drive, Right Hand 694.218
(Assembled)

IMPLANT INSERTER COMPONENTS



Insertor Base 693.001



Torque-Limiting Implant Driver 693.002



Insertor Fork 8mm 693.011



Insertor Fork 10mm 693.012



Insertor Fork 12mm 693.013



Insertor Tube 693.021



Insertor QC Shaft 693.031



*Insertor Base 693.001
Torque-Limiting Implant Driver 693.002
Insertor Fork 10mm 693.012
Insertor Tube 693.021
Insertor QC Shaft 693.031
Torque-Limiting Palm Handle 694.002
(Assembled)*

OTHER INSTRUMENTS



Spanner Wrench 687.509



Slap Hammer 693.141

SURGICAL TECHNIQUE

RISE[®]

Please refer to the package insert, also printed at the back of this manual, for complete description, indications, contraindications, precautions, and warnings.

Minimally Invasive Surgery

Advances in minimally invasive surgery (MIS) are intended to lessen the disruption to the patient's anatomic structures. Without compromising surgical goals, MIS for interbody fusion has been shown to^{1,2}:

- Reduce soft tissue disruption
- Reduce postoperative pain
- Reduce blood loss
- Shorten hospital stay
- Reduce scarring
- Shorten recovery time


STEP

1

TRANSFORAMINAL ACCESS AND APPROACH

Approach

The patient is placed under anesthesia and positioned prone. Lateral C-arm fluoroscopy or other radiographic methods may be utilized throughout the surgery to ensure correct implant placement. In addition to the described interbody fusion technique, posterior stabilization, such as CREO MIS[®], must be used at the appropriate level(s).



The incision can be made 4 to 4.5cm lateral to the midline and the trajectory should be in line with the disc. Finger dissect between the multifidus and longissimus muscles until the facet joint is palpable.

1. Peng CW, Yue WM, et al. Clinical and Radiological Outcomes of Minimally Invasive Versus Open Transforaminal Lumbar Interbody Fusion. *SPINE* 34: 1385-1389, 2009.

2. Kim KT, Lee SH, et al. The Quantitative Analysis of Tissue Injury Markers After Mini-Open Lumbar Fusion. *SPINE* 6: 712-716, 2006.

STEP 2 RETRACTION

MARS™3V dilators may be used to retract soft tissue and to surround the facet. Keep downward pressure on the dilators and rotate as needed when approaching the facet. With the initial dilator in place, a series of cannulas are progressively passed over the initial dilator.

Ensure that the Retractor 3 Blade Frame is in the fully closed position and the blades are securely attached to the frame. Slide the retractor over the cannulas and apply gentle downward pressure on the frame.

Before removing the cannulas, articulate all three blades to one full rotation of the silver knobs. Articulating the blades in this manner will help prevent tissue creep as the cannulas are removed.

Once the retractor has been securely positioned and the Articulating Arm Assembly tightened to stabilize the retractor, remove the cannulas.

Use AP fluoroscopy to verify positioning.



⚙️ TABLE ARM ATTACHMENT


In order to use the MARS™3V Retractor, the Table Arm Attachment must be secured. Attach the Table Clamp onto the bed rail attachment. Insert the Articulating Arm Assembly into the Table Clamp and secure. The opposite end of the Articulating Arm Assembly is then attached to the Retractor 3 Blade Frame.

There are two options for attachment positions on the retractor, as shown at right.

Attaching the arm to point A maintains retractor position relative to the posterior blade position and translates the cephalad and caudad blades laterally when the retractor is opened.

Attaching the arm to point B maintains the retractor position relative to the cephalad and caudad blade position, and translates the posterior blade medially.

Refer to the *LLIF Surgical Technique Guide (GMTGD32)* for additional information.



STEP**3****CREATING TRANSFORAMINAL ACCESS**

Use an osteotome* to remove the inferior facet of the cephalad vertebrae and the superior facet of the caudal vertebrae at the appropriate level(s). This creates a working transforaminal access window to the disc.



Approach using MARS™ 3V
Retractor System



Transforaminal access


*Available in the Posterior Disc Prep Instruments Set II

STEP

4

DISCECTOMY/ENDPLATE PREPARATION

After creating an adequate annular window, remove disc material using rongeurs, rasps, curettes,* and other suitable preparation instruments. **Shavers** may be used to remove superficial layers of the cartilaginous endplates. Insert the smallest shaver into the disc space for further disc removal and endplate preparation, moving to larger shavers as needed. Use caution while using the shavers to avoid damage to the endplate. Careful disc removal and endplate preparation maximizes the potential for a successful fusion.



STEP


5

DISTRACTION AND IMPLANT SIZING




ASSEMBLING THE ADJUSTABLE TRIAL ASSEMBLY

1. Ensure that the Removable Drive, Right Hand is in the unlocked position by verifying that the tan push button is protruding beyond the cylinder.



2. Hold the shaft of the trial while turning the removable drive onto the threaded end until it stops and turn back 1/8th of a rotation to reduce tension. To lock, push the protruding push button on the removable drive until it is flush with the cylinder.



3. Thread the **MIS Handle** into the threaded hole, so the handle can be placed relative to the spine.
4. Attach the **Torque-Limiting Palm Handle**.

Insert the Adjustable Trial Assembly into the disc space at its contracted height. Expand the assembly gradually to the desired height by rotating the Torque-Limiting Palm Handle clockwise. Use caution while expanding to avoid excessive distraction and endplate damage.



Inserting adjustable trial



Trial height shown on side of trial

Rotate clockwise to expand trial to desired height

DISTRACTION AND IMPLANT SIZING (CONT'D)

Determine which height and length best fit the prepared disc space. A secure fit is desirable in order to maintain disc height and stabilize the segment. The final implant size may be confirmed using fluoroscopy.



Implant size confirmed

Note: Alternatively, the shaver/sizers may be used for distraction. Begin with the smallest shaver and use larger sizes until the desired distraction is achieved. Use caution while using shavers to avoid damage to the endplate.

STEP

6

IMPLANT INSERTION


Assemble the inserter as described below.

ASSEMBLING THE IMPLANT INSERTER

Position the instrument as shown below to assemble.

1. Partially thread the **Inserter QC Shaft** into the **Inserter Base**. Do not fully tighten the shaft as it will prevent the **Inserter Fork** from assembling properly in step 5.
2. Slide the **Inserter Tube**, with the flats on the distal end facing upwards, into the Inserter Base until it stops.
Note: The Inserter Tube may also be positioned with the cleaning windows facing upwards.
3. Thread the gold knob of the Inserter Base counterclockwise until it stops.
4. Select the proper size Inserter Fork that matches the implant to be used.
5. Insert the Inserter Fork into the Inserter Tube until it stops, making sure the tabs of the fork are aligned with the cleaning holes in the tube. The end of the shaft should be visible, halfway through the visualization hole in the Inserter Base (circled in red, below).
6. Tighten the Inserter QC Shaft to lock the fork in the holder.
7. Attach the Torque-Limiting Palm Handle, aligning up to the etched line on the Inserter QC Shaft.

Note: Do not torque out the Torque-Limiting Palm Handle during implant inserter assembly.




IMPLANT INSERTION (CONT'D)


Select an appropriate sized implant and pack autogenous bone graft material into the graft chamber.

IMPLANT ATTACHMENT

With the Inserter Tube threaded counterclockwise, attach the implant by placing it between the tips of the Inserter Fork. Rotate the tube clockwise until the implant is secure.



Insert the implant into the disc space using the Implant Inserter Assembly.




STEP

7

IMPLANT EXPANSION

Insert the **Torque-Limiting Implant Driver** through the Implant Inserter Assembly, engaging the drive screw in the implant. Rotate the driver clockwise to expand the implant to the appropriate height.



Expansion of the implant should be determined by the tactile feel of the implant in the disc space as it is expanded. This is determined by gently toggling the implant until the desired fit is achieved.

The overall height can be determined by counting the number of drive screw revolutions of the Torque-Limiting Implant Driver. Depending on the implant, approximately 1.75 revolutions equal 1mm of expansion. The arrow at the end of the driver may be used to help indicate revolutions.

Note: If inserting two RISE® implants, count the number of revolutions of the drive screw and/or confirm implant heights using fluoroscopy to ensure that both implants are expanded to the same height.



**Use arrow on back
end of driver to count
revolutions**

IMPLANT EXPANSION (CONT'D)

The Torque-Limiting Implant Driver is designed to allow the user to identify when either of the following conditions has occurred:

- The implant has reached its maximum height expansion

OR

- The implant is exerting the maximum allowable output of distraction force on the vertebral endplates.

Use caution while expanding the implant to avoid excessive distraction and endplate damage.

The maximum expansion of an implant can be verified by tracking the revolutions of the driver and referencing the expansion range chart printed on the implant module.

Drive Screw Revolutions Required											
Expansion Ranges	Final Height (mm)*										
	7	8	9	10	11	12	13	14	15	16	17
7-13mm	0	1.75	3.25	5.0	6.75	8.25	10.0	-	-	-	-
7-14mm	0	1.75	3.25	5.0	6.75	8.25	10.0	11.75	-	-	-
8-14mm	-	0	1.75	3.25	5.0	6.75	8.25	10.0	-	-	-
8-15mm	-	0	1.75	3.25	5.0	6.75	8.25	10.0	11.75	-	-
9-15mm	-	-	0	1.75	3.25	5.0	6.75	8.25	10.0	-	-
10-17mm	-	-	-	0	1.75	3.25	5.0	6.75	8.25	10.0	11.75
11-17mm	-	-	-	-	0	1.75	3.25	5.0	6.75	8.25	10.0



RISE® Inserter


STEP

8

FINAL POSITIONING

Using fluoroscopy, verify implant position before disengaging. Once the desired position is achieved, disengage the inserter from the implant by first removing the Torque-Limiting Implant Driver. Then rotate the gold knob on the Inserter Base counterclockwise until disengaged. The **Spanner Wrench** can be used to loosen the gold knob if necessary.

If repositioning is needed, the implant must be fully contracted before repositioning.



Removing Inserter

OPTIONAL: IMPLANT REMOVAL

For removal, the implant height may be reduced by inserting the Torque-Limiting Implant Driver into the back of the implant and rotating counterclockwise. Forceps or other manual surgical instruments may then be used to grasp and extract the implant. Alternatively, the implant may be rotated 90° and extracted using forceps or other manual surgical instruments.

SUPPLEMENTAL FIXATION

RISE® implants are intended for use with supplemental fixation, such as posterior stabilization. CREO MIS® may be used for posterior stabilization.

CREO MIS® STABILIZATION SYSTEM

Minimized Soft Tissue Disruption

Extended screw heads provide a small outer diameter to help reduce soft tissue disruption and screw sleeve interference.

Powerful MIS Correction


Deformity Adapters rigidly attach to the extended screw head and allow for screw manipulation and deformity correction.*

Integrated Rod Reduction

Ensures proper thread alignment while capturing, reducing, and locking the rod in one simplified step.



FINAL CONSTRUCTS



Lateral view



Posterior view

RISE®

IMPLANT SETS

RISE® 8mm Implant Set 993.904

Part No.	Description	Qty
193.001	RISE® Spacer 8x22mm, 7-13mm	2
193.051	RISE® Spacer 8x22mm, 9-15mm	2
193.011	RISE® Spacer 8x22mm, 11-17mm	2
193.021	RISE® Spacer 8x22mm, 8-14mm, 10°	2
193.041	RISE® Spacer 8x22mm, 11-17mm, 15°	2
993.004	RISE® 8mm Implant Module	

RISE® 10mm Implant Set 993.905

Part No.	Description	Qty
193.101	RISE® Spacer 10x22mm, 7-13mm	2
193.102	RISE® Spacer 10x26mm, 7-14mm	2
193.103	RISE® Spacer 10x30mm, 7-14mm	2
193.151	RISE® Spacer 10x22mm, 9-15mm	2
193.111	RISE® Spacer 10x22mm, 11-17mm	2
193.112	RISE® Spacer 10x26mm, 10-17mm	2
193.113	RISE® Spacer 10x30mm, 10-17mm	2
193.121	RISE® Spacer 10x22mm, 8-14mm, 10°	2
193.122	RISE® Spacer 10x26mm, 8-15mm, 10°	2
193.123	RISE® Spacer 10x30mm, 8-15mm, 10°	2
193.141	RISE® Spacer 10x22mm, 11-17mm, 15°	2
193.142	RISE® Spacer 10x26mm, 10-17mm, 15°	2
193.143	RISE® Spacer 10x30mm, 10-17mm, 15°	2
993.005	RISE® 10mm Implant Module	
993.907	RISE® Implant Graphic Case	

RISE® 12mm Implant Set 993.906

Part No.	Description	Qty
193.202	RISE® Spacer 12x26mm, 7-14mm	2
193.203	RISE® Spacer 12x30mm, 7-14mm	2
193.212	RISE® Spacer 12x26mm, 10-17mm	2
193.213	RISE® Spacer 12x30mm, 10-17mm	2
193.222	RISE® Spacer 12x26mm, 8-15mm, 10°	2
193.223	RISE® Spacer 12x30mm, 8-15mm, 10°	2
193.242	RISE® Spacer 12x26mm, 10-17mm, 15°	2
193.243	RISE® Spacer 12x30mm, 10-17mm, 15°	2
993.006	RISE® 12mm Implant Module	

RISE® IMPLANT SETS





RISE®

INSTRUMENT SET 993.903

	Instrument	Qty
1	601.800 T-Handle	2
2	668.507 Sizer/Shaver, 7mm	1
3	668.508 Sizer/Shaver, 8mm	1
4	668.509 Sizer/Shaver, 9mm	1
5	668.510 Sizer/Shaver, 10mm	1
6	668.511 Sizer/Shaver, 11mm	1
7	668.512 Sizer/Shaver, 12mm	1
8	668.513 Sizer/Shaver, 13mm	1
9	668.514 Sizer/Shaver, 14mm	1
10	668.515 Sizer/Shaver, 15mm	1
11	668.516 Sizer/Shaver, 16mm	1
12	668.517 Sizer/Shaver, 17mm	1
13	673.003 MIS Handle	1
14	687.509 Spanner Wrench	1
15	693.001 Inserter Base	2
16	693.002 Torque-Limiting Implant Driver	2
17	693.011 Inserter Fork, 8mm	2
18	693.012 Inserter Fork, 10mm	2
19	693.013 Inserter, 12mm	2
20	693.021 Inserter Tube	2
21	693.031 Inserter QC Shaft	2
22	693.141 Slap Hammer	1
23	693.212 Adjustable Trial 10x26, 7-14mm	1
24	693.215 Adjustable Trial 10x26, 8-15mm, 15°	1
25	694.002 Torque-Limiting Palm Handle	1
26	694.218 Removable Drive, Right Hand	1
	993.003 RISE® Instrument Graphic Case	

RISE[®]

INSTRUMENT SET 993.903






MARS™3V

RETRACTOR SET 998.901

Instrument			Qty	Retractor Blades (Cont'd)		Qty
1	623.003	K-Wire Gripper	1	698.510	Blade, CC, 40mm	2
2	698.100	Retractor 3 Blade Frame	1	698.512	Blade, CC, 50mm	2
3	632.102	Retractor 2 Blade Frame	1	698.514	Blade, CC, 60mm	2
4	632.150	10mm Socket Driver	1	698.516	Blade, CC, 70mm	2
5	698.250	Hook and Latch Driver	1	698.518	Blade, CC, 80mm	2
6	675.403	Flouro Modulator	1	698.520	Blade, CC, 90mm	2
7	675.404	Incision Locator	1	698.522	Blade, CC, 100mm	2
8	675.513	8" Suction	1	698.524	Blade, CC, 110mm	2
9	675.800	Radiolucent Initial Dilator Holder	1	698.526	Blade, CC, 120mm	2
10	698.205	Cannula A	1	698.528	Blade, CC, 130mm	2
11	698.210	Cannula B	1	698.530	Blade, CC, 140mm	2
12	698.215	Cannula C	1	698.532	Blade, CC, 150mm	2
13	698.220	Cannula D	1	698.534	Blade, CC, 160mm	2
14	698.230	Frame Handle	1	698.536	Blade, CC, 170mm	2
15	698.240	Shim Tool, CC	1	Disposables		Qty
16	698.260	Docking Pin Tool	1			
17	698.330	Disc Shim Tool	1	632.678S	Bipolar Forceps, 10" Bayo, 1.0mm Tip	1
18	698.350	Docking Pin Sleeve	4	698.600S	MARS™3V Disposable Kit	1
19	Retractor Blades		Qty	698.300S	Lengthening Shim	2
	698.450	Blade, Posterior, 40mm	2	698.305S	Widening Shim	2
	698.452	Blade, Posterior, 50mm	2	698.310S	Docking Pin, 10mm	2
	698.454	Blade, Posterior, 60mm	2	698.315S	Docking Pin, 20mm	2
	698.456	Blade, Posterior, 70mm	2	698.325S	Disc Shim, Aluminum	1
	698.458	Blade, Posterior, 80mm	2	698.326S	Disc Shim, Stainless Steel	
	698.460	Blade, Posterior, 90mm	2			
	698.462	Blade, Posterior, 100mm	2			
	698.464	Blade, Posterior, 110mm	2			
	698.466	Blade, Posterior, 120mm	2			
	698.468	Blade, Posterior, 130mm	2			
	698.470	Blade, Posterior, 140mm	2			
	698.472	Blade, Posterior, 150mm	2			
	698.474	Blade, Posterior, 160mm	2			
	698.476	Blade, Posterior, 170mm	2			

MARS™ 3V

RETRACTOR SET 998.901





MARS™

INSTRUMENT II SET 932.902

	Instrument	Qty
1	632.300 Fiber-Optic Cord	1
2	632.305 Adapter, ACMI	1
3	632.306 Adapter, Wolf	1
4	632.307 Adapter, Olympus	1
5	632.308 Adapter, Storz	1
6	632.310S Light Cable	1
7	632.390 Port Lock, 19mm	1
8	632.391 Port Lock, 22mm	1
9	632.392 Port Lock, 26mm	1
10	632.401 2mm Cannula	1
11	632.402 5mm Cannula	1
12	632.403 8mm Cannula	1
13	632.404 12mm Cannula	1
14	632.405 15mm Cannula	1
15	632.406 18mm Cannula	1
16	632.407 22mm Cannula	1
17	632.408 26mm Port Mount	1
18	632.409 26mm Port Positioner	1
19	632.410 22mm Port Mount	1
20	632.411 22mm Port Positioner	1
21	632.412 19mm Port Mount	1
22	632.413 19mm Port Positioner	1
23	632.500 Table Clamp	1
24	632.750 Articulating Arm Assembly	1
	932.002 MARS™ Instrument II Graphic Case	

MARS™

INSTRUMENT II SET 932.902




MARS™

INSTRUMENT III SET 932.903

Instrument			Qty	Instrument			Qty
1	632.600	Pituitary, 2mm Bayoneted	1	14	632.650	Penfield #2 Push, Bayoneted	1
	632.601	Pituitary, 2mm, Down-Biting, Bayoneted	1	15	632.651	Penfield #2 Pull, Bayoneted	1
	632.602	Pituitary, 2mm, Up-Biting, Bayoneted	1	16	632.652	Penfield #4 Push, Bayoneted	1
	632.605	Pituitary, 4mm, Up-Biting, Bayoneted	1	17	632.653	Penfield #4 Pull, Bayoneted	1
2	632.610	Micro Pituitary, 2mm, Up-Biting, Bayoneted	1	18	632.655	Nerve Hook, Straight, Bayoneted	1
	632.611	Micro Pituitary, 2mm, Bayoneted	1	19	632.656	Nerve Hook, Left, Bayoneted	1
3	632.615	Pituitary, Ring Handle, 2mm	1	20	632.657	Nerve Hook, Right, Bayoneted	1
4	632.616	Scissors, Straight	1	21	632.660	90° Dissector, Straight, Bayoneted	1
	632.618	Scissors, Curved Left	1	22	632.661	90° Dissector, Left, Bayoneted	1
	632.619	Scissors, Curved Right	1	23	632.662	90° Dissector, Right, Bayoneted	1
5	632.620	Kerrison 40°, 3mm, Bayoneted	1	24	632.673	Nerve Root Retractor	1
	632.621	Kerrison 90°, 3mm, Bayoneted	1	25	632.674	Nerve Root Retractor, Wide	1
	632.622	Kerrison 40°, 4mm, Bayoneted	1	26	632.675	Suction Retractor	1
	632.623	Kerrison 90°, 4mm, Bayoneted	1	27	632.676	Bi-Polar Forcep, Straight, Bayoneted, US Connection	1
	632.624	Kerrison 40°, 5mm, Bayoneted	1		632.677	Bi-Polar Forcep, Angled, Bayoneted, US Connection	1
	632.625	Kerrison 90°, 5mm, Bayoneted	1		932.003	MARS™ Instrument Graphic Case	
6	632.630	Bone Curette Straight, 5.2 Cup, Bayoneted	1				
	632.631	Bone Curette Straight, 3.6 Cup, Bayoneted	1				
	632.632	Bone Curette Angled, 5.2 Cup, Bayoneted	1				
	632.633	Bone Curette Angled, 3.6 Cup, Bayoneted	1				
	632.634	Bone Curette Reverse Angled, 5.2 Cup, Bayoneted	1				
	632.635	Bone Curette Reverse Angled, 3.6 Cup, Bayoneted	1				
7	632.640	Woodson Probe	1				
8	632.641	90° Ball Probe Short, Straight, Bayoneted	1				
9	632.642	90° Ball Probe Short, Left, Bayoneted	1				
10	632.643	90° Ball Probe Short, Right, Bayoneted	1				
11	632.644	90° Ball Probe Long, Straight, Bayoneted	1				
12	632.645	90° Ball Probe Long, Left, Bayoneted	1				
13	632.646	90° Ball Probe Long, Right, Bayoneted	1				


MARS™
INSTRUMENT III SET 932.903



POSTERIOR DISC PREP INSTRUMENTS I SET 926.901

	Instrument		Qty
1	626.210	Push Rod Assembly, Bone Funnel	1
2	626.215	Nerve Retractor, 5mm, Suction	1
3	626.220	Nerve Retractor, Corner	1
4	603.061	Nerve Root Retractor, Fine, 5mm	1
5	603.062	Nerve Root Retractor, Medium, 10mm	1
6	679.015	Bone Funnel	1
7	679.015	Bone Funnel - Tube	1
8	626.235	Disc Rongeur, 250x2mm, Straight	1
9	626.236	Disc Rongeur, 250x2mm, Up-Biting	1
10	626.240	Disc Rongeur, 250x4mm, Straight	1
11	626.241	Disc Rongeur, 250x6mm, Straight	1
12	626.242	Disc Rongeur, 250x4mm, Up-Biting	1
13	626.243	Disc Rongeur, 250x4mm, Down-Biting	1
14	626.250	Kerrison, 250x3mm, Straight	1
15	626.252	Kerrison, 250x5mm, Straight	1
16	626.260	Lamina Spreader, Hinged	1
	926.102	Graphic Case	


POSTERIOR DISC PREP INSTRUMENTS I SET 926.901



POSTERIOR DISC PREP INSTRUMENTS II SET 926.902

	Instrument	Qty
1	626.150 Bone Curette, 6.5x9.5mm, Straight	1
2	626.151 Bone Curette, 6.5x9.5mm, Right	1
3	626.152 Bone Curette, 6.5x9.5mm, Left	1
4	626.153 Bone Curette, 6.5x9.5mm, Up-Pushing	1
5	626.154 Bone Curette, 6.5x9.5mm, Down-Pushing	1
6	626.190 Rake, 8mm, Straight	1
7	626.191 Rake, 8mm, Angled	1
8	626.140 Bone Curette, 5.0x7.5mm, Straight	1
9	626.143 Bone Curette, 5.0x7.5mm, Up-Pushing	1
10	626.144 Bone Curette, 5.0x7.5mm, Down-Pushing	1
11	626.160 Bone Curette, 8.0x11.5mm, Straight	1
12	626.161 Bone Curette, 8.0x11.5mm, Right	1
13	626.162 Bone Curette, 8.0x11.5mm, Left	1
14	626.170 Bone Curette, 5.0x10mm, Axial	1
15	626.180 Osteotome, 7mm	1
16	626.185 Rasp, 8x20mm, Knurled, Straight	1
17	626.186 Rasp, 8x20mm, Knurled, Angled	1
18	626.200 Ring Curette, 6mm, Straight	1
19	626.201 Ring Curette, 6mm, Angled Right	1
20	626.202 Ring Curette, 6mm, Angled Left	1
	926.101 Graphic Case II	


POSTERIOR DISC PREP INSTRUMENTS II SET 926.902



MIS LUMBAR DISCECTOMY INSTRUMENT SET 979.901

	Instrument	Qty	Additionally Available
1	679.005 Elevator 6mm	1	673.018 Push Rod, Bone Funnel
2	679.007 Osteotome, 8mm QR	1	679.021 Bone Curette, Angled, 10.7 Serrated Cup
3	679.008 Osteotome, 6mm QR	1	679.022 Bone Curette, Straight, 10.7 Serrated Cup
4	679.010 L-Handle	1	679.023 Bone Curette, Angled, 10.7 Serrated Cup
5	679.011 Rake, Long 8mm, Bayoneted	1	679.024 Bone Curette, Straight, 10.7 Serrated Cup
6	679.012 Rake, Long 9mm, Bayoneted	1	679.061 Bone Curette, 10.0 Rectangle Cup, 75° Up
7	679.013 Rake, Long 10mm, Bayoneted	1	679.062 Bone Curette, 10.0 Rectangle Cup, 75° Down
8	679.015 Bone Funnel	1	679.063 Bone Curette, 12.0 Rectangle Cup, 75° Up
9	679.016 Bone Funnel Clamp	1	
10	679.017 Bone Pusher Rod	1	913.001 MIS Lumbar Discectomy Graphic Case
11	679.025 Bone Curette, 10.0 Serrated Cup	1	
	679.026 Bone Curette, Straight, 10.0 Serrated Cup	1	
	679.027 Bone Curette, Angled, 10.0 Serrated Cup	1	
	679.028 Bone Curette, Straight, 10.0 Serrated Cup	1	
	679.031 Bone Curette, Angled, 12.0 Serrated Cup, Bayoneted, LH	1	
	679.032 Bone Curette, Straight, 12.0 Serrated Cup, Bayoneted, LH	1	
	679.033 Bone Curette, Angled, 12.0 Serrated Cup, Bayoneted, RH	1	
	679.034 Bone Curette, Straight, 12.0 Serrated Cup, Bayoneted, RH	1	
12	679.041 Bone Curette, 10.7 Serrated Cup, 90° Up	1	
13	679.042 Bone Curette, 10.7 Serrated Cup, 90° Down	1	
14	679.051 Ring Curette, 6mm	1	
	979.001 MIS Discectomy Instruments Graphic Case		


MIS LUMBAR DISCECTOMY INSTRUMENT SET 979.901



PRESERVE® POSTERIOR UNILATERAL INSTRUMENT SET 904.907

Instrument	Qty	Instrument	Qty
1 601.001 Implant Holder	1	13 604.317 Scraper, Oblique, 17mm	1
2 601.800 T-Handle	2	14 604.808 Paddle Distractor, 8mm	1
3 604.001 Holder, Straight	1	15 604.809 Paddle Distractor, 9mm	1
4 604.002 Holder, Angled	1	16 604.810 Paddle Distractor, 10mm	1
5 604.108 Trial, SUSTAIN®-R Oblique, 26mm length, 8mm	1	17 604.811 Paddle Distractor, 11mm	1
604.109 Trial, SUSTAIN®-R Oblique, 26mm length, 9mm	1	18 604.812 Paddle Distractor, 12mm	1
604.110 Trial, SUSTAIN®-R Oblique, 26mm length, 10mm	1	19 604.813 Paddle Distractor, 13mm	1
604.111 Trial, SUSTAIN®-R Oblique, 26mm length, 11mm	1	20 604.815 Paddle Distractor, 15mm	1
604.112 Trial, SUSTAIN®-R Oblique, 26mm length, 12mm	1	21 604.817 Paddle Distractor, 17mm	1
604.113 Trial, SUSTAIN®-R Oblique, 26mm length, 13mm	1	22 673.017 Slide Hammer, Quick Disconnect	1
604.115 Trial, SUSTAIN®-R Oblique, 26mm length, 15mm	1	23 673.108 Trial Shaft, 8x22mm wide, 8mm, SUSTAIN® Oblique, Small	1
604.117 Trial, SUSTAIN®-R Oblique, 26mm length, 17mm	1	673.109 Trial Shaft, 8mm wide, 9mm, SUSTAIN® Oblique, Small	1
604.208 Trial, SUSTAIN®-R Oblique, 30mm length, 8mm	1	673.110 Trial Shaft, 8mm wide, 10mm, SUSTAIN® Oblique, Small	1
604.209 Trial, SUSTAIN®-R Oblique, 30mm length, 9mm	1	673.111 Trial Shaft, 8x22mm wide, 11mm, SUSTAIN® Oblique, Small	1
604.210 Trial, SUSTAIN®-R Oblique, 30mm length, 10mm	1	673.112 Trial Shaft, 8x22mm wide, 12mm, SUSTAIN® Oblique, Small	1
604.211 Trial, SUSTAIN®-R Oblique, 30mm length, 11mm	1	673.113 Trial Shaft, 8x22mm wide, 13mm, SUSTAIN® Oblique, Small	1
604.212 Trial, SUSTAIN®-R Oblique, 30mm length, 12mm	1	673.115 Trial Shaft, 8x22mm wide, 15mm, SUSTAIN® Oblique, Small	1
604.213 Trial, SUSTAIN®-R Oblique, 30mm length, 13mm	1	673.117 Trial Shaft, 8x22mm wide, 17mm, SUSTAIN® Oblique, Small	1
604.215 Trial, SUSTAIN®-R Oblique, 30mm length, 15mm	1	673.208 Trial Shaft, 10x22mm wide, 8mm, SUSTAIN® Oblique, Small	1
604.217 Trial, SUSTAIN®-R Oblique, 30mm length, 17mm	1	673.209 Trial Shaft, 10x22mm wide, 9mm, SUSTAIN® Oblique, Small	1
6 604.308 Scraper, Oblique, 8mm	1	673.210 Trial Shaft, 10x22mm wide, 10mm, SUSTAIN® Oblique, Small	1
7 604.309 Scraper, Oblique, 9mm	1	673.211 Trial Shaft, 10x22mm wide, 11mm, SUSTAIN® Oblique, Small	1
8 604.310 Scraper, Oblique, 10mm	1	673.212 Trial Shaft, 10x22mm wide, 12mm, SUSTAIN® Oblique, Small	1
9 604.311 Scraper, Oblique, 11mm	1	673.213 Trial Shaft, 10x22mm wide, 13mm, SUSTAIN® Oblique, Small	1
10 604.312 Scraper, Oblique, 12mm	1	673.215 Trial Shaft, 10x22mm wide, 15mm, SUSTAIN® Oblique, Small	1
11 604.313 Scraper, Oblique, 13mm	1	673.217 Trial Shaft, 10x22mm wide, 17mm, SUSTAIN® Oblique, Small	1
12 604.315 Scraper, Oblique, 15mm	1	904.009 PRESERVE® Posterior Unilateral Instruments	

PRESERVE® POSTERIOR UNILATERAL INSTRUMENT SET 904.907



IMPORTANT INFORMATION ON THE RISE® SPACER

DESCRIPTION

RISE® Spacers are lumbar interbody fusion devices used to provide structural stability in skeletally mature individuals following discectomy. RISE® Spacers are provided in different shapes to accommodate various surgical approaches to the lumbar spine (posterior, transforaminal [posterolateral] or lateral) and can expand to the desired height. The implants are available in various heights and geometric options to fit the anatomical needs of a wide variety of patients. This device is to be filled with autograft bone and/or allogenic bone graft composed of cancellous and/or corticocancellous bone. Protrusions on the superior and inferior surfaces of each device grip the endplates of the adjacent vertebrae to resist expulsion.

RISE® Spacers are manufactured from titanium alloy, as specified in ASTM F136 and F1295. An internal component is manufactured from radiolucent PEEK polymer, as specified in ASTM F2026.

INDICATIONS

The RISE® Spacer is an interbody fusion device intended for use at one or more levels of the thoracic spine (T1-T12), thoracolumbar junction (T12-L1), or lumbosacral spine (L1-S1) as an adjunct to fusion in patients with the following indications: degenerative disc disease (DDD), disc herniation (with myelopathy and/or radiculopathy), spondylolisthesis, deformity (degenerative scoliosis or kyphosis), spinal stenosis, and failed previous fusion (pseudarthrosis). DDD is defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies. These patients should be skeletally mature and have had at least six (6) months of non-operative treatment.

The RISE® Spacer is to be filled with autograft bone and/or allogenic bone graft composed of cancellous and/or corticocancellous bone. This device is intended to be used with supplemental fixation systems that have been cleared for use in the thoracolumbosacral spine (e.g., posterior pedicle screw and rod systems, anterior plate systems, and anterior screw and rod systems).

WARNINGS

One of the potential risks identified with this system is death. Other potential risks which may require additional surgery, include:

- device component fracture,
- loss of fixation,
- non-union,
- fracture of the vertebrae,
- neurological injury, and
- vascular or visceral injury.

Interbody fusion devices for the treatment of degenerative conditions are designed to withstand both full load bearing and the loads associated with long-term use which could result from the presence of non-union or delayed union.

Certain degenerative diseases or underlying physiological conditions such as diabetes, rheumatoid arthritis, or osteoporosis may alter the healing process, thereby increasing the risk of implant breakage or spinal fracture.

Possible adverse effects which may occur include: failed fusion or pseudarthrosis leading to implant breakage; allergic reaction to implant materials; device fracture or failure; device migration or loosening; decrease in bone density; pain, discomfort, or abnormal sensations due to the presence of the device; injury to nerves, vessels, and organs; venous thrombosis, lung embolism and cardiac arrest; and death.

Components of this system are manufactured from titanium alloy. Dissimilar metals in contact with each other can accelerate the corrosion process due to galvanic corrosion effects. Mixing of implant components with different materials is not recommended, for metallurgical, mechanical, and functional reasons.

These warnings do not include all adverse effects which could occur with surgery in general, but are important considerations particular to orthopedic implants. General surgical risks should be explained to the patient prior to surgery.

Patients with previous spinal surgery at the level(s) to be treated may have different clinical outcomes compared to those without previous surgery.

PRECAUTIONS

The implantation of intervertebral fusion devices should be performed only by experienced spinal surgeons with specific training in the use of this system because this is a technically demanding procedure presenting a risk of serious injury to the patient. Preoperative planning and patient anatomy should be considered when selecting implant size.

Surgical implants must never be reused. An explanted implant must never be reimplanted. Even though the device appears undamaged, it may have small defects and internal stress patterns which could lead to breakage.

Adequately instruct the patient. Mental or physical impairment which compromises or prevents a patient's ability to comply with necessary limitations or precautions may place that patient at a particular risk during postoperative rehabilitation.

Metallic implants can loosen, fracture, corrode, migrate, cause pain, or stress shield bone even after a fracture has healed, particularly in young, active patients. While

the surgeon must have the final decision on implant removal, we recommend that whenever possible and practical for the individual patient, fixation devices should be removed once their service as an aid to healing is accomplished. Implant removal should be followed by adequate postoperative management.

Factors such as the patient's weight, activity level, and adherence to weight bearing or load bearing instructions have an effect on the stresses to which the implant is subjected.

For optimal implant performance, the surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact the performance of the system.

MRI SAFETY INFORMATION



The RISE® Spacers are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla only
- Maximum spatial field gradient of 3,000 gauss/cm (30 T/m) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 1 W/kg

Under the scan conditions defined above, the RISE® Spacers are expected to produce a maximum temperature rise of less than or equal to 3.9°C after 15 minutes of continuous scanning.

The image artifact caused by the device is not expected to extend beyond 35mm from the device when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

CONTRAINDICATIONS

Use of the RISE® spacer is contraindicated in patients with the following conditions:

1. Active systemic infection, infection localized to the site of the proposed implantation, or when the patient has demonstrated allergy or foreign body sensitivity to any of the implant materials
2. Prior fusion at the level(s) to be treated
3. Severe osteoporosis, which may prevent adequate fixation
4. Conditions that may place excessive stresses on bone and implants, such as severe obesity or degenerative diseases, are relative contraindications. The decision whether to use these devices in such conditions must be made by the physician taking into account the risks versus the benefits to the patient.
5. Patients whose activity, mental capacity, mental illness, alcoholism, drug abuse, occupation, or lifestyle may interfere with their ability to follow postoperative restrictions and who may place undue stresses on the implant during bony healing and may be at a higher risk of implant failure.
6. Any condition not described in the indications for use
7. Signs of local inflammation
8. Fever or leukocytosis
9. Morbid obesity
10. Pregnancy
11. Mental illness
12. Any other condition which would preclude the potential benefit of spinal implant surgery, such as the presence of tumors or congenital abnormalities, fracture local to the operating site, elevation of sedimentation rate unexplained by other diseases, elevation of the white blood count (WBC), or a marked left shift in the WBC differential count
13. Suspected or documented allergy or intolerance to composite materials
14. Any case not needing a fusion
15. Any patient not willing to cooperate with postoperative instruction
16. Patients with a known hereditary or acquired bone friability or calcification problem should not be considered for this type of surgery
17. These devices must not be used for pediatric cases, nor where the patient still has general skeletal growth
18. Spondylolisthesis unable to be reduced to Grade 1
19. Any case where the implant components selected for use would be too large or too small to achieve a successful result
20. Any case that requires the mixing of metals from two different components or systems
21. Any patient having inadequate tissue coverage at the operative site or inadequate bone stock or quality
22. Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance

COMPLICATIONS AND POSSIBLE ADVERSE EVENTS

Prior to surgery, patients should be made aware of the following possible adverse effects in addition to the potential need for additional surgery to correct these effects:

- Loosening, bending or breakage of components
- Displacement/migration of device components
- Tissue sensitivity to implant material
- Potential for skin breakdown and/or wound complications
- Non-union or delayed union or mal-union
- Infection

IMPORTANT INFORMATION ON THE RISE® SPACER

- Nerve damage, including loss of neurological function (sensory and/or motor), paralysis, dysesthesia, hyperesthesia, paresthesia, radiculopathy, reflex deficit, cauda equina syndrome
- Dural tears, cerebral spinal fluid leakage
- Fracture of vertebrae
- Foreign body reaction (allergic) to components or debris
- Vascular or visceral injury
- Change in spinal curvature, loss of correction, height and/or reduction
- Urinary retention or loss of bladder control or other types of disorders of the urogenital system
- Ileus, gastritis, bowel obstruction or other types of gastrointestinal system compromise
- Reproductive system compromise including impotence, sterility, loss of consortium and sexual dysfunction.
- Pain or discomfort
- Bursitis
- Decrease in bone density due to stress shielding
- Loss of bone or fracture of bone above or below the level of surgery
- Bone graft donor site pain, fracture, and/or delayed wound healing
- Restriction of activities
- Lack of effective treatment of symptoms for which surgery was intended
- Need for additional surgical intervention
- Death

PACKAGING

These implants and instruments may be supplied pre-packaged and sterile, using gamma irradiation. The integrity of the sterile packaging should be checked to ensure that sterility of the contents is not compromised. Packaging should be carefully checked for completeness and all components should be carefully checked to ensure that there is no damage prior to use. Damaged packages or products should not be used, and should be returned to Globus Medical. During surgery, after the correct size has been determined, remove the products from the packaging using aseptic technique.

The instrument sets are provided nonsterile and are steam sterilized prior to use, as described in the STERILIZATION section below. Following use or exposure to soil, instruments must be cleaned, as described in the CLEANING section below.

HANDLING AND USE

All instruments and implants should be treated with care. Improper use or handling may lead to damage and/or possible malfunction. Products should be checked to ensure that they are in working order prior to surgery. All products should be inspected prior to use to ensure that there is no unacceptable deterioration such as corrosion, discoloration, pitting, cracked seals, etc. Non-working or damaged instruments should not be used, and should be returned to Globus Medical.

Implants are single use devices and should not be cleaned. Re-cleaning of single use implants might lead to mechanical failure and/or material degradation. Discard any implants that may have been accidentally contaminated.

CLEANING

All instruments that can be disassembled must be disassembled for cleaning. All handles must be detached. Instruments may be reassembled following sterilization. The instruments should be cleaned using neutral cleaners before sterilization and introduction into a sterile surgical field or (if applicable) return of the product to Globus Medical.

Cleaning and disinfecting of instruments can be performed with aldehyde-free solvents at higher temperatures. Cleaning and decontamination must include the use of neutral cleaners followed by a deionized water rinse. Note: certain cleaning solutions such as those containing formalin, glutaraldehyde, bleach and/or other alkaline cleaners may damage some devices, particularly instruments; these solutions should not be used.

The following cleaning methods should be observed when cleaning instruments after use or exposure to soil, and prior to sterilization:

1. Immediately following use, ensure that the instruments are wiped down to remove all visible soil and kept from drying by submerging or covering with a wet towel.
2. Disassemble all instruments that can be disassembled.
3. Rinse the instruments under running tap water to remove all visible soil. Flush the lumens a minimum of 3 times, until the lumens flush clean.
4. Prepare Enzo® (or a similar enzymatic detergent) per manufacturer's recommendations.
5. Immerse the instruments in the detergent and allow them to soak for a minimum of 2 minutes.
6. Use a soft bristled brush to thoroughly clean the instruments. Use a pipe cleaner for any lumens. Pay close attention to hard to reach areas.
7. Using a sterile syringe, draw up the enzymatic detergent solution. Flush any lumens and hard to reach areas until no soil is seen exiting the area.
8. Remove the instruments from the detergent and rinse them in running warm tap water.
9. Prepare Enzo® (or a similar enzymatic detergent) per manufacturer's recommendations in an ultrasonic cleaner.

10. Completely immerse the instruments in the ultrasonic cleaner and ensure detergent is in lumens by flushing the lumens. Sonicate for a minimum of 3 minutes.
11. Remove the instruments from the detergent and rinse them in running deionized water or reverse osmosis water for a minimum of 2 minutes.
12. Dry instruments using a clean soft cloth and filtered pressurized air.
13. Visually inspect each instrument for visible soil. If visible soil is present, then repeat cleaning process starting with Step 3.

CONTACT INFORMATION

Globus Medical may be contacted at 1-866-GLOBUS1 (456-2871). A surgical technique manual may be obtained by contacting Globus Medical.

STERILIZATION

These implants and instruments may be available sterile or nonsterile.

Sterile implants and instruments are sterilized by gamma radiation, validated to ensure a Sterility Assurance Level (SAL) of 10⁻⁶. Sterile products are packaged in a heat sealed container/pouch. The expiration date is provided in the package label. These products are considered sterile unless the packaging has been opened or damaged. Sterile implants and instruments that become nonsterile or have expired packaging are considered nonsterile and may be sterilized according to instructions for nonsterile implants and instruments below. Sterile implants meet pyrogen limit specifications.

Nonsterile implants and instruments have been validated to ensure an SAL of 10⁻⁶. The use of an FDA-cleared wrap is recommended, per the Association for the Advancement of Medical Instrumentation (AAMI) ST79, *Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities*. It is the end user's responsibility to use only sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization cassettes) that have been cleared by the FDA for the selected sterilization cycle specifications (time and temperature).

When using a rigid sterilization container, the following must be taken into consideration for proper sterilization of Globus devices and loaded graphic cases:










- Recommended sterilization parameters are listed in the table below.
- Only FDA-cleared rigid sterilization containers for use with pre-vacuum steam sterilization may be used.
- When selecting a rigid sterilization container, it must have a minimum filter area of 176 in² total, or a minimum of four (4) 7.5in diameter filters.
- No more than one (1) loaded graphic case or its contents can be placed directly into a rigid sterilization container.
- Stand-alone modules/racks or single devices must be placed, without stacking, in a container basket to ensure optimal ventilation.
- The rigid sterilization container manufacturer's instructions for use are to be followed; if questions arise, contact the manufacturer of the specific container for guidance.
- Refer to AAMI ST79 for additional information concerning the use of rigid sterilization containers.

For implants and instruments provided NONSTERILE, sterilization is recommended (wrapped or containerized) as follows:

Method	Cycle Type	Temperature	Exposure Time	Drying Time
Steam	Pre-vacuum	132°C (270°F)	4 Minutes	30 Minutes
Steam	Pre-vacuum	134°C (273°F)	3 Minutes	30 Minutes

These parameters are validated to sterilize only this device. If other products are added to the sterilizer, the recommended parameters are not valid and new cycle parameters must be established by the user. The sterilizer must be properly installed, maintained, and calibrated. Ongoing testing must be performed to confirm inactivation of all forms of viable microorganisms.

CAUTION: Federal Law (USA) Restricts this Device to Sale by or on the order of a Physician.

SYMBOL TRANSLATION			
	CATALOGUE NUMBER		STERILIZED BY IRRADIATION
	LOT NUMBER		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
	CAUTION		MANUFACTURER
	SINGLE USE ONLY		USE BY (YYYY-MM-DD)
	QUANTITY		



Globus Medical
Valley Forge Business Center
2560 General Armistead Avenue
Audubon, PA 19403
www.globusmedical.com

Customer Service:

Phone 1-866-GLOBUS1 (or 1-866-456-2871)

Fax 1-866-GLOBUS3 (or 1-866-456-2873)

©2023 Globus Medical. All rights reserved. Patent www.globusmedical.com/patents.
Life moves us is a registered trademark of Globus Medical. Please refer to package insert for
description, indications, contraindications, warnings, precautions and other important information.

CE 0297

GMTGD85
08.23 Rev D