
M.flow[®] Instruments

- DE Gebrauchsanweisung | EN Instructions for use | FR Mode d'emploi |
ES Instrucciones de manejo | PT Instruções de utilização | IT Istruzioni per l'uso |
NL Gebruiksaanwijzing |

 www.miethke.com

US This Instructions for Use is NOT intended for United States users. Please discard.
The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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info@miethke.com

CONTENTS

1	PREFACE AND IMPORTANT INFORMATION	18
2	INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE	18
2.1	EXPLANATION OF THE WARNINGS	18
2.2	DISPLAY CONVENTIONS	18
2.3	OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION	18
2.4	FEEDBACK ON THE INSTRUCTIONS FOR USE	18
2.5	COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION	18
3	DESCRIPTION <i>M.flow Instruments</i>	19
3.1	MEDICAL PURPOSE	19
3.2	CLINICAL BENEFITS	19
3.3	INDICATIONS	19
3.4	CONTRAINDICATIONS	19
3.5	INTENDED PATIENT GROUPS	19
3.6	INTENDED USERS	19
3.7	INTENDED USE ENVIRONMENT	19
3.8	TECHNICAL DESCRIPTION	19
3.9	SYSTEM COMPONENTS	20
3.10	FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	20
4	PROPERTIES <i>M.flow Instruments</i>	20
4.1	PRODUCT DESCRIPTION	20
4.2	IMPORTANT SAFETY INFORMATION	21
4.3	PATIENT EDUCATION	21
4.4	TRANSPORT AND STORAGE	22
4.5	USE OF THE PRODUCT	22
4.6	DISPOSAL	28
4.7	TROUBLESHOOTING	28
4.8	TECHNICAL INFORMATION	29
4.9	SYMBOLS USED FOR LABELLING	30
5	MEDICAL DEVICE CONSULTANTS	30

1 PREFACE AND IMPORTANT INFORMATION

Preface

Thank you for purchasing the medical device *M.flow Instruments*. Please contact us if you have any questions about the contents of these instructions for use or the use of the product. Your team at Christoph Miethke GmbH & Co. KG

Relevance of the instructions for use

**WARNING**

Improper handling and non-intended use of this product can cause risks and damages. Therefore, please read and closely follow these instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.

Scope

M.flow Instruments include the following components:

- ▶ *M.flow Compass*
- ▶ *M.blue plus Adjustment Ring*

Additional options:

- ▶ *M.blue plus Adjustment Assistant*
- ▶ *M.flow Checkmate*

2 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE

2.1 EXPLANATION OF THE WARNINGS

**DANGER**

Indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.

**WARNING**

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

**CAUTION**

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE

Indicates information considered important, but not hazard-related (e.g. messages relating to property damage).

The symbols for danger, warning and caution are yellow warning triangles with a black edge and a black exclamation mark.

2.2 DISPLAY CONVENTIONS

Display	Description
<i>Italics</i>	Indicates <i>product names</i>

2.3 OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION

The latest version of the instructions for use as well as translations into additional languages can be found on our website:
<https://miethke.com/downloads/>
If you still need additional information despite carefully reading the instructions for use and the additional information, please contact us or your authorised distributor.

2.4 FEEDBACK ON THE INSTRUCTIONS FOR USE

Your opinion is important to us. Please let us know if you have any requests and criticisms about these instructions for use. We will analyse your feedback and take it into account for the next version of the instructions for use where appropriate.

2.5 COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION

Christoph Miethke GmbH & Co. KG guarantees a faultless product that is free of material and manufacturing defects upon delivery.
No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in this document, if it is combined with products by another manufacturer or if it is used in any way other than for the intended purpose and the intended use.
Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights.

3 DESCRIPTION *M.flow Instruments*

3.1 MEDICAL PURPOSE

The *M.flow Instruments* are intended for adjusting and verifying the adjustable flow-reduction levels of the flow-reducing device *M.flow*, as well as for adjusting and verifying the adjustable pressure levels of the valve *M.blue*.

3.2 CLINICAL BENEFITS

The *M.flow Instruments* are auxiliary instruments with the indirect clinical benefit of optimising hydrocephalus shunt therapy by enabling the valve adjustment function. They must only be used together with the respective adjustable hydrocephalus shunt. Accordingly, they are to be regarded as accessories to the MIETHKE shunt system.

3.3 INDICATIONS

The following indications apply to *M.flow Instruments*:

- ▶ Hydrocephalus

3.4 CONTRAINDICATIONS

The following contraindications apply to *M.flow Instruments*:

- ▶ Intolerance to materials of the instruments

3.5 INTENDED PATIENT GROUPS

- ▶ Patients who are treated with a CSF shunt system on account of their clinical picture

3.6 INTENDED USERS

In order to avoid risks due to false diagnoses, incorrect treatments and delays, the product must only be used by users with the following qualifications:

- ▶ Medical professionals, e. g. neurosurgeons
- ▶ Knowledge of mode of operation and intended use of the product. These can be acquired, e.g., in training offers organised by the Christoph Miethke GmbH & Co. KG (see chap. 5).

3.7 INTENDED USE ENVIRONMENT

Professional Healthcare Facilities

- ▶ Medical facilities, such as hospitals, doctor's surgeries and/or rehabilitation facilities

3.8 TECHNICAL DESCRIPTION

M.flow Instruments can be used to determine, change and monitor the flow-reduction level of the adjustable flow-reducing unit of the *M.flow* valve as well as the pressure level of the adjustable gravitational unit of the *M.blue* valve. The *M.flow* (Fig. 1) is a valve with an adjustable flow-reducing unit and an additional shunt switch-off function that is used for shunting cerebrospinal fluid (CSF) in the treatment of hydrocephalus. The flow-reducing unit allows the flow resistance to be adjusted and thus the shunted cerebrospinal fluid volume per unit of time to be controlled. The *M.flow*'s shunt switch-off function allows cerebrospinal fluid shunting to be fully halted.

Fig. 1: *M.flow*

The *M.blue* (Fig. 2) is a posture-dependent valve for the shunting of cerebrospinal fluid (CSF) in the treatment of hydrocephalus. It consists of a differential pressure unit and an adjustable gravitational unit.

Fig. 2: *M.blue*

The *M.flow Compass* (Fig. 3) is used for localising and reading the adjustable flow-reducing unit of the *M.flow* as well as the adjustable gravitational unit of the *M.blue*.


Fig. 3: M.flow Compass
a) open, 1. Scale ring
b) closed, 2. Float gauge

The *M.blue plus* Adjustment Ring (Fig. 4) can be used to set the flow-reduction level of the *M.flow* from 0 to 10 as well as the shunt switch-off range 'M.flow off' (⊗).
The *M.blue* Adjustment Ring can also be used to set the opening pressure of the gravitational unit of the *M.blue* from 0 to 40 cmH₂O.

Fig. 4: M.blue plus Adjustment Ring

The flow-reduction level for the adjustable flow-reducing unit of the *M.flow* can be changed before or after implantation. It is pre-set to 0 (minimum flow-reducing effect) by the manufacturer.
The opening pressure of the adjustable gravitational unit of the *M.blue* can be changed before or after implantation. It is pre-set by the manufacturer to 20 cmH₂O.

3.9 SYSTEM COMPONENTS

The *M.flow Instruments* can be used to set the following valves:

Product name	Figure
<i>M.flow</i>	
<i>M.blue</i>	

3.10 FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

The product *M.flow Instruments* is MR Unsafe.

CAUTION

The *M.flow Instruments* must not be taken into an MRI facility as this may pose a safety risk to the patient and/or the user.

4 PROPERTIES M.flow Instruments

4.1 PRODUCT DESCRIPTION

4.1.1 APPROVED COMPONENTS

The *M.flow Instruments* comprise the following components:

- *M.flow* Compass (1) and
- *M.blue plus* Adjustment Ring (2)

The *M.flow Instruments* can be supplemented by the components:

- *M.blue plus* Adjustment Assistant (3) and
- *M.flow* Checkmate (4)

	Item no.	Figure
1	FX901T	
2	FX892T	
3	FX893T	
4	FX904T	

4.1.2 SCOPE OF DELIVERY

Box content (FX900T)	Number
<i>M.flow Compass</i>	1
<i>M.blue plus Adjustment Ring</i>	1
Instructions for use for <i>M.flow Instruments</i>	1

The components listed in Ch. 4.1.1 can also be ordered individually.

4.1.3 STERILITY

WARNING

- ▶ If the sterile packaging of the *Checkmate*, or the product itself, is damaged, the product must not be used.
- ▶ The *Checkmate* must not be used after the expiry date and must be re-sterilised as described in chapter 4.1.4 .

The *Checkmate* is sterilised with steam under strictly controlled conditions. The respective expiry date is printed on the packaging.

With the exception of the *M.flow Checkmate*, the *M.flow Instruments* cannot be sterilised, only disinfected.

4.1.4 REPEATED USE AND RESTERILISATION

The *M.flow Checkmate* must be sterilised in double sterile packaging in an autoclave (steam sterilisation, fractionated vacuum method) at 134 °C and 5 minutes retention time.

4.1.5 PRODUCT LIFETIME

These medical devices are designed to operate accurately and reliably over long periods of time. The expected product life of *M.flow Instruments* is three years after their first use, provided that the product is exposed to normal usage conditions and is properly maintained (see chapter 4.5 Use of the product).

The instruments may be used beyond this time if they are in perfect working condition. However, we cannot guarantee that these medical devices will require replacement for medical or technical reasons.

4.1.6 PRODUCT CONFORMITY

The product meets regulatory requirements.

4.2 IMPORTANT SAFETY INFORMATION

4.2.1 SAFETY INSTRUCTIONS

Important! Read all safety instructions carefully before using the product. Follow the safety instructions in order to avoid injuries and life-threatening situations.

WARNING

- ▶ The products must not be used if the packaging or the products are damaged.
- ▶ The *Checkmate* must not be used after the expiry date and must be re-sterilised as described in chapter 4.1.4 .
- ▶ Due to the risk of injury resulting from incorrect use of the product, the instructions for use must be carefully read and understood before the product is used for the first time.
- ▶ Prior to use, it is essential to check the product for completeness and integrity.

4.2.2 COMPLICATIONS, SIDE EFFECTS, PRECAUTIONS AND RESIDUAL RISKS

The following complications and side effects can occur in conjunction with the product *M.flow Instruments*:

- ▶ Allergic reaction/intolerance to the materials of the instruments

As a precaution, a physician must be consulted immediately if the patient suffers from skin rashes and tightness, severe headaches, dizzy spells or similar.

4.2.3 REPORTING OBLIGATION

All serious incidents (damage, injuries, infections, etc.) occurring in relation to the product must be reported to the manufacturer and the responsible state authority.

4.3 PATIENT EDUCATION

The attending physician is responsible for informing the patient and/or his/her proxy in advance. The patient is to be informed about warnings, precautions, contraindications, precautionary measures to be taken as well as restrictions on use in relation to the product. (see Ch. 4.1.5, 4.2, 4.5, 4.8).

4.4 TRANSPORT AND STORAGE

The medical devices must always be transported and stored in a clean and dry place.

4.4.1 TRANSPORT

Transport conditions

Temperature range for transport	-30 °C ... +40 °C
---------------------------------	-------------------

4.4.2 STORAGE

Storage conditions

Temperature range for storage	0 °C ... +40 °C for <i>Checkmate</i> room temperature
-------------------------------	---

4.5 USE OF THE PRODUCT

4.5.1 INTRODUCTION

The *M.flow Instruments* can be used to determine, change and monitor the flow-reduction level of the adjustable flow-reducing unit of the *M.flow* valve as well as the pressure level of the adjustable gravitational unit of the *M.blue* valve.

4.5.2 SAFETY NOTICES AND WARNINGS

WARNING

- ▶ The position '*M.flow* off' (⊗), represented by a circle with a cross, completely stops the flow through the valve. Danger of under-drainage!
- ▶ Because of the magnets inside the *M.flow Instruments*, *M.flow Instruments* must not be used in the vicinity of active implants that are influenced by magnetism, such as cardiac pacemakers. When adjusting with the *M.flow Instruments*, the safety distance specified by the manufacturer of the active implant should be observed.

CAUTION

- ▶ The *M.flow Instruments* must not be taken into an MRI facility as this may pose a safety risk to the patient and/or the user.
- ▶ Only *M.flow Instruments* approved for that purpose must be used to determine, change and monitor the flow-reduction level of the *M.flow*.

CAUTION

- ▶ The opening pressure of the gravitational unit of *M.blue* can be determined, changed and monitored with the *M.flow Instruments* or with the *M.blue plus Instruments* approved for that purpose.

NOTICE

The *M.blue plus Adjustment Ring* emits a magnetic field. Metallic objects and magnetic storage media should be placed at a sufficient safety distance.

4.5.3 REQUIRED MATERIALS

Adjustment of the adjustable flow-reduction level of the *M.flow* or the adjustable gravitational unit of the *M.blue* after implantation requires the *M.flow Compass*, the *M.blue plus Adjustment Ring* and possibly the *M.blue plus Adjustment Assistant*.

To adjust the adjustable flow-reduction level of the *M.flow* before or during valve implantation, the re-sterilisable *M.flow Checkmate* can be used.

The product *M.flow Instruments* is designed so that it can be safely used with the shunt components described in chapter 3.9.

4.5.4 PREPARING FOR USE

Checking the sterile packaging

Immediately before using the *Checkmate*, the sterile packaging must be visually inspected in order to check the integrity of the sterile barrier system. The products should only be removed from the packaging immediately prior to use.

Checking the integrity of the *M.flow Instruments*

Prior to use, the instruments must be checked for integrity. To this end, all instruments must be visually inspected.

4.5.5 APPLICATION OF *M.flow Instruments*

The following steps must be performed in order to set the flow-reduction level of the *M.flow* or the pressure level of the gravitational unit of the *M.blue*:

1. Localisation

CAUTION

- ▶ The *M.flow Compass* reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the *M.blue plus Adjustment Ring* should not be placed in the immediate vicinity of the *M.flow Compass* when determining the flow-reduction level or opening pressure. We recommend a minimum distance of 30 cm.
- ▶ Swelling of the skin may make adjustment difficult for a few days after surgery. If the valve setting cannot be checked conclusively using the *M.flow Compass*, we recommend checking it with an imaging method.
- ▶ The *M.flow Compass* should be placed as centrally as possible over the valve, as the determined flow-reduction level or pressure level may otherwise be incorrect.
- ▶ Possible air pockets in the *M.flow Compass* do not affect the compass function.

When the *M.flow Compass* is opened, a circular cut-out becomes visible, which is used to localise the valve in the patient's head as centrally as possible using the index finger (Fig. 5).

Fig. 5: Localising the valve

The directional markings (proximal and distal) must be aligned according to the inlet and outlet connector of the valve. The directional markings show the flow direction.

2. Verification procedure

In order to determine the selected flow-reduction level or pressure level, the *M.flow Compass* is then closed again. The float gauge should now be centred within the designated circular marking by moving the *M.flow Compass*

(Fig. 6). Once the float gauge is centred, the currently set flow-reduction level of the *M.flow* or the currently set pressure level of the gravitational unit of the *M.blue* can be read off via the line mark on the float gauge (Fig. 7).


Fig. 6: Determining the flow-reduction level or pressure level with the *M.flow Compass*

The scale ring is marked with two scales (Fig. 7). The grey marked adjustment range from 0 (lowest flow-reduction level) to 10 (highest flow-reduction level) on the outer scale applies to the flow-reduction level of the *M.flow*. The setting '*M.flow off*' (⊗) (shunt switch-off range) causes the flow to be blocked. No cerebrospinal fluid shunting is possible with this setting. The adjustment range from 0 to 40 cmH₂O marked '*M.blue*' on the inner scale applies to the pressure level of the *M.blue* gravitational unit.


Fig. 7: M.flow Compass scale ring

1. Outer: Scale for the M.flow from 0 to 10 with additional shunt switch-off range 'M.flow off' (⊗) (flow-reduction level in the illustrated example: 0).

2. Inner: Scale for the M.blue gravitational unit from 0 to 40 cmH₂O (the opening pressure of the M.blue gravitational unit in the illustrated example: 28 cmH₂O)

3. Adjustment process

WARNING

The position 'M.flow off' (⊗), represented by a circle with a cross, completely stops the flow through the valve. Danger of underdrainage!

CAUTION

Between flow-reduction level 0 (minimal flow-reducing effect) and the setting 'M.flow off' (⊗) (valve impermeable/cerebrospinal fluid shunting prevented) there is a stop that prevents the extreme positions from being accidentally set incorrectly.

When adjusting, care must be taken to change the valve by a maximum of 4 levels per adjustment process; otherwise, errors may occur.

Example: The flow-reduction level is to be changed from 0 to 8. The correct method is an adjustment in two stages: an initial adjustment from 0 to 4 and subsequently from 4 to 8.

Example: The flow-reduction level is to be changed from 'M.flow off' (⊗) to 0. The correct method is an adjustment in three stages: an initial adjustment from 'M.flow off' (⊗) to 8, from 8 to 4 and subsequently from 4 to 0.

CAUTION

When adjusting the gravitational unit of the M.blue, care must be taken to change the opening pressure by a maximum of 16 cmH₂O per adjustment process; otherwise, errors may occur.

Example: The opening pressure is to be changed from 6 to 36 cmH₂O. The correct method is an adjustment in two stages: an initial adjustment from 6 to 22 cmH₂O and subsequently from 22 to 36 cmH₂O.

3a. Adjusting with the M.blue plus Adjustment Ring

In order to adjust the flow-reduction level of the M.flow or the pressure level of the M.blue, the M.flow Compass is opened, but without changing the position of the scale ring. The M.blue plus Adjustment Ring is now inserted into the scale ring in such a manner that the line marking points to the desired value on the scale of the scale ring. (Fig. 8)


Fig. 8: Inserting the M.blue plus Adjustment Ring:

1. M.blue plus Adjustment Ring, 2. Scale ring


Fig. 9: Setting the flow-reduction level of the M.flow in the example illustrated to 3 or setting the gravitational unit of the M.blue in the example illustrated to 40 cmH₂O.

By applying slight pressure with the index finger to the valve membrane located in the centre of the M.blue plus Adjustment Ring under the skin, the rotor brake is released and the flow-reduction level of the M.flow or the pressure level of the gravitational unit of the M.blue is changed to the desired value (Fig. 10).


Fig. 10: Adjustment using the M.blue plus Adjustment Ring

Both the flow-reducing unit of the M.flow and the gravitational unit of the M.blue are equipped with a feedback mechanism.

Due to the valve housing design, targeted pressure exerted on the valve membrane produces an audible acoustic signal (a clicking sound) and/or palpable resistance as soon as the rotor brake has been released. The valve thus shows both acoustically and haptically when the pressure is sufficient for uncoupling. Once this pres-

sure has been released, the rotor is once again adjustment-proof. Although the click caused by releasing the rotor brake is easily audible before implantation, it may be considerably reduced after implantation and the filling of the valve depending on its position and the condition of the implant surroundings. Normally, however, it should be audible to the patient or by using a stethoscope.

3b. Adjustment using the M.blue plus Adjustment Assistant

Alternatively, the M.blue plus Adjustment Assistant can be used to adjust the flow-reduction level or the pressure level. To do this, insert the M.blue plus Adjustment Assistant into the M.blue plus Adjustment Ring aligned to the desired value and press it with your index finger (Fig. 11).


Fig. 11: M.blue plus Adjustment Assistant inserted in the M.blue plus Adjustment Ring

4. Checking after adjustment

After adjusting the flow-reduction level or pressure level, it is recommended to check the setting. To do this, proceed as in Points 1 and 2. Should the measured value not correspond with the required flow-reduction level or pressure level, the adjustment process should be repeated. To do this, start again at Point 3.

5. Checking and adjusting using the M.flow Checkmate C€₀₂₉₇

The M.flow Checkmate (Fig. 12) is supplied sterile and can be re-sterilised. The M.flow Checkmate can be used to change and check the flow-reduction level adjustment before and during valve implantation directly on the M.flow. To determine the flow-reduction level, the

M.flow Checkmate is centred over the *M.flow*. The *M.flow Checkmate* automatically aligns itself over the valve. The flow-reduction level can be read from the direction of the proximal catheter (at the inlet connector). If the flow-reduction level is to be adjusted, the *M.flow Checkmate* is placed centrally onto the *M.flow*. When doing so, the desired flow-reduction level must point in the direction of the proximal catheter (leading towards the valve). By slightly pressing the *M.flow Checkmate* onto the valve, the rotor brake in the *M.flow* is released and the flow-reduction level set.

When adjusting, make sure that the valve is changed by a maximum of 4 flow-reduction levels per adjustment; otherwise, errors may occur (see Point '3. Adjustment process').

Fig. 12: *M.flow Checkmate*

4.5.6 CLEANING RECOMMENDATION FOR NON-STERILISABLE *M.flow Instruments*

CAUTION

***M.flow Instruments* are made from thermolabile components that are sensitive to heat and moisture as well as chemically reactive. Do not immerse *M.flow Instruments* in cleaning solutions or allow liquids to enter the housing, as moisture, corrosion and impurities could impair proper functioning.**

Surface contamination on the *M.flow Instruments* should always be removed immediately after use by wiping with alcohol-based cleaning agents (minimum 75% alcohol).

The contact time should be at least 60 seconds and should be adapted to the degree of contamination. The instruments should then be wiped with a dry cloth.

The following cleaning methods are not suitable for processing the *M.flow Instruments* (except for the *Checkmate*): Irradiation, ultrasound, sterilisation, mechanical processing, immersion in cleaning fluids.

4.5.7 CLEANING AND DISINFECTION OF THE *Checkmate*

Avoid damaging the product due to unsuitable cleaning agents/disinfectants and/or too high temperatures!

- ▶ Use cleaning agents/disinfectants approved for surgical steels in accordance with the manufacturer's instructions.
- ▶ Observe specifications regarding concentration, temperature and contact time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 55°C.
- ▶ Perform ultrasonic cleaning:
 - ▶ as an effective mechanical support for manual cleaning/disinfection
 - ▶ for precleaning products with dried residues prior to mechanical cleaning/disinfection
 - ▶ as integral mechanical support during mechanical cleaning/disinfection
 - ▶ for the subsequent cleaning of products with residues still present after mechanical cleaning/disinfection
- ▶ Provided the instruments can be fastened securely or safely for cleaning in machines or on storage aids, clean and disinfect the instruments mechanically.

Manual cleaning/disinfection

- Check visible surfaces for residues after manual cleaning/disinfection.
- If necessary, repeat the cleaning process.

Mechanical cleaning/disinfection

- Place product on suitable screen basket.

Mechanical cleaning/disinfection with manual precleaning

- Manual precleaning with ultrasound

Phase	Step	T (°C/°F)	t (min)	Conc. (%)	Water quality	Chemical specifications
I	Disinfected ultrasonic cleaning	RT (cold)	15	2	D-W	B. Braun Stabimed®; free of aldehyde, phenol and quaternary ammonium compounds; pH = 9
II	Intermediate rinsing	RT (cold)	1	-	D-W	

D-W: drinking water; RT: room temperature

Mechanical alkaline cleaning and thermal disinfection

Phase I

- Clean product in the ultrasonic cleaning bath (frequency 35 kHz). Ensure that all accessible surfaces are wetted.

Phase II

- Fully rinse/flush product (all accessible surfaces) under running water.

Device type: Single chamber cleaning/disinfection device without ultrasound

- Place product on suitable screen basket (avoid missed areas)

Phase	Step	T (°C/°F)	t (min)	Water quality	Chemical specifications
I	Pre-rinsing	< 25/77	3	D-W	-
II	Cleaning	55/131	10	P-W	- Concentration, alkaline: pH = 10.9; < 5 % anionic surfactants 1% solution pH = 10.5
III	Neutralisation	20/68	2	P-W	- Concentration, acidic: pH = 2.6; Base: Citric acid - 1% solution pH = 3.0

D-W: Drinking water, P-W: Purified water demineralised

Phase	Step	T (°C/°F)	t (min)	Water quality	Chemical specifications
IV	Inter-mediate inspection	70/158	1	P-W	-
V	Thermal disinfection	94/201	10	P-W	-
VI	Drying	90/194	40	-	-

D-W: Drinking water, P-W: Purified water demineralised

Checking, maintenance, testing

- ▶ Allow product to cool down to room temperature.
- ▶ Check product after each cleaning and disinfection for: Cleanliness, functioning and damage, e.g. insulation, loose, bent, broken, cracked, worn and detached parts.
- ▶ Immediately discard damaged products. After the cleaning process, it is recommended to place the instrument inside suitable double sterile packaging. This will ensure reuse in the best way possible.

4.6 DISPOSAL

Adhere to national regulations when disposing of or recycling the product, its components or its packaging.

The products are to be disposed of correctly as potentially infectious material in accordance with medical practice as well as the respective regional laws and regulations.

4.7 TROUBLESHOOTING

- ▶ The *M.flow Compass* should be placed as centrally as possible over the valve, as the determined flow-reduction level or pressure level may otherwise be incorrect.
- ▶ It must be ensured that the arrow markings on the compass housing (proximal, distal) correspond with the flow direction in the shunt system.
- ▶ To ensure a correct reading of the set pressure level or flow-reduction level, the float gauge in the *M.flow Compass* should be centred by slightly moving the instrument within the designated circular marking.
- ▶ The *M.flow Compass* is marked with two scale rings:
 - ▶ For the flow-reduction level of the *M.flow*, the adjustment range marked grey from 0 (lowest flow-reduction level) to 10 (highest flow-reduction level) as well as the shunt switch-off range (⊗) on the outer scale applies.
 - ▶ For the gravitational unit of the *M.blue*, the pressure level setting range marked with 'M.blue' from 0 to 40 cmH₂O on the inner scale applies.
- ▶ When adjusting the *M.flow*, care must be taken to change the valve by a maximum of 4 levels per adjustment process; otherwise, errors may occur.
Example: The flow-reduction level must be changed from 0 to 8. The correct method is an adjustment in two stages: An initial adjustment from 0 to 4 and subsequently from 4 to 8.
Example: The flow-reduction level must be changed from 'M.flow off' (⊗) to 0. The correct method is an adjustment in three stages: An initial adjustment from 'M.flow off' (⊗) to 8, from 8 to 4 and subsequently from 4 to 0.
- ▶ When adjusting the gravitational unit of the *M.blue*, care must be taken to change the opening pressure by a maximum of

16 cmH₂O per adjustment process; otherwise, errors may occur.

Example: The opening pressure must be changed from 6 to 36 cmH₂O. The correct method is an adjustment in two stages: an initial adjustment from 6 to 22 cmH₂O and subsequently from 22 to 36 cmH₂O.

- ▶ The *M.flow Compass* reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the *M.blue plus Adjustment Ring* should not be placed in the immediate vicinity of the *M.flow Compass* when determining the flow-reduction level or pressure level. We recommend a minimum distance of 30 cm.
- ▶ Swelling of the skin may make adjustment difficult for a few days after surgery. If

the valve setting cannot be checked conclusively using the *M.flow Compass*, we recommend checking it with an imaging method.

- ▶ When implanting the *M.flow* and/or *M.blue*, the distance between the surface of the valve and the skin should not exceed 10 mm. If the skin and tissue thickness exceeds 10 mm, it may be difficult or impossible to locate the valve and read and set the flow-reduction level or pressure level.
- ▶ The adjustable MIETHKE valve should not be implanted in an area that makes the detection or palpation of the valve difficult (e. g. underneath heavily scarred tissue).

4.8 TECHNICAL INFORMATION

4.8.1 MATERIALS THAT COME INTO CONTACT WITH BODILY TISSUE/FLUIDS

<i>M.flow Compass</i>	Aluminium (AlMgSi0,5, AlMgSi1)
<i>M.blue plus Adjustment Ring</i>	Aluminium (AlMgSi0,5, AlMgSi1)
<i>M.blue plus Adjustment Assistant</i>	Aluminium (AlMgSi0,5, AlMgSi1), polyurethane
<i>M.flow Checkmate</i>	Titanium alloy (TiAl6V4)

4.8.2 TECHNICAL DATA

Manufacturer	Christoph Miethke GmbH & Co. KG
Product designation	<i>M.flow Instruments</i>
Medical purpose	Adjusting and checking the adjustable flow-reduction levels of the <i>M.flow</i> valve as well as for adjusting and checking the adjustable pressure levels of the <i>M.blue</i> valve
Suitability for sterilisation/disinfection	<i>M.flow Checkmate</i> suitable for re-sterilisation <i>M.flow Compass</i> suitable for disinfection <i>M.blue plus Adjustment Ring</i> suitable for disinfection <i>M.blue plus Adjustment Assistant</i> suitable for disinfection
Storage	Store in a clean and dry place
Dimensions	<i>M.flow Compass</i> : ø 58.7 mm, H: 16.4 mm <i>M.blue plus Adjustment Ring</i> : ø 43 mm; H: 11.5 mm <i>M.blue plus Adjustment Assistant</i> : ø 23.5 mm; H: 28.5 mm <i>M.flow Checkmate</i> : ø 20 mm; H: 27.5 mm

4.9 SYMBOLS USED FOR LABELLING

Symbol	Explanation
	EU conformity marking, xxxx indicates the identifier of the responsible notified body.
	Medical device
	Manufacturer
	Date of manufacture
	Use-by date
	Batch code
	Catalogue number
	Quantity
	Unique Device Identifier
UDI-DI	UDI-DI number
	Sterilised using steam
	Double sterile barrier system
	Do not use if package is damaged and consult instructions for use.
	Keep dry
	Consult instructions for use / electronic instructions for use
	Caution
	Non-pyrogenic
	Free of natural rubber latex, latex-free
Rx _{only}	Indicates that in the USA, the product may only be issued to physicians.
	MR Unsafe

Symbol	Explanation
	Product contains strong magnets

5 MEDICAL DEVICE CONSULTANTS

In compliance with regulatory requirements, Christoph Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions.

You can contact our medical device consultants at:

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany