


SAXXONY®

Posterior Cervical Thoracic System



Connector Options

Variety of adjustable rod-to-rod, head-to-head, offset, and cross connectors for surgeon need.




System Compatibility

Pair with Inertia® Connexx™ using Ø3.5-5.5mm transition rods and connectors.



Screw Options

Shanks available in Ø3.5-5.0mm and with full or partial threads.



Standard Set Contents



Full Thread Screws*

Ø3.5 x 10mm	8
Ø3.5 x 12mm	8
Ø3.5 x 14mm	8
Ø3.5 x 16mm	8

Ø4.0 x 10mm	2
Ø4.0 x 12mm	2
Ø4.0 x 14mm	2
Ø4.0 x 16mm	2
Ø4.0 x 18mm	4
Ø4.0 x 20mm	4
Ø4.0 x 22mm	4
Ø4.0 x 24mm	4
Ø4.0 x 26mm	4
Ø4.0 x 28mm	2
Ø4.0 x 30mm	2
Ø4.0 x 35mm	2
Ø4.0 x 40mm	2

Prebent Rods

Ø3.5 x 20mm	3
Ø3.5 x 25mm	3
Ø3.5 x 30mm	3
Ø3.5 x 35mm	3
Ø3.5 x 40mm	3
Ø3.5 x 50mm	3
Ø3.5 x 60mm	3
Ø3.5 x 70mm	3
Ø3.5 x 80mm	3

Straight Rods

Ø3.5 x 120mm	3
Ø3.5 x 240mm	3

Transition Rod QTY

Ø3.5-5.5x400mm	3
----------------	---

Set Screws QTY

T20 Set Screw	20
---------------	----

Head to Head Conn.


H2H 25-30mm	2
H2H 30-35mm	2
H2H 35-40mm	2
H2H 40-45mm	2
H2H 45-50mm	2
H2H 50-55mm	2
H2H Set Screw	8
H2H Lock Nut	8

Cross Connectors

R2R 20-32mm	2
R2R 32-38mm	2
R2R 38-44mm	2
R2R 40-52mm	2

*Screws available
Ø3.5-5.0mm
by request

Optional Sets



Highlighted Systems



Struxxure® MCS
Single Level Plates



Struxxure® -C
Single & Multi-Level Plates



NEXXT MATRIXX® Cervical
Standalone, Standard, Convexx™

Nexxt Spine, LLC
14425 Bergen Blvd, Suite B
Noblesville, IN 46060
(317) 436-7801
Info@NexxtSpine.com
NexxtSpine.com
70-054, Rev C



Contact customer service at Info@NexxtSpine.com for full SKU offering. For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert or contact your local representative; visit NexxtSpine.com for additional product information. All rights reserved. All content herein is protected by copyright, trademarks and other intellectual property rights owned by Nexxt Spine, LLC and must not be redistributed, duplicated or disclosed, in whole or in part, without the expressed written consent of Nexxt Spine, LLC. This material is intended for healthcare professionals, the Nexxt Spine sales force and authorized representatives. Distribution to any other recipient is prohibited.