

INERTIA®

Pedicle Screw System

- Degenerative Vertebral Body Fusion
- Deformity Correction System




Nexxt Spine, LLC
14425 Bergen Blvd, Suite B
Noblesville, IN 46060
(317) 436-7801
Info@NexxtSpine.com

70-060, Rev C

HOUSING / SCREW OPTIONS

Screws available $\varnothing 4.5-10.5\text{mm}$ x 20-120mm long.

Some implants offered By Request. Contact Info@NexxtSpine.com for implant availability.



ROD FEATURES

$\varnothing 5.5\text{mm}$ - Ti-6Al-4V Alloy

Optional Hex End and Laser Lined

Straight rods are available up to 600mm

Curved, pre-cut rods are available up to 175mm

CONNECTOR OPTIONS


Ti-6Al-4V ELI per ASTM F136



Single Open
Wedding Band



Double Open
Wedding Band




Double
Closed Domino



Double Closed
Wedding Band



Inline
2 Screw



Cross Connectors
23-94mm



Offset Connectors
25-60mm

HIGHLIGHTED INSTRUMENTS

Pistol Grip Reducer
20mm of Reduction



Turbo Tower Reducer
25mm of Reduction



Screw Inserter
Auto-Locking



COMPATIBLE NEXXT MATRIX® SYSTEMS



Lateral + STRUXXURE® L



ALIF + STRUXXURE® A



SA ALIF



Corpectomy



TLIF



TLIF Oblique



Nexxt Spine, LLC
14425 Bergen Blvd, Suite B
Noblesville, IN 46060
(317) 436-7801
Info@NexxtSpine.com
NexxtSpine.com



For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert or contact your local representative; visit NexxtSpine.com for additional product information.

All rights reserved. All content herein is protected by copyright, trademarks and other intellectual property rights owned by Nexxt Spine, LLC and must not be redistributed, duplicated or disclosed, in whole or in part, without the expressed written consent of Nexxt Spine, LLC. This material is intended for healthcare professionals, the Nexxt Spine sales force, and authorized representatives. Distribution to any other recipient is prohibited.