


CastleLoc-C®



ANTERIOR CERVICAL INTERBODY FUSION CAGE SYSTEM

CastleLoc is a Registered Trademark.



CastleLoc-C Anterior Cervical Interbody Fusion Cage System

Overview


The *CastleLoc-C* Cervical Interbody Fusion Cage System is indicated for Anterior Cervical Interbody Fusion procedures in skeletally mature patients with cervical disc disease at one level from the C3 disc to the C7 disc.

Features & Benefits

The *CastleLoc-C* Cervical Interbody Fusion Cage System is available in a variety of different sizes and configurations to accommodate anatomical variation in different vertebral levels and/or patient anatomy.


Design

Anatomical shape of Wedge Cage which allows the surgeon to choose the size suited to the patient's anatomy.



Radiopaque Markers

Allow visualization of cage in the disc space in radiographic image. Marker Pins made of Unalloyed Tantalum for better radiographic image without artifact.



Simple instruments for easier and faster surgery



Indications

CastleLoc-C Cervical Interbody Fusion Cage System is indicated for use in skeletally mature patients with degenerative disc disease (DDD) of the cervical spine with accompanying radicular symptoms at one disc level. DDD is defined as disogenic pain with degeneration of the disc confirmed by patient history and radiographic studies. CastleLoc-C Cervical Interbody Fusion Cage System is used to facilitate Interbody Fusion in the cervical spine at the C3 to C7 disc levels using autograft bone. CastleLoc-C Cervical Interbody Fusion Cage System is to be used with supplemental fixation. Patients should have at least six (6) weeks of non-operative treatment prior to treatment with an intervertebral cage.

Ordering Information

Ordering Information – Cage

Wedge Type



Cervical Cage, 7°, Wedge

CAT. NO.	HEIGHT (mm)	WIDTH (mm)	DEPTH (mm)	ANGLE (°)
5205-1214	5	12	14	7
5206-1214	6	12	14	7
5207-1214	7	12	14	7
5208-1214	8	12	14	7
5209-1214	9	12	14	7
5210-1214	10	12	14	7
5205-1417	5	14	17	7
5206-1417	6	14	17	7
5207-1417	7	14	17	7
5208-1417	8	14	17	7
5209-1417	9	14	17	7
5210-1417	10	14	17	7

Ordering Information – Instrument

Cat No.	Description	QTY (PER SET)
CC02-4705	CERVICAL CAGE TRIAL 12X14X5 / 7°	1
CC02-4706	CERVICAL CAGE TRIAL 12X14X6 / 7°	1
CC02-4707	CERVICAL CAGE TRIAL 12X14X7 / 7°	1
CC02-4708	CERVICAL CAGE TRIAL 12X14X8 / 7°	1
CC02-4709	CERVICAL CAGE TRIAL 12X14X9 / 7°	1
CC02-4710	CERVICAL CAGE TRIAL 12X14X10 / 7°	1
CC02-5705	CERVICAL CAGE TRIAL 14X17X5 / 7°	1
CC02-5706	CERVICAL CAGE TRIAL 14X17X6 / 7°	1
CC02-5707	CERVICAL CAGE TRIAL 14X17X7 / 7°	1
CC02-5708	CERVICAL CAGE TRIAL 14X17X8 / 7°	1
CC02-5709	CERVICAL CAGE TRIAL 14X17X9 / 7°	1
CC02-5710	CERVICAL CAGE TRIAL 14X17X10 / 7°	1
CC02-1001	CAGE HOLDER	2
CC02-2001	BONE GRAFT HOLDER	1
CC02-2002	BONE GRAFT IMPACTOR	1
CC02-3001	CAGE FINAL IMPACTOR	1
CC02-0001	CONTAINER COVER	1
CC02-0002	CONTAINER BASE	1
CC02-0003	UPPER TRAY	1
CC02-0005	LOWER TRAY	1
CC02-0008	IMPLANT CADDY CASE COVER	1
CC02-0009	IMPLANT CADDY CASE	1

Option

Ordering Information – Cage

Convex Type



Cervical Cage, 0°, Convex

CAT. NO.	HEIGHT (mm)	WIDTH (mm)	DEPTH (mm)	ANGLE (°)
5105-1214	5	12	14	0
5106-1214	6	12	14	0
5107-1214	7	12	14	0
5108-1214	8	12	14	0
5109-1214	9	12	14	0
5110-1214	10	12	14	0
5105-1417	5	14	17	0
5106-1417	6	14	17	0
5107-1417	7	14	17	0
5108-1417	8	14	17	0
5109-1417	9	14	17	0
5110-1417	10	14	17	0

Ordering Information – Instrument

Cat No.	Description	QTY (PER SET)
CC02-6705	CERVICAL CAGE TRIAL 12X14X5 / 0°	1
CC02-6706	CERVICAL CAGE TRIAL 12X14X6 / 0°	1
CC02-6707	CERVICAL CAGE TRIAL 12X14X7 / 0°	1
CC02-6708	CERVICAL CAGE TRIAL 12X14X8 / 0°	1
CC02-6709	CERVICAL CAGE TRIAL 12X14X9 / 0°	1
CC02-6710	CERVICAL CAGE TRIAL 12X14X10 / 0°	1
CC02-7705	CERVICAL CAGE TRIAL 14X17X5 / 0°	1
CC02-7706	CERVICAL CAGE TRIAL 14X17X6 / 0°	1
CC02-7707	CERVICAL CAGE TRIAL 14X17X7 / 0°	1
CC02-7708	CERVICAL CAGE TRIAL 14X17X8 / 0°	1
CC02-7709	CERVICAL CAGE TRIAL 14X17X9 / 0°	1
CC02-7710	CERVICAL CAGE TRIAL 14X17X10 / 0°	1
CC02-1001	CAGE HOLDER	2
CC02-2001	BONE GRAFT HOLDER	1
CC02-2002	BONE GRAFT IMPACTOR	1
CC02-3001	CAGE FINAL IMPACTOR	1
CC02-0001	CONTAINER COVER	1
CC02-0002	CONTAINER BASE	1
CC02-0003	UPPER TRAY	1
CC02-0005	LOWER TRAY	1
CC02-0008	IMPLANT CADDY CASE COVER	1
CC02-0009	IMPLANT CADDY CASE	1

WORLD-LEADING ENTERPRISE
IN THE BIOTECH INDUSTRY



Head Office : #201, 16-25, Dongbaekjungang-ro, Giheung-gu, Yongin-si, Gyeonggi-do, 17015, Korea
Seoul Office : 17F, 159-1, Mokdongseo-ro, Yangcheon-gu, Seoul, 07997, Korea
Tel. 82-1600-0841 Fax. 82-1600-0843 www.lnkbimed.com



BR-01(32_1711)