


ENTASIS®

SACROILIAC JOINT FUSION



JOINT COMPRESSION

Dual-pitch threading provides up to 4.5mm of compression, designed to facilitate fusion across the joint

EASE OF INSERTION

Self-drilling, dual-lead threads allow for quick insertion

Stackable Guide Wires provide stability while streamlining technique

THE ENTASIS ADVANTAGE

The Entasis SI Joint Fusion System was designed to provide a faster and potentially safer SI joint procedure with joint compression that the surgeon can feel.

Self-drilling, dual-lead, and dual-pitch threads allow for quick insertion and compression across the joint. The system features rigid, stackable Guide Wires to allow constant proximal control and stability.

PRODUCT SPECIFICATIONS

Screws:


- Diameters: 7.0mm, 9.5mm, 11.5mm
- Lengths: 30mm – 70mm (5mm increments)
- Optional 13mm washer for 7.0mm implant

Guide Wires:

- Blunt
- Trocar Tip



FEATURES



For product information, including indications for use, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the product Instructions For Use (IFU) and information on the CoreLink website. To obtain a copy of the current Instructions for Use (IFU) for full prescribing and risk information, please call CoreLink Customer Service at (888) 349-7808 or visit corelinksurgical.com/ifu. © 2021 CoreLink, LLC. All rights reserved.



Designed, developed, and manufactured
in the United States of America.

St. Louis, MO 63026
CoreLinkSurgical.com

P: (888) 349-7808

CL-MKTG-0054, Rev. 3