


ParamountTM

Anterior Cervical Cage



Controlled
Plate
Deployment

Zero Profile

Less Invasive

SPINE*WAVE*



Controlled Anchor Plate Deployment

- Threaded deployment avoids aggressive impaction of Anchor Plates in the cervical spine
- Both Anchor Plates deploy simultaneously, simplifying the procedure

Zero Profile

- Anchor Plates provide engagement without instrumentation protruding anterior of the vertebral bodies




Less Invasive

- Self-contained Anchor Plates are deployed in the same plane as the disc space, minimizing the exposure

Pre-deployed



Post-deployed



SPINEWAVE

3 Enterprise Drive, Suite 210, Shelton, CT 06484, U.S.A
1-877-844-4225
www.spinewave.com