

prima**LOK**[™]**FF**

Facet Fixation System

The Solution for Posterior Facet Fixation.

primaLOK™ FF


Facet Fixation, Simplified

Minimally invasive posterior fixation achieved by simple placement of self-tapping screws with a polyaxial, anti-backout washer.

Design Rationale


primaLOK FF is a next-generation Facet Fixation device. Designed with a patented, polyaxial anti-migration washer that conforms to the patient anatomy, the implant minimizes tissue disruption and irritation. These unique features, combined with a simplified instrument kit, result in an easy-to-use device.

- Polyaxial washer for secure contact with the bone
- Streamlined instruments allow for easy device placement
- Complements the VariLift®-LX System



Broad Flexibility in a Small Footprint

With a full range of screw lengths, primaLOK FF provides ideal supplemental fixation for Anterior Lumbar Interbody Fusion (ALIF), Transforaminal Lumbar Interbody Fusion (TLIF), Lateral Lumbar Interbody Fusion (LLIF), and Posterior Lumbar Interbody Fusion (PLIF).





Wenzel Spine, Inc.

1130 Rutherford Lane • Suite 200 • Austin, TX 78753

BRO-0009C© • Wenzel Spine 2018

1-855-WENZEL1 • www.wenzelspine.com