



# Cantonale™ [MAX] TLIF Ti

TITANIUM TLIF CAGE SYSTEM

## System Features

Trabecular Structure for Multi-Axial New Bone In-Growth

Porous Titanium Providing Excellent Long-Term Fixation

Superior and Inferior Tooth Pattern Prevents Implant Migration


Multiple Lengths Accommodate Varying Patient Anatomy

Bullet Nose for Ease of Insertion


Large Graft Space for Adequate Grafting Material

Side Windows for Additional Bone Growth in the Axial Plane

Non-Lordotic



Lordotic



Footprint



### Cage Sizes (Offered in Non-Lordotic & 5° Lordotic)

| Heights (mm) | Lengths (mm) x Widths (mm)         |
|--------------|------------------------------------|
| 7 – 15       | 26 x 11, 28 x 11, 30 x 11, 32 x 11 |

[neurostructures.com](http://neurostructures.com)

199 Technology Drive, Suite 110 | Irvine, CA 92618 | Ph/Fax: 800.352.6103

All content herein is protected by copyright, trademarks, and other intellectual property rights owned by or licensed to NeuroStructures, Inc. or one of its affiliates, unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of NeuroStructures, Inc. This material is intended for the NeuroStructures, Inc. sales force and healthcare professionals. Distribution to any other recipient is prohibited.

© 2020 NeuroStructures, Inc. All Rights Reserved.



09/20 Rev. 001 NS-100-202